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Nazwa produktu: SARS-CoV-2 Antigen Rapid Test Kit (Immunochromatography)

Numer katalogowy: LX-401301

Producent: Labnovation Technologies, Inc.

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SARS-CoV-2 Antigen Rapid Test Kit (Immunochromatography)

INTENDED USE

This rapid test kit is intended for the qualitative detection of SARS-CoV-2 infection from patients. It is for professional use only. It is an aid in the diagnosis of the patients with suspected SARS-CoV-2 infection in conjunction with clinical presentation and results of other laboratory tests. Results from this test kit should not be used as the sole basis for diagnosis.

The test provides preliminary test results. Negative results do not preclude SARSCoV-2 infection and should not be used as the sole basis for treatment or other management decision.

PRINCIPLE

This kit is an immunochromatography assay. According to the gold immunochromatographic test principle, double antibody sandwich method was used to detect SARS-CoV-2 nucleocapsid antigen in the samples. When there is virus antigen presence in the sample, the antigen binds with the corresponding colloidal gold monoclonal antibody and the coated monoclonal antibody at the detection line to form a compound and then condenses into a red band, indicating a positive result. If there is no antigen in the sample, complex cannot be formed at the detection line, and no red band is shown, indicating negative result. Whether the sample contains antigen or not, the gold monoclonal antibody will bind to the enveloped antibody at the quality control line, form a compound and condense into a red band.

KIT COMPONENTS

- 20 Test cassettes
- 20 Sample tubes with prefilled sample extraction buffer
- 20 Single packaged Swabs (Nylon tip)
- 1 Positive control swab and 1 negative control swab
- 1 Tube stand
- 1 Instruction for use

STORAGE AND EXPIRY

Store as packaged in the sealed pouch at 2-30°C, avoid hot and sunshine, dry place, valid for 24 months. Do not freeze. Some protective measures should be taken in hot summer and cold winter to avoid high temperature or freeze-thaw. Do not open the inner packaging until ready, it must be used in one hour if opened (Humidity≤60%, Temp: 20°C-30°C). Please use immediately when the humidity > 60%.

SAMPLE PRESEVATION

Sample of human nasopharyngeal, oropharyngeal and anterior nasal swabs should be processed as soon as possible after sample collection. If the test cannot be performed immediately, the sample should be stored in a sealed state, stored at 2-8°C for 8 hours, and stored below -20°C for 1 month. Long-term storage is not recommended.

QUALITY CONTROL

A quality control should be performed prior to using the test kit. Use the positive and negative control swabs provided in accordance with the instructions " SAMPLE TREATMENT" and "TEST STEPS". If the quality control result is neither negative nor positive, the test kit should not be used.

SAMPLE COLLECTION


By nasopharyngeal swab: The sampler gently holds the head of the person to be collected with one hand, holds the swab with the other hand, sticks the swab to the nostril to enter, and slowly penetrates backwards along the bottom of the lower nasal passage, so as not to exert too much force to avoid traumatic hemorrhage. When the tip of the swab reaches the posterior wall of the nasopharynx cavity, gently rotate it for several times (in case of reflex cough, stop for a minute), and then slowly take out the swab.

By oropharyngeal swab: The head of the person to be collected is slightly tilted and his mouth is wide open, exposing the pharyngeal tonsils on both sides. Wipe the swab across the root of the tongue. Wipe the swab over the root of the tongue and then wipe the palatine tonsils on both sides at least three times with a little force. Finally, wipe the back of the throat wall up and down at least three times.

By anterior nasal swab: Gently insert the entire absorbent tip of the swab (around 1.5cm) into your nostril. Firmly and slowly rotate the swab in a circular path against the inside of your nostril at least 5 times. Be sure to collect any nasal drainage that maybe present on the swab. Gently take out the swab from the nostril. Use the same swab to repeat steps in the other nostril. Then slowly take out the swab.

SAMPLE TREATMENT


1. Take a sample tube (with prefilled sample extraction solution). Remove the aluminum foil sealing of the prefilled sample extraction tube and then put the tube into the tube stand.
2. Dip the swab into the sample extraction liquid after sample collection, fully immerse the tip of the swab, rotate and squeeze the swab 10 times, then take the stranded liquid as much as possible.
3. Break off the swab shaft from the breaking point. Discard the upper end of the swab.
4. Close the cap of the sample tube.

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1. Take out the test cassette from the sealed pouch, place it on a clean and level surface with the sample port well up.
2. Apply 3 full drops of the treated sample (60µl-70µl) vertically into the sample well of the test cassette.
3. Observe the test results immediately within 15-20 minutes, the result is invalid over 20 minutes.

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
INTERPRETATION OF RESULTS

POSITIVE: Two (2) distinct colored lines appear. One line should be in the control region (C) and the other line should be in the test region (T).

NEGATIVE: One (1) colored line appears in the control region(C). No apparent colored line appears in the test region (T). The negative result does not indicate the absence of analytes in the sample, it only indicates the level of tested analytes in the sample is less than the minimum detection limit.

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PERFORMANCE CHARACTERISTICS

1. Sensitivity and Specificity

For the clinical study, 102 oro- / nasopharyngeal samples and 204 anterior nasal samples within 7 days of onset of symptoms were collected. In addition, 400 oro- / nasopharyngeal and 200 anterior nasal samples which showed no symptoms and no specific exposure risk for SARS CoV-2, were collected to confirm the specificity. The verification is carried out with enterprise reference and the results meet the following requirements:

The National Reference Panel for SARS-Cov-2 Antigen Detection Kit of the National Institutes for Food and Drug Control was used for verification, and the results met the following requirements:

Coincidence rate of negative reference products: National negative reference materials (N1-N20) were tested, and the results were all negative, with a coincidence rate of 100%.

Coincidence rate of positive reference products: tested with national positive reference products (P1-P8), the results are all positive, and the coincidence rate is 100%.

Repeatability: The 10 test results of R1 and R2 should be positive, and the color rendering is uniform.

Clinical Sensitivity and Specificity of naso-oropharyngeal samples:

| SARS-CoV-2 Antigen Rapid Test | Nucleic Acid Test | | Total |
|-------------------------------|-------------------|----------|-------|
| | Positive | Negative | |
| Positive | 100 | 0 | 100 |
| Negative | 2 | 400 | 402 |
| TOTAL | 102 | 400 | 502 |

Diagnostic sensitivity: 98.04% (93.1%-99.76%)*
 Diagnostic specificity: 100% (99.08%-100.00%)*
 Total coincidence rate: 99.6% (98.57%-99.95%) *
 *95% CI (Confidence Interval)

| Test Results | Sensitivity | 95% CI (%) | |
|--------------|-------------|------------|-----|
| CT≤30 | 100% | 96.27 | 100 |
| CT≤27 | 100% | 95.98 | 100 |
| CT≤24 | 100% | 94.72 | 100 |

Clinical Sensitivity and Specificity of anterior nasal samples:

| SARS-CoV-2 Antigen Rapid Test | Nucleic Acid Test | | Total |
|-------------------------------|-------------------|----------|-------|
| | Positive | Negative | |
| Positive | 201 | 0 | 201 |
| Negative | 3 | 200 | 203 |
| TOTAL | 204 | 200 | 404 |

Diagnostic sensitivity: 98.5% (95.8% - 99.9%)*
 Diagnostic specificity: 100% (98.1% - 100%)*
 Total coincidence rate: 99.3% (97.8% - 99.8%)*
 *95% CI (Confidence Interval)

| Test Results | Sensitivity | 95% CI (%) | |
|--------------|-------------|------------|------|
| CT≤30 | 99.1% | 94.9 | 99.8 |
| CT≤32 | 99.4% | 96.5 | 99.9 |
| CT≤34 | 99.5% | 97.0 | 99.9 |

2. Limit of Detection

| | |
|-------------------|---------------------------|
| LOD concentration | 30 TCID ₅₀ /mL |
|-------------------|---------------------------|

3. Cross-reactivity

With human coronavirus 229E, human coronavirus OC43, human coronavirus NL63, adenovirus, human metapneumovirus, parainfluenza virus 1, parainfluenza virus 2, parainfluenza virus 3, parainfluenza virus 4, influenza A, type B Influenza, enterovirus, respiratory syncytial virus, rhinovirus, SARS coronavirus, MERS coronavirus, Haemophilus influenzae, Streptococcus pneumoniae, Streptococcus pyogenes, Candida albicans, Bordetella pertussis, Mycoplasma pneumonia,

pneumonia Chlamydia, Legionella pneumophila, etc. have no cross reaction.

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PRODUCENT

LABNOVATION TECHNOLOGIES, INC.

Adres: 101 and 5th Floor, Building 1, No. 68, 18th Road, Guangming Hi-Tech Park, Tangjia Community, Fenghuang Street, Guangming District, Shenzhen 518107, Guangdong, Chiny.

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

SARS-CoV-2 Antigen & Influenza A/B Rapid Test Kit

(Immunochromatography)

For professional and in vitro diagnostic use only.

PRODUCT NAME

SARS-CoV-2 Antigen & Influenza A/B Rapid Test Kit (Immunochromatography)

INTENDED USE

The SARS-CoV-2 Antigen & Influenza A/B Rapid Test Kit

(Immunochromatography) is a lateral flow chromatography immunoassay intended for the qualitative detection and differentiation of nucleocapsid protein antigens from SARS-CoV-2, influenza type A and influenza type B in the same time directly from nasal/nasopharyngeal (NS/NP) swab, oropharyngeal (OP) swab specimens from patients with signs and symptoms of respiratory viral infection. It is for professional use only. It is an aid in the diagnosis of the patients with suspected SARS-CoV-2 or influenza infection in conjunction with clinical presentation and results of other laboratory tests. Results from this test kit should not be used as the sole basis for diagnosis.

Negative results should be treated as presumptive and confirmed by viral culture or a molecular assay. Negative results do not rule out COVID-19 and influenza virus infection and should not be used as the sole basis for treatment or patient management decisions, including infection control decisions. Negative results should be considered in the context of a patient's recent exposures, history and the presence of clinical signs and symptoms consistent with COVID-19.

The kit is intended for use by trained clinical laboratory personnel and individuals trained in point of care settings.

For prescription use only. For in vitro diagnostic use only.

PRINCIPLE

The SARS-CoV-2 Antigen & Influenza A/B Rapid Test Kit employs a lateral flow chromatographic technology to detect the presence of the nucleocapsid protein antigen directly from influenza A, influenza B, and SARS-CoV-2.

After the patient sample is collected and treated with the antigen extraction buffer, the viral nucleoproteins antigen will be exposed. Add the extracted specimen into the test cassette, the specimen will migrate forward along with the test strips through capillary effect. If Influenza A, Influenza B, SARS-CoV or SARS-CoV-2 viral antigen is present, they will be captured and detected on the T or A/B line, respectively on each test strips within 15 minutes from the addition of the samples, resulting in purplish red band on the test region, indicating a positive result. If the nucleocapsid protein antigen is not present or present at very low levels in the

sample, there is no red line appears in "T" or "A/B" positions. The "Control Line" (C) is used for procedural control.

Control line should always appear if the test procedure is performed properly and the test reagents of control line are working.

KIT COMPONENTS

20 Test cassettes

2 Sample extraction

20 Sample tubes, 20 Swabs

1 Instruction Manual

STORAGE AND EXPIRY

Store as packaged in the sealed pouch at 2-30°C, avoid hot and sunshine, dry place, valid for 12 months. Do not freeze. Some protective measures should be taken in hot summer and cold winter to avoid high temperature or freeze-thaw. Do not open the inner packaging until ready, it must be used in one hour if opened (Humidity≤60%, Temp: 20°C -30°C). Please use immediately when the humidity>60%.

QUALITY CONTROL

Good laboratory practice suggests the use of positive and negative controls to ensure that test reagents are working and that the test is correctly performed. Labnovation SARS-CoV-2 Antigen & Influenza A/B Rapid Test Kit (Immunochromatography) contains three Positive Control Swabs and a Negative Control Swab.

SARS-CoV-2 Positive Control Swab: Non-infectious recombinant SARS-CoV-2 nucleocapsid antigen dried onto a swab

Positive Control A+/B-: Non-infectious recombinant Influenza A nucleoprotein antigens dried onto a swab

Positive Control A-/B+: Non-infectious recombinant Influenza B nucleoprotein antigens dried onto a swab.

Negative Control Swab: Sterile swab

1. Procedural Controls

A. The distinct colored lines at the "Control" position is an internal procedural control. If the test flows and the reagents work, this line will always appear.

B. The clearing of background color from the result window is a negative background control. The background color in the window should be light pink to white within 15minutes. Background color should not hinder reading of the test.

2. External Positive and Negative Controls

Positive and Negative control swabs are supplied with each kit. Prepare

kit control swabs and test using the same procedure as used for patient specimens. These swabs will monitor the entire assay. Test these swabs once with each new shipment received and once for each untrained operator. Further controls may be tested in order to conform with local, state and/or federal regulations, accrediting groups, or your lab's standard Quality Control procedures.

If the correct control results are not obtained, do not perform patient tests or report patient results. Contact Technical Support during normal business hours before testing patient specimens.

Sample Preservation

Sample of human nasopharyngeal swabs and oropharyngeal swabs should be processed as soon as possible after sample collection. If the test cannot be performed immediately, the sample should be stored in a sealed state, stored at 2-8°C for 8 hours, and stored below -20°C for 1 month. Long-term storage is not recommended.


SAMPLE COLLECTION

By nasopharyngeal swab: The sampler gently holds the head of the person to be collected with one hand, holds the swab with the other hand, sticks the swab to the nostril to enter, and slowly penetrates backwards along the bottom of the lower nasal passage, so as not to exert too much force to avoid traumatic hemorrhage. When the tip of the swab reaches the posterior wall of the nasopharynx cavity, gently rotate it for one week (in case of reflex cough, stop for a minute), and then slowly take out the swab.

By oropharyngeal swab: The head of the person to be collected is slightly tilted and his mouth is wide open, exposing the pharyngeal tonsils on both sides. Wipe the swab across the root of the tongue. Wipe the pharyngeal tonsils on both sides of the person to be collected back and forth with a little force for at least 3 times, and then wipe up and down the posterior pharyngeal wall for at least 3 times.

SAMPLE TREATMENT


Add 550µl sample extraction(about 22-24 drops) into the sample tube, dip the swab after collecting the sample into the sample extraction liquid, fully immerse the tip of the swab, rotate and squeeze the swab 10 times, then pull out the swab, and take the stranded liquid as the sample to test cassette.

TEST STEPS

Instructions must be read entirely before taking the test. Place the test and conduct the testing at room temperature. Do not open the inner packing until the sample preparation is ready. Use it as soon as possible after opening the inner packing.

1. Take out the test cassette from the sealed pouch, place it on a clean and level surface with the sample port well up.
2. Apply 2 full drops of the treated sample (60µl-70µl) vertically into each of the two sample wells of the test cassette.
3. Observe the test results immediately within 15-20 minutes, the result is invalid over 20 minutes.


Note: When using viral transport medium (VTM), it is important to ensure that the VTM containing the sample is warmed to room temperature. Cold samples will not flow correctly and can lead to erroneous or invalid results. Several minutes will be required to bring a cold sample to room temperature.

INTERPRETATION OF RESULTS

POSITIVE: The presence of T (T on nCoV /A or B on Flu) and C lines within the reaction window indicate a positive result on SARS-CoV-2 or Flu A and/or B infection or co-infection.

NEGATIVE: One colored line appears in the control region(C). No apparent colored line appears in the test region (T on nCoV /A or B on Flu). The negative result does not indicate the absence of analytes in the sample, it only indicates the level of tested analytes in the sample is less than the minimum detection limit.

INVALID: No colored lines appear, or control line fails to appear, indicating that the operator error or reagent failure. Verify the test procedure and repeat the test with a new testing device.

LIMITATIONS

1. This reagent is a qualitative detection reagent, which cannot determine the exact content of antigen.
2. The test results of this reagent are only for the reference of clinicians and should not be taken as the sole basis for clinical diagnosis and treatment. Clinical management of patients should be considered in the light of their symptoms/signs, medical history, other laboratory tests and treatment responses.
3. Restricted by antigen detection reagent method, the lowest detection limit (sensitivity analysis) is generally lower than that of nucleic acid detection, so the researchers deal with negative result to give more attention, should be combined with other test results comprehensive judgment, advice to doubt the negative result of nucleic acid detection or virus isolation culture identification method for review.
4. False negative results may be caused by unreasonable sample collection, transport and treatment, and low viral load in samples.

PERFORMANCE CHARACTERISTICS

1. Clinical Sensitivity and Specificity

SARS-CoV-2 test strip performance against PCR

| SARS-CoV-2 Antigen | PCR | | Total |
|--------------------|------------|------------|------------|
| | Positive | Negative | |
| Positive | 100 | 0 | 100 |
| Negative | 2 | 400 | 402 |
| TOTAL | 102 | 400 | 502 |

*95% CI (Confidence Interval)

Clinical Sensitivity: 98.04%, (93.10%-99.76%)

Clinical Specificity: 100%, (99.08%-100.00%)

Total consistent: 99.60%, (98.57% to 99.95%)

Influenza A test strip performance against PCR

| Influenza A | PCR | | Total |
|-------------|----------|----------|-------|
| | Positive | Negative | |
| Positive | 100 | 3 | 103 |
| Negative | 4 | 426 | 430 |

| TOTAL | 104 | 429 | 533 |
|-------|-----|-----|-----|
|-------|-----|-----|-----|

*95% CI (Confidence Interval)

Clinical Sensitivity: 96.15%(90.44%~98.94%)

Clinical Specificity: 99.30% (97.97%~99.86%)

Total consistent: 98.69% (97.31%~99.47%)

Influenza B test strip performance against PCR

| Influenza B | PCR | | Total |
|--------------|------------|------------|------------|
| | Positive | Negative | |
| Positive | 108 | 3 | 111 |
| Negative | 3 | 419 | 422 |
| TOTAL | 111 | 422 | 533 |

*95% CI (Confidence Interval)

Clinical Sensitivity: 97.30% (92.30%~99.44%)

Clinical Specificity: 99.29% (97.94%~99.85%)

Total consistent: 98.87% (97.57%~99.59%)

Limit of Detection

| Name | Subtype | LOD concentration TCID ₅₀ /mL |
|-------------|----------|--|
| Influenza A | H1N1 | 3.54 x 10 ² |
| | H3N2 | 2.28 x 10 ² |
| Influenza B | Victoria | 6.25 x 10 ² |
| | Yamagata | 1.56 x 10 ³ |
| SARS-CoV-2 | / | 1.6 x 10 ² |

3. Cross-reactivity

Cross reactivity of the SARS-CoV-2 antigen & Influenza A/B Rapid Test Kit was evaluated by testing normal respiratory track pathogenic microorganisms (bacteria, viruses, yeast) and a pooled human nasal wash that may be present in the nasal cavity. Each of the bacteria, viruses, and yeast were tested in duplicate in the absence or presence of heat inactivated SARS-CoV-2 virus (2.4 x 10² TCID₅₀/mL) or Influenza A/B. No cross-reactivity or interference came out of these tests.

Cross reactivity of the SARS-CoV-2 test strip:

With human coronavirus 229E, human coronavirus OC43, human coronavirus NL63, adenovirus, human metapneumovirus, parainfluenza virus 1, parainfluenza virus 2, parainfluenza virus 3, parainfluenza virus 4, influenza A, type B Influenza,

enterovirus, respiratory syncytial virus, rhinovirus, SARS coronavirus, MERS coronavirus, Haemophilus influenzae, Streptococcus pneumoniae, Streptococcus pyogenes, Candida albicans, Bordetella pertussis, Mycoplasma pneumonia, pneumonia Chlamydia, Legionella pneumophila, etc. have no cross reaction.

Cross reactivity of the Influenza A/B test strip:

With respiratory Syncytial Virus A, measles Virus, Mumps Virus, adenovirus Type 3, parainfluenza Type 2, partial Pulmonary Virus, human coronavirus OC43, human coronavirus 229E, human herpes virus 6, cytomegalovirus, epstein Barr Virus, enterovirus CA16, rhinovirus, staphylococcus aureus, streptococcus pneumoniae, mycoplasma pneumoniae, parapertussis, chlamydia pneumoniae, escherichia coli, pooled human nasal wash, candida albicans, etc. have no cross reaction.

4. Interfering

Common interfering substances in the sample, such as blood, mucin, and pus, have no effect on the test results.

PRECAUTIONS

1. For IN VITRO diagnostic use only.
2. Reagents should be used as soon as possible after opened. This reagent cannot be reused for disposable.
3. The test device should remain in the sealed pouches until use. If sealing problem happens, do not test. Don't use after the expiration date.
4. All specimens and reagents should be considered potentially hazardous and handled in the same manner as an infectious agent after use.

MANUFACTURER

LABNOVATION TECHNOLOGIES, INC.
 Add.: 101 and 5th Floor, Building 1, No. 68, 18th Road, Guangming Hi-Tech Park, Tangjia Community, Fenghuang Street, Guangming District, Shenzhen 518107, Guangdong, China
 Tel: 0086-755-86368398 Fax: 0086-755-86368318
 Web: www.labnovation.com
 E-mail: export@labnovation.com

Shandong Weigao Group Medical Polymer Co.,Ltd.
 No.20 Xingshan Road, Torch Hi-tech Science Park 264210 Weihai,
 Shandong Province PEOPLE'S REPUBLIC OF CHINA

INSTRUCTIONS OF SYMBOL

| | | | |
|--|-----------------------------|--|------------------------------------|
| | Consult instruction for use | | Keep dry |
| | Temperature | | Batch number |
| | For single use | | In vitro diagnostic medical device |
| | Manufacturer | | Date of manufacture |
| | Expire date | | Contains sufficient for <n> tests |
| | European representative | | CE Marked |

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PRODUCENT

LABNOVATION TECHNOLOGIES, INC.

Add.: 101 and 5th Floor, Building 1, No. 68, 18th Road, Guangming Hi-Tech Park, Tangjia Community, Fenghuang Street, Guangming District, Shenzhen

518107, Guangdong, Chiny

Tel: 0086-755-86368398 Fax: 0086-755-86368318

Web: www.labnovation.com

E-mail: export@labnovation.com

Shangdong Weigao Group Medical Polymer Co.,Ltd.

No.20 Xingshan Road, Torch Hi-tech Science Park 264210 Weihai,

Shangdong Province PEOPLE'S REPUBLIC OF CHINA

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FOR CLINICAL LABORATORY PROFESSIONAL USE ONLY

PRODUCT NAME

COVID-19 (SARS-CoV-2) IgM/IgG Antibody Test Kit
(Immunochromatography/Colloidal Gold)

INTENDED USE

The kit is a rapid test intended for the qualitative detection of the IgM and IgG antibodies to SARS-CoV-2 in human serum, plasma or whole blood sample from patients suspected of COVID-19 infection. It is an aid in the diagnosis of the patients with suspected SARS-CoV-2 infection in conjunction with clinical presentation and results of other laboratory tests. Results from this test kit should not be used as the sole basis for diagnosis. IgM antibodies to SARS-CoV-2 are generally detectable in blood several days after initial infection, Positive results for both IgG and IgM could occur after infection and can be indicative of acute or recent infection. This test is intended for professional and laboratory use.

The test provides preliminary test results. Negative results don't preclude SARS-CoV-2 infection and should not be used as the sole basis for treatment or other management decision. IgM antibodies may not be detected in the first few days of infection; the sensitivity of the qSARS-CoV-2 IgG/IgM Rapid Test early after infection is unknown.

TEST PRINCIPLE

This kit is an immunochromatography assay, using capture method for rapid, qualitative detection of severe acute respiratory syndrome coronavirus 2 (SARS-CoV-2) IgM/IgG antibody in human serum, plasma or whole blood sample. When the sample contains the SARS-CoV-2 IgM antibody, it forms a complex with the gold label antigen (SARS-CoV-2 recombinant antigen). The complex moves forward under the action of chromatography and combines with the coated antibody (Mouse anti-human IgM monoclonal antibody) at the IgM line to form a complex and develop color, which is a positive result. When the sample does not contain the SARS-CoV-2 IgM antibody, no complex can be formed and no red band appears in IgM line, which is a negative result.

When the sample contains the SARS-CoV-2 IgG antibody, it forms a complex with the gold label antigen (SARS-CoV-2 recombinant antigen) . The complex moves forward under the action of chromatography and combines with the coated antibody(Mouse anti-human IgG monoclonal antibody). at the IgG line to form a complex and develop color, which is a positive result. When the sample does not contain the SARS-CoV-2 IgG antibody, no complex can be formed and no red band appears in IgG line, which is a negative result. Regardless of whether the

SARS-CoV-2 IgG or IgM antibody is contained in the sample, the gold label quality control antibody will bind with the coated antibody at the C line to form a complex and develop color (C line).

MAIN COMPONENTS

The test cassette.

Sample dilution: composed of 20 mM phosphate buffer solution (PBS).

STORAGE AND EXPIRY

Store as packaged in the sealed pouch at 2-30°C, avoid hot and sunshine, dry place, valid for 12 months. DO NOT FREEZE. Some protective measures should be taken in hot summer and cold winter to avoid high temperature or freeze-thaw. Do not open the inner packaging until ready, it must be used in one hour if opened (Humidity ≤60%, Temp: 20°C-30°C). Please use immediately when the humidity >60%.

SAMPLE REQUIREMENT

1. The reagent can be used for the serum, plasma and whole blood samples.
2. A serum / plasma / whole blood sample must be collected in a clean and dry container. EDTA, sodium citrate, heparin can be used as anticoagulants in plasma / whole blood samples. Detect immediately after collecting blood.
3. Serum and plasma samples may be stored at 2-8°C for 3 days prior to assay. If testing is delayed more than 3 days, the sample should be frozen (-20°C or colder). Repeat freeze and thaw for no more than 3 times. Whole blood samples with anticoagulant can be stored at 2-8°C for 3 days, and should not be frozen; whole blood samples without anticoagulant should be used immediately (if the sample has agglutination, it can be detected by serum) .

TEST METHODS

Instructions must be read entirely before taking the test. Allow the test device controls to equilibrate to room temperature for 30 minutes (20°C -30°C) prior to testing. Do not open the inner packaging until ready, it must be used in one hour if opened (Humidity ≤60%, Temp: 20°C -30°C). Please use immediately when the humidity >60%.

For Serum/Plasma

1. Remove the test device from the sealed pouch, place it on a clean and level surface with the sample well up.
2. Add serum or plasma (10µl) vertically into the sample well .
3. Add two (2) drops (80-100µl) of sample buffer into the sample well .
4. Observe the test results immediately within 15~20 minutes, the

result is invalid over 20 minutes.

For Whole Blood

1. Remove the test cassette from the sealed pouch, place it on a clean and level surface with the sample well up.
2. Add whole blood (20µl) vertically into the sample well .
3. Add two (2) drops (80-100µl) of sample buffer into the sample well .
4. Observe the test results immediately within 15~20 minutes, the result is invalid over 20 minutes.


INTERPRETATION OF RESULTS

POSITIVE: Two or three distinct red lines appear. One line should be in the control region (C) and the other line would be in the test region (IgG line, IgM line or both).

NEGATIVE: One red line appears in the control region(C). No red or pink line appears in the test region (IgG line, IgM line).

INVALID: No red lines appear or control line fails to appear, indicating that the operator error or reagent failure. Verify the test procedure and repeat the test with a new testing device.

LIMITATIONS

1. This reagent is designed for the qualitative screening test. Concentration of SARS-CoV-2 IgM/IgG antibody cannot be determined by this qualitative test. The depth of the T-line(IgG/IgM) color is not necessarily related to the concentration of the antibody in the sample.

2. The results of the reagent are only for clinical reference, which is not the only basis for clinical diagnosis and treatment. A confirmed diagnosis and treatment should only be made by a physician after all clinical and laboratory findings have been evaluated.

PERFORMANCE CHARACTERISTICS

1. Sensitivity and Specificity

| SARS-CoV-2 IgM Ab Rapid Test | PCR Test | | Total |
|---------------------------------|----------|----------|-------|
| | Positive | Negative | |
| Positive | 246 | 40 | 286 |
| Negative | 54 | 960 | 1014 |
| TOTAL | 300 | 1000 | 1300 |

Analysis of coincidence rate of SARS-CoV-2 IgM Ab rapid test and PCR Test in serum samples:

Positive coincidence rate= $246 / (246+54) \times 100\% = 82\%$,

Negative coincidence rate= $960 / (40+960) \times 100\% = 96\%$,

Total coincidence rate= $(246+960) / (246+54+40+960) \times 100\% = 92.8\%$.

| SARS-CoV-2 IgG Ab Rapid Test | PCR Test | | Total |
|---------------------------------|----------|----------|-------|
| | Positive | Negative | |
| Positive | 279 | 25 | 304 |
| Negative | 21 | 975 | 996 |
| TOTAL | 300 | 1000 | 1300 |

Analysis of coincidence rate of SARS-CoV-2 IgG Ab rapid test and PCR Test in serum samples:

Positive coincidence rate= $279 / (279+21) \times 100\% = 93\%$,

Negative coincidence rate= $975 / (25+975) \times 100\% = 97.5\%$,

Total coincidence rate= $(279+975) / (279+21+25+975) \times 100\% = 96.5\%$.

2. Cross-reactivity

Specimens which contain high titered antibodies to the pathogens listed below were investigated with SARS-CoV-2 IgM/IgG Ab Rapid Test. The results showed no

cross reactivity.

| | |
|-----------------------------|------------------------------|
| Human coronavirus 229E | Rhinovirus |
| Human coronavirus OC43 | Streptococcus pneumoniae |
| Human coronavirus NL63 | Streptococcus pyogenes |
| Adenovirus type 3 | Epstein-Barr virus |
| Human Metapneumovirus | Human Immunodeficiency virus |
| Parainfluenza virus 1 | Plasmodium falciparum |
| Parainfluenza virus 2 | Plasmodium ovale |
| Parainfluenza virus 3 | Dengue virus |
| Parainfluenza virus 4 | Bordetella pertussis |
| Influenza A | Mycoplasma pneumoniae |
| Influenza B | Chlamydia pneumoniae |
| Human Metapneumovirus | Legionella pneumophila |
| Respiratory syncytial virus | |

3. Interferences

The test results of SARS-CoV-2 IgM/IgG Ab Rapid Test do not be interfered with the substance at the following concentration:

| | |
|-----------------------------|------------|
| Hemoglobin | 10 mg/ml |
| Bilirubin Conjugated | 0.4 mg/ml |
| Bilirubin Unconjugated | 0.4 mg/ml |
| Human serum protein | 60 mg/ml |
| Triglycerides | 15 mg/ml |
| Cholesterol | 4 mg/ml |
| Human anti-mouse antibody | 100 IU/ml |
| Rheumatoid Factor | 1500 IU/ml |
| ANA anti-nuclear antibodies | 100 IU/ml |

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BIBLIOGRAPHY

- [1] Nanshan Chen, Min Zhou, Xuan Dong, et al. (2020). Epidemiological and clinical characteristics of 99 cases of 2019 novel coronavirus pneumonia in Wuhan, China: a descriptive study. The Lancet, 395(10223), 507-513.
- [2] Chaolin Huang, Yeming, et al. (2020). Clinical features of patients infected with 2019 novel coronavirus in Wuhan, China. The Lancet, 395(10223), 497-506.
- [3] Chen Wang, Peter W Horby, Frederick G Hayden, George F Gao. (2020). A novel coronavirus outbreak of global health concern. The Lancet, 395(10223), 470-473.
- [4] Na Zhu, Ph.D., Dingyu Zhang, M.D., Wenling Wang, Ph.D., et al. (2020). A Novel Coronavirus from Patients with Pneumonia in China, 2019. The New England Journal of Medicine.
- [5] World Health Organization: Clinical management of severe acute respiratory infection when Novel coronavirus (nCoV) infection is suspected: Interim Guidance. 12 January, 2020.

MANUFACTURER

LABNOVATION TECHNOLOGIES, INC.

Add.: 101 and 5th Floor, Building 1, No. 68, 18th Road, Guangming Hi-Tech Park, Tangjia Community, Fenghuang Street, Guangming District, Shenzhen 518107, Guangdong, China

Tel: 0086-755-86368398 Fax: 0086-755-86368318

Web: www.labnovation.com



E-mail: export@labnovation.com

EUROPEAN REPRESENTATIVE

MedNet EC-REP GmbH

Borkstrasse 10, 48163 Muenster, Germany

INSTRUCTIONS OF SYMBOL

| | | | |
|---|-----------------------------|---|------------------------------------|
|  | Consult instruction for use |  | Keep dry |
|  | Temperature | LOT | Batch number |
|  | For single use | IVD | In vitro diagnostic medical device |
|  | Manufacturer |  | Date of manufacture |
|  | Expire date |  | Contains sufficient for <n> tests |
| EC REP | European representative | CE | CE Marked |

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BIBLIOGRAFIA

(1) Nanshan Chen, Min Zhou, Xuan Dong, et al. (2020). Epidemiological and clinical characteristics of 99 cases of 2019 novel coronavirus pneumonia in Wuhan, China: a descriptive study. The Lancet, 395(10223), 507-513.

(2) Chaolin Huang, Yeming, et al. (2020). Clinical features of patients infected with 2019 novel coronavirus in Wuhan, China. The Lancet, 395(10223), 497-506.

(3) Chen Wang, Peter W Horby, Frederick G Hayden, George F Gao. (2020). A novel coronavirus outbreak of global health concern. The Lancet, 395(10223), 470-473.

(4) Na Zhu, Ph.D., Dingyu Zhang, M.D., Wenling Wang, Ph.D., et al. (2020). A Novel Coronavirus from Patients with Pneumonia in China, 2019. The New England Journal of Medicine.

(5) World Health Organization: Clinical management of severe acute respiratory infection when Novel coronavirus (nCoV) infection is suspected: Interim Guidance. 12 Styczeń, 2020.

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Park, Tangjia Community, Fenghuang Street, Guangming District, Shenzhen 518107, Guangdong, Chiny

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