Product Catalogue Number FP-014-2 For clinical diagnosis and scientific research.

CLL Chromosome and Gene Anomaly Probe Detection Kit

[Product Name] CLL Chromosome and Gene Anomaly Probe Detection Kit (Fluorescence In Situ Hybridization Method)

[Intended use]

This product is mainly used to detect chromosome and gene abnormalities in chronic lymphoblastic leukemia; The applicable population is patients suspected or diagnosed with chronic lymphocytic leukemia by routine clinical examination.

Chronic lymphocytic leukemia (CLL) is a mature B-lymphocyte clonal proliferative tumor characterized by lymphocyte aggregation in peripheral blood, bone marrow, spleen and lymph nodes. Peripheral blood B lymphocytes with clinical diagnosis of persistence (3 months) $\geq 5 \times 10^9/L$ (e.g. peripheral blood B cells $< 5 \times 10^9/L$), Chronic lymphocytic leukemia is also diagnosed in patients with hematopenia caused by bone marrow infiltration or disease-related symptoms. Fluorescence in situ hybridization showed that about 80% of patients with chronic lymphoblastic leukemia had chromosome abnormalities, and the most common deletion was in the long arm del of chromosome 13 (13q14.1); Chromosome 12 deletion or trisomy; Chromosome 17 short arm deletion del (17p). These abnormalities are of great significance for the diagnosis, differential diagnosis, treatment and prognosis of chronic lymphocytic leukemia.

This kit is not clinically verified in combination with gene targeted therapeutic drugs, but only for gene detection performance. This kit is only applicable to the detection of chronic lymphocytic leukemia and provides doctors with auxiliary information for diagnosis.

[Detection principle]

Based on the fluorescence in situ hybridization technology, a nucleotide of the nucleic acid probe is labeled with fluorescein. The detected target gene is homologous and complementary with the used nucleic acid probe. After denaturation, annealing and renaturation, they can form a hybrid between the target gene and the nucleic acid probe. The detected gene is analyzed qualitatively, quantitatively or relatively under the microscope by the fluorescence detection system. The kit adopts orange probe labeled with orange fluorescein and green probe labeled with green fluorescein. The two probes can be combined with the target detection site by in situ hybridization.

Under normal conditions (without gene deletion and chromosome abnormality), it is displayed as two orange signals and two green signals under fluorescence microscope. When there is gene deletion, there will be a lack of green or orange signal. When there is chromosome polysomy, the probe signal will increase. Gene deletion and chromosome abnormalities were detected by this method, so as to provide reference basis for clinical identification, prognosis judgment and medication of leukemia patients.

[Product Introduction]


The kit uses orange fluorescein-labeled RB1 probe and green fluorescein-labeled ATM probe to bind the RB1 and ATM probe to the target detection site by in situ hybridization.

[Product Composition]

The kit consists of P53/CEP17 dual-color probes, as shown in Table 1.

Table 1 Kit composition

Package	Component name	Specifications	Quantity	Main components
10 Tests/box	P53/CEP17 dual-color probe	100 μ L/Tube	1	P53 Orange probe, CEP17 Green

HealthCare
Fast FISH Technology Leader

Wuhan HealthCare Biotechnology Co., Ltd.

Building #8, Optics Valley Precision Medicine Industry Base,
#9 Gaokeyuan 3rd Road, East Lake High-Tech Zone. 430206, Wuhan, Hubei, China.

Tel#: +86(027)8757-0662 Ext. #814 E-mail: cs@healthcare-bio.com | Website: www.healthcare-biotech.com

[Storage conditions & Validity]

Keep sealed away from light at -20°C±5°C. The product is valid for 20 months. Avoid unnecessary repeated freezing and thawing that should not exceed 10 times. After opening, within 24 hours for short-term preservation, keep sealed at 2-8°C in dark. For long-term preservation after opening, keep the lid sealed at -20°C±5°C away from light. The kit is transported below 0°C.

[Applicable Instruments]

1. Fluorescence microscopy imaging system includes fluorescence microscopy and filter sets for DAPI, Green, and Orange.

[Sample Requirements]

1. Applicable specimen types: Fresh bone marrow specimens stored at 4°C for no more than 24 hours and used for FISH testing after culture.
2. If the sample is stored at a too high or too low (such as freezing) temperature, the sample will not be used for testing and should be discarded.
3. Bone marrow cell suspension for karyotyping should be stored at -20°C for FISH detection.
4. If the cells suspension is too volatile or contaminated during storage, the sample should be discarded.

[Test method]

1.Related Reagents

The following reagents are required for the experiment but not provided in this kit

① 20×SSC, pH 5.3±0.2

Weigh 176g of sodium chloride and 88g of sodium citrate, dissolve in 800mL of deionized water, adjust the pH to 5.3±0.2 at room temperature, and complete to 1 L with deionized water. High-pressure steam sterilization, stored at 2-8°C, the solution shelf life is of 6 months. Discard if the reagent appears cloudy (turbid) or contaminated.

② 2×SSC, pH 7.0±0.2

Take 100mL of the above 20×SSC, dilute with 800mL deionized water, mix, adjust the pH to 7.0±0.2 at room temperature, complete to 1L with deionized water, stored at 2-8°C, the shelf life is of 6 months. Discard if the reagent appears cloudy (turbid) or contaminated.

③ 0.3% NP-40/0.4×SSC solution, pH 7.0-7.5

Take 0.6mL NP-40 and 4mL 20×SSC, add 150mL deionized water, mix, adjust the pH to 7.0-7.5 at room temperature, with deionized water complete to a volume of 200mL. Stored at 2-8°C, the shelf life is of 6 months. Discard if the reagent appears cloudy (turbid) or contaminated.

④ Fixation solution (methanol: glacial acetic acid = 3:1)

Prepare a ready to use fixation solution by mixing thoroughly 30ml of methanol and 10ml of glacial acetic acid.

⑤ PBS buffer, pH 7.4 ± 0.2

Prepare a ready to use fixation solution by mixing thoroughly 30ml of methanol and 10ml of glacial acetic acid.

Sodium chloride	8g
Potassium chloride	0.2g
Sodium hydrogen phosphate	3.58g
Potassium dihydrogen phosphate	0.27g

Dissolve the above reagents in 800mL deionized water, mix, adjust the pH to 7.4 ± 0.2 at room temperature, with deionized water complete to a volume of 1000mL. Stored at room temperature for 6 months. Discard if the reagent appears cloudy (turbid) or contaminated.

⑥ Ethanol Solution: 70% ethanol, 85% ethanol

Dilute 700ml, 850ml of ethanol with deionized water to 1L. The shelf life is of 6 months. Discard if the reagent appears cloudy (turbid) or contaminated.

Wuhan HealthCare Biotechnology Co., Ltd. | Fast[®]robe[®]

HealthCare
Fast FISH Technology Leader

Your Authentic Partner In Molecular Cytogenetics

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Wuhan HealthCare Biotechnology Co., Ltd.

Building #8, Optics Valley Precision Medicine Industry Base,

#9 Gaokeyuan 3rd Road, East Lake High-Tech Zone. 430206, Wuhan, Hubei, China. **Tel#:** +86(027)8757-0662

Wew. #812, **E-mail:** cs@healthcare-bio.com | **Website:** www.healthcare-biotech.com

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TP53/ATM FISH Probe Kit

Introduction

The TP53/ATM FISH Probe Kit is designed to detect rearrangements involving the human *TP53* and *ATM* genes, located on chromosome bands 17p13.1 and 11q22.3, respectively. Abnormalities in both genes, the *TP53* gene – also known as *P53*, *BCC7*, *LFS1* or *TRP53* – and the *ATM* gene – also called *AT1*, *ATA*, *ATC*, *ATD*, *ATE*, *ATDC*, *TEL1* or *TELO1*, have been reported in a number of leukemia types and other malignancies.

Intended Use

To measure the copy number of the human *TP53* and *ATM* genes located on chromosome bands 17p13.1 and 11q22.3, respectively.


Cont.

LSP TP53 FISH Probe
LSP ATM FISH Probe

Color

CytoOrange
CytoGreen

Probe Design

LSP TP53 FISH Probe covers a chromosomal region which includes the entire *TP53* gene. LSP ATM FISH Probe covers a chromosomal region including the entire *ATM* gene.

Cat. No.

CT-PAC106-10-OG

Volume

10 Tests (100 µL)

Signal Pattern Interpretation

Normal Pattern

2O + 2G

Abnormal Pattern

Other Patterns

- 1) Varley JM, et al. *Br J Cancer*. 76(1):1-14 (1997).
- 2) Vogelstein B, et al. *Nature*. 408(6810):307-10 (2000).
- 3) Soussi T & Bérout C. *Nat Rev Cancer*. 1(3):233-40 (2001).
- 4) Olivier M, et al. *Cancer Res*. 63(20):6643-50 (2003).
- 5) Petitjean A, et al. *Oncogene*. 26(15):2157-65 (2007).

CytoTest Inc.
9430 Key West Ave., Suite 210
Rockville, MD 20850, USA

* CE IVD only available in certain countries. All other countries are either ASR or RUO. Please contact your local dealer or our headquarters for more information.

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www.cytotest.com

numer wersji

Product Catalogue Number FP-197 For clinical diagnosis and scientific research.

1q21 and 1p32 Gene Anomaly Probe Detection Kit

[Product Name] 1q21 and 1p32 Gene Anomaly Probe Detection Kit (Fluorescence In-Situ Hybridization Method).

[Product Introduction]


The kit uses orange fluorescein-labeled 1p32(CDKN2C), and green fluorescein-labeled 1q21(CKS1B), and bind 1p32/1q21gene with the target site by in situ hybridization.

[Product Composition]

The kit consists of 1q21 orange-red probe as shown in Table 1.

Table 1: Kit composition

Package Specifications	Component name	Specifications	Quantity	Main components
10 Tests/box	1p32/1q21 orange-red probe	100μL/Tube	1	1p32 Orange probe, 1q21Green probe

[Storage conditions & Validity]

Keep sealed away from light at -20°C±5°C. The product is valid for 12 months. Avoid unnecessary repeated freezing and thawing that should not exceed 10 times. After opening, within 24 hours for short-term preservation, keep sealed at 2-8°C in dark. For long-term preservation after opening, keep the lid sealed at -20°C±5°C away from light. The kit is transported below 0°C.

[Applicable Instruments]

Fluorescence microscopy imaging systems, including fluorescence microscopy and filter sets suitable for DAPI (367/452), Green (495/517), and Orange (547/565).

[Sample Requirements]

1. Applicable specimen types: Fresh bone marrow samples that have not been fixed are stored at 4°C for less than 24 hours; bone marrow cells suspension after fixation is stored at -20°C for less than 6 months; the prepared bone marrow cell slides can be set at -20°C for less than 1 month.
2. When specimens are stored at too high or too low a temperature (eg, frozen), the specimen will not be used for testing and should be discarded.
3. If the cells suspension is excessively volatile or contaminated during storage, the sample should be discarded.

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CBFB-MYH11 Fusion/Translocation FISH Probe Kit

Introduction

The CBFB-MYH11 Fusion/Translocation FISH Probe Kit is designed to detect rearrangements involving the human *CBFB* and *MYH11* genes located on chromosome bands 16q22.1 and 16p13.11, respectively. Rearrangements between the two genes, the *CBFB* gene – also known as *CBFb* or *PEBP2B* – and the *MYH11* gene – also called *AAT4*, *FAA4*, *SMHC* or *SMMHC*, have been observed in acute myeloid leukemia (AML) and other hematological malignancies.

Intended Use

To detect rearrangements involving the human *CBFB* and *MYH11* genes located on chromosome bands 16q22.1 and 16p13.11, respectively.


Cont.

LSP CBFB 5'-3' FISH Probe
LSP MYH11 5'-3' FISH Probe

Color

CytoOrange
CytoGreen

Probe Design

Not to Scale

LSP CBFB 5'-3' FISH Probe covers the 5' (start) and 3' (end) portion of the *CBFB* gene and some genomic sequences adjacent to the two ends of the gene. LSP MYH11 5'-3' FISH Probe covers about the entire *MYH11* gene as well as sequences upstream (5' start) and downstream (3' end) of the gene. The probe set is optimized to reveal translocations between the two genes.

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DEK-NUP214 Fusion/Translocation FISH Probe Kit

Introduction

The DEK-NUP214 Fusion/Translocation FISH Probe Kit is designed to detect rearrangements involving the human *DEK* and *NUP214* genes located on chromosome bands 6p22.3 and 9q34.13, respectively. Rearrangements involving portions of these two genes, the *DEK* gene – also known as *D6S231E* – and the *NUP214* gene – also called *CAN*, *CAIN*, or *D9S46E*, have been observed in acute myeloid leukemia (AML), myelodysplastic syndrome (MDS) and many other hematological malignancies.

Intended Use

To detect rearrangements involving the human *DEK* and *NUP214* genes located on chromosome bands 6p22.3 and 9q34.13, respectively.


Cont.

LSP DEK FISH Probe
LSP NUP214 FISH Probe

Color

CytoGreen
CytoOrange

Probe Design

LSP DEK FISH Probe covers the entire *DEK* gene along with some upstream (5') and downstream (3') flanking genomic sequences. LSP NUP214 FISH Probe spans the complete *NUP214* gene, as well as adjacent 5' and 3' genomic sequences. The probe set is optimized to reveal the typical t(6;9)(p22;q34) translocation as well as other translocations between the two genes.

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NUP98 Break Apart FISH Probe Kit


Introduction

The NUP98 Break Apart FISH Probe Kit is designed to detect rearrangements in the human *NUP98* gene located on chromosome band 11p15.4. In addition to revealing breaks, which can lead to translocation of parts of the gene, inversion, or its fusion to other genes, the probe set can also be used to identify other *NUP98* aberrations such as deletions or amplifications. Rearrangements and abnormal expression of the *NUP98* gene – also known as *NUP96*, *NUP196* or *ADIR2* – have been observed in acute lymphoblastic leukemia (ALL), acute myeloid leukemia (AML), myelodysplastic syndrome (MDS) and other malignancies.

Intended Use
To detect rearrangements in the human <i>NUP98</i> gene located on chromosome band 11p15.4.

Cont.	Color
LSP NUP98 5' FISH Probe	CytoOrange
LSP NUP98 3' FISH Probe	CytoGreen

Probe Design

Cat. No.	Volume
CT-PAC180-10-OG	10 Tests (100 µL)

Signal Pattern Interpretation	
<u>Normal Patterns</u>	<u>Abnormal Patterns</u>
2F*	Other Patterns

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ZBTB16-RARA Fusion/Translocation FISH Probe Kit

Introduction

The ZBTB16-RARA Fusion/Translocation FISH Probe Kit is designed to detect rearrangements involving the human *ZBTB16* and *RARA* genes located on chromosome bands 11q23.3 and 17q21.2, respectively. Rearrangements between the two genes, the *ZBTB16* gene – also known as *PLZF* or *ZNF145* – and the *RARA* gene – also called *RAR* or *NR1B1*, have been observed in acute promyelocytic leukemia (APL, the M3 subtype of acute myelogenous leukemia) and other malignancies.

Intended Use

To detect rearrangements involving the human *ZBTB16* and *RARA* genes located on chromosome bands 11q23.2 and 17q21.2, respectively.


Cont.

Color

LSP ZBTB16 FISH Probe
LSP RARA FISH Probe

CytoOrange
CytoGreen

Probe Design

Not to Scale

LSP ZBTB16 FISH Probe covers a chromosomal region which includes the entire *ZBTB16* gene. LSP RARA FISH Probe covers a chromosomal region which includes the entire *RARA* gene. The probe set is optimized to reveal translocations between the two genes.

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JAK2 Break Apart FISH Probe Kit

Introduction

The JAK2 Break Apart FISH Probe Kit is designed to detect rearrangements in the human *JAK2* gene located on chromosome band 9p24.1. In addition to revealing breaks, which can lead to translocation of parts of the gene, inversion, or its fusion to other genes, the probe set can also be used to identify other *JAK2* aberrations such as deletions or amplifications. Rearrangements and abnormal expression of the *JAK2* gene – also known as *THCYT3* or *JTK10* – have been observed in acute myeloid and lymphoid leukemias and other malignancies.

Intended Use

To detect rearrangements in the human *JAK2* gene located on chromosome band 9p24.1.


Cont.

LSP JAK2 5' FISH Probe
LSP JAK2 3' FISH Probe

Color

CytoOrange
CytoGreen

Probe Design

LSP JAK2 5' FISH Probe covers the 5' (start) portion of the *JAK2* gene and some adjacent genomic sequences. LSP JAK2 3' FISH Probe covers the 3' (end) part as well as sequences downstream of the gene. The two probes are flanking sequences across the *JAK2* gene in which variable breakpoints have been observed.

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MYC Break Apart FISH Probe Kit

Introduction

The MYC Break Apart FISH Probe Kit is designed to detect rearrangements in the human *MYC* gene located on chromosome band 8q24.21. In addition to revealing breaks, which can lead to translocation of parts of the gene, inversion, or its fusion to other genes, the probe set can also be used to identify other *MYC* aberrations such as deletions or amplifications. Rearrangements and abnormal expression of the *MYC* gene – also known as *MRTL*, *MYCC*, *c-Myc* or *bHLHe39* – have been observed in Burkitt's Lymphoma and other hematological malignancies, myeloma, as well as breast, cervical, colon, ovarian and other tumor types.

Intended Use

To detect rearrangements in the human *MYC* gene located on chromosome band 8q24.21.


Cont.

LSP MYC 5' FISH Probe
LSP MYC 3' FISH Probe

Color

CytoOrange
CytoGreen

Probe Design

LSP MYC 5' FISH Probe covers genomic sequences adjacent to the 5' (start) portion of the *MYC* gene. LSP MYC 3' FISH Probe covers sequences downstream the 3' end of the gene. The two probes are flanking sequences across the *MYC* gene in which variable breakpoints have been observed.

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IGK Break Apart FISH Probe Kit

Introduction

The IGK Break Apart FISH Probe Kit is designed to detect rearrangements in the human IGK locus mapping to chromosome band 2p11.2. In addition to revealing breaks, which can lead to translocation of parts of the gene, inversion, or its fusion to other genes, the probe set can also be used to identify other IGK aberrations such as deletions or amplifications. Rearrangements and abnormal expression of the IGK gene – also known as IGK@ – have been observed in various B-cell lymphoma subtypes and other malignancies.

Intended Use

To detect rearrangements in the human *IGK* locus situated on chromosome band 2p11.2.


Cont.

LSP IGK 5' FISH Probe
LSP IGK 3' FISH Probe

Color

CytoOrange
CytoGreen

Probe Design

Not to Scale

LSP IGK 5' FISH Probe covers the 5' (start) portion of the *IGK* locus and some adjacent genomic sequences. LSP IGK 3' FISH Probe covers sequences at the 3' (end) of the gene. The two probes are flanking sequences across the *IGK* locus in which variable breakpoints have been observed.

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BCL2 Break Apart FISH Probe Kit

Introduction

The BCL2 Break Apart FISH Probe Kit is designed to detect rearrangements in the human *BCL2* gene located on chromosome band 18q21.33. In addition to revealing breaks, which can lead to translocation of parts of the gene, inversion, or its fusion to other genes, the probe set can also be used to identify other *BCL2* aberrations such as deletions or amplifications. Rearrangements and abnormal expression of the *BCL2* gene – also known as *Bcl-2* or *PPP1R50* – is routinely found in follicular lymphoma but also occurs in many other hematological and solid cancer types. One particularly common rearrangement is a reciprocal translocation to the *IGH@* locus on chromosome 14.

Intended Use

To detect rearrangements in the human *BCL2* gene located on chromosome band 18q21.33.


Cont.

LSP BCL2 5' FISH Probe
LSP BCL2 3' FISH Probe

Color

CytoOrange
CytoGreen

Probe Design

LSP BCL2 5' FISH Probe covers sequences adjacent to the 5' (start) portion of the *BCL2* gene. LSP BCL2 3' FISH Probe covers some sequences downstream (3' end) of the gene. The two probes are flanking sequences across the *BCL2* gene in which variable breakpoints have been observed.

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BCL6 Break Apart FISH Probe Kit

Introduction

The BCL6 Break Apart FISH Probe Kit is designed to detect rearrangements in the human *BCL6* gene located on chromosome band 3q27.3. In addition to revealing breaks, which can lead to translocation of parts of the gene, inversion, or its fusion to other genes, the probe set can also be used to identify other *BCL6* aberrations such as deletions or amplifications. Rearrangements of the *BCL6* gene – also known as *BCL5*, *LAZ3*, *BCL6A*, *ZNF51* or *ZBTB27* – have been observed in B-cell lymphomas and leukemias. *BCL6* is also dysregulated in multiple myeloma cases and several solid tumor types. More than 30 different translocation partner genes have been described.

Intended Use

To detect rearrangements in the human *BCL6* gene located on chromosome band 3q27.3.


Cont.

LSP BCL6 5' FISH Probe
LSP BCL6 3' FISH Probe

Color

CytoOrange
CytoGreen

Probe Design

LSP BCL6 5' FISH Probe covers some genomic sequences adjacent to the 5' end of the *BCL6* gene. LSP BCL6 3' FISH Probe covers some sequence downstream of the 3' end of the gene. The two probes are flanking sequences across the *BCL6* gene in which variable breakpoints have been observed.

Not to Scale

Cat. No.

CT-PAC207-10-OG

Volume

10 Tests (100 µL)

Signal Pattern Interpretation

Normal Pattern

2OG*

Abnormal Pattern

Other Patterns

*Overlapping orange and green signals can appear as yellow.

- 1) Chao DT & Korsmeyer SJ. *Annu Rev Immunol.* 16:395-419 (1998).
- 2) Shaffer AL, et al. *Immunity.* 13(2):199-212 (2000).
- 3) Migliazza A, et al. *Proc Natl Acad Sci U S A.* 92(26):12520-4 (1995).
- 4) Horn H, et al. *Blood.* 121(12):2253-63 (2013).
- 5) Duy C, et al. *Nature.* 473(7347):384-8 (2011).

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KAT6A-CREBBP Dual Fusion/Translocation FISH Probe Kit

Introduction

The KAT6A-CREBBP Dual Fusion/Translocation FISH Probe Kit is designed to detect rearrangements involving the human KAT6A and CREBBP genes, located on chromosome bands 8p11.21 and 16p13.3, respectively. KAT6A is also known as MOZ, MRD32, MYST3, MYST-3, ZNF220, RUNXBP2 or ZC2HC6A. CREBBP is also known as CBP, RSTS, KAT3A or RSTS1. Rearrangements involving portions of these two genes have been observed in various non-Ewing sarcomas, acute non-lymphocytic (AML) and monocytic leukemias, myelodysplastic syndrome and other tumor types.

Intended Use

To detect rearrangements involving the human KAT6A and CREBBP genes, located on chromosome bands 8p11.21 and 16p13.3, respectively.


Cont.

LSP KAT6A 5'-3' FISH Probe
LSP CREBBP 5'-3' FISH Probe

Color

CytoGreen
CytoOrange

Probe Design

LSP KAT6A 5'-3' FISH Probe and LSP CREBBP 5'-3' FISH Probe cover the 5' (start) and 3' (end) portion and some adjacent genomic sequences of the KAT6A and CREBBP genes, respectively. The two probes are flanking sequences across the KAT6A and CREBBP genes in which variable breakpoints have been observed.

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Trypsin (0.25 %) in DPBS (1x)
 Cat. No. TRY-2B (100 ml)

General Information

Trypsin solutions are used to detach adherent cells from culture surfaces. They are composed of natural porcine pancreas-derived trypsin. The concentration of trypsin necessary to dislodge cells from their substrate is dependent primarily on the cell type and the age of the culture. Various formulations should be tested to determine the best product for a specific application.

Product Specifications

Appearance	Clear frozen liquid
Storage and shelf life	Store at $\leq -15^{\circ}\text{C}$. Avoid repeated freeze-thaw cycles. Preparation of aliquots recommended. Once opened, store at 4°C and use within 2-4 weeks.
Shipping conditions	Frozen (Dry ice)
Thawing	$+37^{\circ}\text{C}$ water bath or overnight at $+2^{\circ}\text{C}$ to $+8^{\circ}\text{C}$. Swirl gently to homogenize.

Formulation

Components	Concentration mg/L
KCl	200.00
KH_2PO_4	200.00
NaCl	8000.00
Na_2HPO_4	1150.00
Trypsin	2500.00

Instructions for Use

Detachment of adherent cells using Trypsin:

Trypsin (0.25 %) in DPBS (1x) solution is supplied as a sterile, ready-to-use, frozen liquid. This entire procedure should be done in a laminar flow hood using proper aseptic technique.

1. The product can either be thawed in a $+37^{\circ}\text{C}$ water bath or overnight at $+2^{\circ}\text{C}$ to $+8^{\circ}\text{C}$.
2. Carefully aspirate all of the media from the cell culture flask.
3. Rinse cells with Ca^{2+} and Mg^{2+} -free salt solution (see related products), aspirate, and discard.
4. Prewarm the trypsin solution in a $+37^{\circ}\text{C}$ water bath. Add enough trypsin solution to completely cover the cells.
5. Incubate the flask at $+37^{\circ}\text{C}$, or for more sensitive cultures, at room temperature or $+2^{\circ}\text{C}$ to $+8^{\circ}\text{C}$.
6. When the trypsinization process is complete, cells will appear rounded upon microscopic examination and the solution in the flask will appear cloudy. Check the flask often to avoid overexposure. Trypsin can cause cellular damage and time of exposure should be kept to a minimum.
The time required to detach cells from the culture surface is dependent on the cell type, the age of the culture, population density, serum concentration in the growth medium and time since last subculture.
7. Neutralize trypsin either with serum containing medium or trypsin inhibitor. Gently centrifuge the cell suspension and discard the trypsin-containing supernatant.
8. Resuspend the cell pellet with fresh medium and count or culture as desired.

Product Information

Trypsin (0.25 %) in DPBS (1x)
Cat. No. TRY-2B (100 ml)

Related Products

Product	Cat. No.
Dulbecco's PBS (1x), w/o Ca & Mg, w/o Phenol Red	PBS-1A
Hank's Balanced Salts, w/o Ca & Mg, w/o Phenol Red	HBSS-2A

Precautions and Disclaimer

The product is for research and further manufacturing only.

Help Needed?

If you have any further questions regarding this product, please do not hesitate to contact our cell culture experts by email (techservice@capricorn-scientific.com) or phone (+49 6424 944640).

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©Cytotest Inc