

CareFusion Surgical Clippers general purpose blade

Provides a versatile choice for preoperative hair removal

The general purpose blade provides an effective solution for most procedural hair removal. As part of the CareFusion Surgical Clippers system, the general purpose blade quickly removes all types of body hair with a cut that's similar to a razor, while leaving the skin intact. You and your patients will feel the difference.



The general purpose blade offers the following features:

- Wide blade for efficient cutting
- Clipping closer to the skin without disrupting it
 - .23 mm with the 4406 general purpose blade compared to .6 mm with the 3M 9660 blade
- Existing clipper handle
- Disposable blade
- Wet/Dry capability
- Drop-in charging so clippers are always ready





1. Attach the blade

With gloves on position a new blade on the top of the clipper head. Align the mark on the clipper blade to the arrow. Slide the blade toward the dot to secure it in place until you hear a click.



2. Remove the blade

Turn the clipper off, position the clipper blade facing downward over a sharps trash receptacle and use your thumb to push the blade forward along the blade frame.

Tips for clipping with the general purpose blade:

- Hold the clippers like a pen or pencil with the CareFusion logo facing up
- Rest the base of the blade flat on the surface of the skin
- Clip against the grain of the hair
- Clip only where necessary



General purpose blade

Model no.	4406
Blade width	37.2 mm
Cutting ability (no. strands/second)	12,000
Blade cut height	.23 mm
Facility uses	For all types of body hair

Ordering information

Cat. no.	Description	Qty. cs.
4406	General purpose blade	50
4412A	Neuro blade (for head hair and other very thick hair)	20
4403A	SensiClip blade (for sensitive areas)	20
5513E	Rechargeable surgical clipper	1
5514A	Charging adapter	1

For more information, including a video highlighting the features of the CareFusion Surgical Clippers, visit carefusion.com/surgicalclippers.

CareFusion
Vernon Hills, IL

carefusion.com

