[www.thermofisher.com](http://www.thermofisher.com)

## AnaeroGen™

**[REF] AN0010C, AN0020C, AN0010W,  
AN0020D, AN0025A, & AN0035A**

### Intended Use

AnaeroGen (AN0010C, AN0020C, AN0010W, AN0020D, AN0025A, AN0035A) are culture media accessory devices intended to be used with the in vitro diagnostic devices part of the AnaeroGen products range, for the generation of optimal atmospheric conditions for the growth of anaerobic microorganisms.

The devices accessories are for professional use only, are not automated and nor are they companion diagnostics.

### Summary and Explanation

Anaerobic bacteria are important pathogens that can cause a variety of infections in humans. The site of anaerobic infection is commonly the site of normal colonization. Anaerobic means "life without air." Anaerobic bacteria grow in places which completely, or almost completely, lack oxygen. The spectrum of infections ranges from local abscesses to life-threatening infections. Anaerobic microorganisms, including pathogens such as *Clostridium perfringens*, *Bacteroides fragilis*, *Fusobacterium necrophorum* and *Fusobacterium nucleatum*, can be highly infectious and can cause life-threatening diseases. For example, *Fusobacterium necrophorum* is one of the most commonly isolated species of *Fusobacterium*, and is an uncommon cause of bacteremia, which has a high mortality rate of 10%<sup>1</sup>, *Fusobacterium nucleatum* is another member of the *Fusobacterium* species and has been linked to infections such as necrobacillosis and Lemierre's disease<sup>2</sup>.

Due to the severity of bacterial infections caused by anaerobic microorganisms in a vast number of patients, it is therefore highly important to grow and isolate anaerobic microorganisms from clinical samples. Early diagnosis is vital in the prevention and reduction of bacterial infections.

### Principle of Method

The AnaeroGen (AN0010C, AN0020C, AN0010W, AN0020D, AN0025A, AN0035A) can generate the optimal atmospheric conditions for the growth of anaerobic microorganisms through a very simple method.

The system consists of a paper gas generating sachet and a plastic pouch or a jar/box. The paper sachet contains ascorbic acid and activated carbon which react on contact with air. Oxygen is rapidly adsorbed and carbon dioxide is produced. When the paper sachet is placed in a sealed plastic pouch or a jar, this reaction will create ideal atmospheric conditions for the growth of anaerobes. It proceeds with no evolution of hydrogen and therefore does not require a catalyst. No addition of water is necessary to activate the reaction. When used as directed the AnaeroGen sachet will reduce the oxygen content in the pouch to below 1% within 30 minutes. The resulting carbon dioxide content will be between 8% and 14%. The level of carbon dioxide will depend on how many plates are placed in the pouch or the jar.

The system is designed to support the anaerobic incubation of culture media or identification panels.

Sachets can be provided singularly, or as part of a wider anaerobic gas generator system where 10 plastic pouches,

# Thermo

## SCIENTIFIC

and 10 W-Zip seal pouches (with integral seal) respectively, are provided together with the sachets.

AnaeroGen 2.5L Sachet (AN0025A) and AnaeroGen 3.5L Sachet (AN0035A) function according to the same principle, with the only difference being the sizes of the sachets that are intended to be used with a 2.5L jar or box and a 3.5L jar, respectively. Jar and box are sold separately.

### Description of Components

AnaeroGen Compact (AN0010C) is supplied in packages containing 10 sachets and 10 sealed plastic pouches; AnaeroGen W-Zip Compact Gas Generator System (AN0010W) includes 10 sachets and 10 W-Zip seal pouches.

AnaeroGen Compact (AN0020C) contains 20 sachets only, like the AnaeroGen Compact Sachet (AN0020D) device, that contains 10 sachets: for both these products, pouches are supplied separately.

AnaeroGen 2.5L Sachet (AN0025A) and AnaeroGen 3.5L Sachet (AN0035A) are sold in packages containing 10 sachets each, with sachets sizes of 2.5L and 3.5L respectively, intended to be used with the corresponding size 2.5L jar, 2.5L box and 3.5L jar, all sold separately.

### Materials Required but Not Supplied

- Inoculated media
- Incubators
- Quality control organisms
- AnaeroGen Compact Sealing Clip (AN0005C) A pack of 5 sealing clips intended to be used with the pouches included in AnaeroGen Compact (AN0010C and AN0020C).
- AnaeroGen Compact Plastic Pouch (AG0020C) A pack of 20 plastic pouches (without clips).
- AnaeroGen W-Zip Seal Pouches (AG0060C) A pack of 20 W-Zip seal pouches (with integral seal).
- AnaeroJar™ 2.5L Jar (AG0025A) A 2.5L jar supplied with 2.5L capacity polycarbonate jar base and lid secured to base by four self-venting clips, intended to be used with AnaeroGen 2.5L Sachet (AN0025A) and CampyGen™ 2.5L Sachet (CN0025A).
- AnaeroBox™ 2.5L Box with Lid and Seal (AB0025A) A 2.5L box intended to be used with AnaeroGen 2.5L Sachet (AN0025A) and CampyGen 2.5L Sachet (CN0025A).
- AnaeroBox 3.5L Box with Lid and Seal (AB0035A) A 3.5L box intended to be used with AnaeroGen 3.5L Sachet (AN0035A) and CampyGen 3.5L Sachet (CN0035A).
- Resazurin Anaerobic Indicator (BR0055B) A cotton strip impregnated with resazurin used as a redox indicator.

### Storage

- Store product at 2 - 25 ± 2°C.
- The product may be used until the expiry date stated on the label.

### Warnings and Precautions

- For *in vitro* diagnostic use only.
- For professional use only.
- Inspect the product packaging before first use.
- Do not use the product if there is any visible damage to the packaging or product.
- Do not use the product beyond the stated expiry date.
- Do not use the device if signs of contamination are present.

- Do not use the device if there are signs of deterioration.
- It is the responsibility of each laboratory to manage waste produced according to their nature and degree of hazard and to have them treated or disposed of in accordance with any federal, state and local applicable regulations. Directions should be read and followed carefully. This includes the disposal of used or unused reagents as well as any other contaminated disposable material following procedures for infectious or potentially infectious products.
- As soon as the AnaeroGen paper sachet is exposed to the air, the reaction will start. It is therefore essential that the paper sachet is placed in the pouch and sealed within one minute.
- The reaction of the ascorbic acid with oxygen is exothermic. However, the temperature of the AnaeroGen paper sachet will not exceed 65°C.

Refer to the Safety Data Sheet (SDS) for safe handling and disposal of the product ([www.thermofisher.com](http://www.thermofisher.com)).

### Serious Incidents

Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the relevant regulatory authority in which the user and/or the patient is established.

### Procedure

- Place the inoculated media plates or identification panel in the plastic pouch provided. Disposable plastic Petri dishes should be of the vented variety to aid gas transfer between the interior and exterior of the plates.
- Tear open an AnaeroGen foil sachet at the tear-nick indicated and remove the Anaerogen paper sachet from within.
- Immediately place the AnaeroGen paper sachet in the plastic pouch.  
**NOTE:** The AnaeroGen paper sachet will become warm to touch on exposed to air.
- Expel excess air from the plastic pouch. Seal the plastic pouch immediately with the AnaeroGen Compact Clip (AN0005).  
**NOTE:** The time taken between opening the foil sachet and sealing the plastic pouch should not exceed one minute. Extended exposure will result in loss of reactivity, and full anaerobic conditions may not be achieved in the jar.
- Incubate appropriately.
- After the incubation period remove the plates or ID panel and examine for the presence of colonies or biochemical reaction. If the plates require re-incubation, then a fresh AnaeroGen Compact sachet must be used following steps 2-5 described above.  
**NOTE:** The plates may be initially inspected through the transparent plastic pouch. If the bag is opened, a fresh AnaeroGen Compact sachet is required for re-incubation.
- After incubation, the exhausted AnaeroGen Compact paper sachet and plastic pouch should be sterilised and discarded with the non-hazardous laboratory waste.

### Quality Control

It is the responsibility of the user to perform Quality Control testing taking into account the intended use of the product, and in accordance with any local applicable regulations (frequency, number of strains, incubation temperature etc.).

It is recommended that an OXOID™ Anaerobic Indicator (BR0055) is used in the jar as a visual check that anaerobic conditions have been achieved and maintained.

The user should periodically check their anaerobic system for its ability to provide adequate conditions for the growth of appropriate bacteria.

The performance of this product can be verified by testing the following reference strains.

Incubation Conditions: 24 - 48 h @ 37° ± 2°C anaerobic

Positive Controls	
<i>Clostridium novyi</i> ATCC®9690™	Growth
Negative Controls	
<i>Kocuria rhizophila</i> ATCC®9341™	No Growth

### Performance Characteristics

Evidence of the repeatability and reproducibility of the AnaeroGen (AN0010C, AN0020C, AN0010W, AN0020D, AN0025A, AN0035A) are demonstrated by the performance data generated from In-Process QC testing performed by our third party supplier during the manufacturing of the anaerobic gas generating sachets and by the historical performance achieved through routine QC assessed by the manufacturer.

Performance Specification:

- Oxygen Removal = 0.1% or below after 2.5 hours.
- Carbon Dioxide Production = 11% +/- 4% after 24 hours.

### Bibliography

- Centers for Disease Control and Prevention. 2021a. 'Covid-19-Associated Fusobacterium nucleatum Bacteremia, Belgium'. 27(3). [https://wwwnc.cdc.gov/eid/article/27/3/20-2284\\_article](https://wwwnc.cdc.gov/eid/article/27/3/20-2284_article).
- Public Health England. 2015b. 'Identification of Anaerobic Gram-Negative Rods'. UK Standards for Microbiology Investigations. ID 25 (2). <https://www.gov.uk/government/publications/smi-id-25-identification-of-anaerobic-gram-negative-rods>.

### Symbol Legend

Symbol	Definition
	Catalogue number
	In Vitro Diagnostic Medical Device
	Batch code
	Temperature limit
	Use-by date
	Consult instructions for use or consult electronic instructions for use
	Do not use if packaging damaged and consult instructions for use

	Not for near patient testing
	Manufacturer
	Authorized representative in the European Community/ European Union
	European Conformity Assessment
	UK Conformity Assessment
	Unique device identifier

©2022 Thermo Fisher Scientific Inc. All rights reserved.  
ATCC and ATCC catalogue marks are a trademark of American Type Culture Collection.  
All other trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries.

Oxoid Limited, Wade Road, Basingstoke,  
RG24 8PW, UK

For technical assistance please contact your local distributor.

#### Revision Information

Version	Date of issue and modifications introduced
1.0	2022-11-02. New document. (LIVE)

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© 2022 Thermo Fisher Scientific Inc. Med ensamrätt.
Alla övriga varumärken tillhör Thermo Fisher Scientific Inc. och dess dotterbolag.