<https://www.thermofisher.com>

CE 2797

UK CA

IVD

Rx only

2°C

8°C

EU
US
CA
ROW

+800 135 79 135
1 855 236 0910
1 855 805 8539
+31 20 794 7071

Thermo Scientific™ Oxoid™ X and V Factor Diagnostic Discs REF (DD0003B, DD0003T, DD0004B, DD0004T, DD0005B and DD0005T)

Intended Use

The X and V Factor Diagnostic Discs are impregnated with X and/or V factor and are used in a diagnostic workflow to aid microbiology professionals in identifying potential pathogenic bacteria suspected in patients having a *Haemophilus* spp. infection. These discs enable to growth of *Haemophilus* spp. and provide a qualitative result. The device is for use with an agar grown culture, for professional use only, and is neither automated nor a companion diagnostic. This device is an *in vitro* diagnostic.

Principle of the Test

X and V Factor Diagnostic Discs are used in the qualitative agar diffusion test method for *in vitro* susceptibility testing. A suitable therapeutic agent for *in vivo* use can be determined using filter paper discs impregnated with specified concentrations of growth factor agents placed on the surface of a suitable test medium. Pure cultures of clinical isolates are inoculated onto the test medium and the Diagnostic Disc placed on the surface. The growth factors within the disc diffuses through the agar to form a gradient. After incubation, the zones of inhibition around the discs are measured and compared against recognised zone diameter ranges for the specific growth factor combination of agents/organisms under test.

Components

X and V Factor Diagnostic Discs consist of 6mm diameter paper discs impregnated with concentrations of haemin and/or nicotinamide adenine dinucleotide (NAD). The discs are clearly marked on both sides with the corresponding codes: X DD, V DD or XV DD.

Product	Disc Label	Disc Ingredient
X Factor Discs (DD0003B and DD0003T)	X DD	Haemin A dark brown powder when in dry form.
V Factor Discs (DD0004B and DD0004T)	V DD	Nicotinamide adenine dinucleotide (NAD) A white powder when in dry form.
X and V Factors Discs (DD0005B and DD0005T)	XV DD	Haemin and NAD A white powder when in dry form.

X and V Factor Diagnostic Discs are supplied in cartridges of 50 discs, either as a single cartridge per pack (DD0003T/DD0004T/DD0005T) or five cartridges per pack (DD0003B/DD0004B/DD0005B).

Cartridges are individually packed in a foil-sealed blister pack with a desiccant.

X	Factor	Diagnostic	Discs
DD0003T			50
DD0003B			250
V	Factor	Diagnostic	Discs
DD0004T			50
DD0004B			250
X and V	Factor	Diagnostic	Discs
DD0005T			50
DD0005B			250

Materials Provided for:

DD0003T/DD0004T/DD0005T:

1 x Spring cap and plunger
1 x Desiccant tablet in each cartridge
1 x Cartridge per pack
50 x discs in each cartridge
Foil sealing each individual cartridge

DD0003B/DD0004B/DD0005B:

1 x Spring cap and plunger per cartridge
1 x Desiccant tablet in each cartridge
50 x Discs in each cartridge
5 x Cartridges per pack
Foil sealing each individual cartridge

Materials Required but Not Supplied

- Agar plates with appropriate media, inoculum suspension medium, sterile loops and swabs, sterile forceps, McFarland turbidity standards, incubator, modified atmosphere environments.
- Antibiotic disc dispensers.
- Quality control (QC) strains, apparatus to measure zone sizes, and interpretative criteria for local standard methods.

Warnings and Precautions

- Not intended for susceptibility testing.
- This product is for *in vitro diagnostic* use only.
- For professional use only.
- Failure to observe instructions for use may lead to erroneous results.
- Observe aseptic techniques and established precautions against all microbiological hazards throughout all procedures.
- Cultures, containers, and other contaminated materials must be sterilized after use in accordance with guidelines for the handling and disposal of biohazardous waste.
- Refer to the Material Safety Data Sheet (MSDS) for safe handling and disposal of the product (found on www.thermofisher.com).
- Inspect the product packaging before first use. Do not use product if there is any visible damage to the foil seal.
- Do not use beyond the stated expiry date.
- Once the cartridge is open, ensure it is stored in an opaque desiccated environment to prevent degradation.
- If the discs do not produce the expected reaction with recommended control microorganisms, review the entire procedure.
- In the event of malfunction do not use the device.
- The device does not contain any carcinogenic, mutagenic or toxic to reproduction, endocrine disrupting substances or materials likely to cause sensitisation or an allergic reaction under normal use.

Storage and Handling

Unopened cartridges must be stored at 2° to 8°C until required. Allow cartridges to reach room temperature before removing them from the packaging to minimize condensation. The expiry date is valid only for the unopened blister packs stored under correct conditions. Once opened, cartridges should be stored within a dispenser (sold separately) in the container provided with its unsaturated (orange) desiccant, or other suitable opaque airtight container with a desiccant, to protect the discs from light and moisture. Once opened from their desiccant-containing packaging, discs should be used within 7 days and only if stored as described in this IFU. Dispensers should be stored within the container at 2° to 8°C and be allowed to come to room temperature before opening to prevent the formation of condensation.

Method¹

- Specimens should be collected and handled following recommended guidelines.
- Inoculate a blood agar plate using a loop or swab with a pure culture of the test organism exhibiting the appropriate haemolysis from the primary isolation plate.
- Using either sterile forceps, single disc ejector or disc dispenser; apply the disc/s to the surface of the pre-inoculated agar.
- Lightly touch the disc with the sterile forceps to ensure complete contact between the disc and agar surface.
- Invert and incubate plates at 35-37°C for 18-24 hours.
- Inspect for a zone of inhibition around the discs. If suspect colonies are isolated from the primary plating medium, repeat the described procedure with a pure culture of the suspect colonies.

Interpretation

Organisms that require only X factor – only grow around X and XV Factor discs.
Organisms that require only V factor – only grow around V and XV Factor discs.
If both X and V are required – organism grows only around XV Factor discs.

Table 1. Summary of Growth Patterns

Haemophilus species	Growth around discs		
	X	V	X + V
<i>H. influenzae</i>	-	-	+
<i>H. parainfluenzae</i>	-	+	+
<i>H. haemoglobinophilus</i>	+	-	+
<i>H. aegyptius</i> *	-	-	+
<i>H. haemolyticus</i>	-	-	+
<i>H. pittmaniae</i>	-	+	+
<i>H. parahaemolyticus</i>	-	+	+
<i>H. paraphrohaemolyticus</i>	-	+	+
<i>H. ducreyi</i>	+	-	+
<i>H. sputorum</i>	-	+	+
* <i>H. aegyptius</i> is indistinguishable from <i>H. influenzae</i> biotype III in normal laboratory tests. ²			

User Quality Control Procedures

The use of positive and negative control organisms is recommended.¹

Each batch of X or V Factor discs should be checked with a positive and a negative control.

Each batch of XV Factor discs should be checked with a positive control.

Testing of control organisms should be performed in accordance with established laboratory quality control procedures. If aberrant quality control results are noted, patient results should not be reported.

X and V Factor	<i>Haemophilus influenzae</i> NCTC 11931 or NCTC 12975
V Factor only	<i>Haemophilus parainfluenzae</i> NCTC 10665
X Factor only	<i>Haemophilus haemoglobinophilus</i> NCTC 8540

Analytical Performance Characteristics

Repeatability and reproducibility were measured internally by testing *Haemophilus influenzae* (ATCC® 19418, ATCC® 9334, and ATCC® 49247) and *Haemophilus parainfluenzae* (ATCC® 33392) on X and V Factor Diagnostic Discs. Refer to Table 2.

Results were deemed comparable when either no deviation, or a deviation no greater than 2 mm, between the test discs was observed. (Data on file.)

Table 2. Batch Pass Rate

Product	Batches tested	Time Period	Pass Rate
X Factor DD0003B	176	2001 – 2021	99.43%
X Factor DD0003T	148	1998 – 2021	95.27%
V Factor DD0004B	163	2001 – 2021	95.09%
V Factor DD0004T	118	1998 – 2021	98.31%
X and V Factor DD0005B	113	2001 – 2021	95.58%
X and V Factor DD0005T	170	1998 – 2021	97.06%

Clinical Performance Characteristics

A total of 97 isolates (88 throat swabs; nine reference strains) were collected from 15 healthy men isolated at the British Antarctic Survey Base, Antarctica. Swabs were inoculated and cultured.⁴ *Haemophilus* spp. like bacterial colonies were sub-cultured and confirmed using Oxoid V Factor Diagnostic Discs. Further speciation of *Haemophilus* spp. was determined using Oxoid X Factor Discs and/or Oxoid Factor V Diagnostic

amongst the 15 healthy men and exhibited transmission within the isolated community. The study demonstrated that Oxoid X Factor and V factor Diagnostic Discs are effective in the identification of *H. parainfluenzae*.

Cleary et al. used Oxoid X and V Factor Diagnostic Discs to identify *H. influenzae* pathogens from nasopharyngeal samples of young children in the UK. 1,569 nasopharyngeal swabs were collected.⁵ *H. influenzae* was identified using Oxoid X and V Factor Diagnostic Discs in 275 of the isolates; of that group, 273 (99.3%) were confirmed as non-typable *H. influenzae* (NTHi) by PCR. The carriage rate of NTHi among the children in 2008 to 2009 was 14.6%, and in 2009 to 2010, it was 14.8%. The PCV13 vaccine was introduced in 2010. The NTHi carriage rate from 2010 to 2011 was 22.7%, which was significantly higher than the previous years ($P < 0.05$). The carriage rate from 2011 to 2012 was 16.6% and from 2012 to 2013 was 21.52%. This study demonstrated that the introduction of the PCV13 vaccine in young children altered the epidemiology of NTHi and it demonstrated that indirect negative and positive associations between NTHi and *S. pneumoniae* carriage exists. In particular, at the time of vaccine introduction, the carriage rate of NTHi significantly increased. This study did not report any inconsistent or out of specification results for the Oxoid X and V.

Serious Incidents

Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the relevant regulatory authority in which the user and/or the patient is established.

Discs. Tests identified 67 of the isolates as *H. parainfluenzae* with up to five strains being isolated from one swab. This study demonstrated a high prevalence of *H. parainfluenzae*

Contact Information

International Telephone: +44 (0) 1256 841144

Email: microbiology@thermofisher.com

Website: www.thermofisher.com

References

- Public Health England. 2019. "UK Standards for Microbiology Investigations – X and F Factor Test." https://assets.publishing.service.gov.uk/government/uploads/system/uploads/attachment_data/file/800302/TP_38i4.pdf
- MacFaddin JF. X and V Factors. Biochemical Tests for Identification of Medical Bacteria. 3rd ed. Philadelphia: Lippincott Williams and Wilkins; 2000. p. 451-3. B, III.
- Jones AM. Haemophilus influenzae and H. parainfluenzae: the influence of media and CO2 on differentiation using X, V and XV discs. *MedLab Sci* 1982;39:189-91. B, III.
- Kerr, G, K Forbes, A Williams, and T Pennington. 1993. "An Analysis of the Diversity of Haemophilus Parainfluenzae in the Adult Human Respiratory Tract by Genomic DNA Fingerprinting." *Epidemiology and Infection* 111 (1): 89–98. <https://doi.org/10.1017/S0950268800056715>.
- Cleary, D, V Devine, D Morris, K Oxman, R Gladstone, S Bentley, S Faust, and S Clarke. 2018. "Pneumococcal Vaccine Impacts on the Population Genomics of Non-Typeable Haemophilus Influenzae." *Microbial Genomics* 4 (9). <https://doi.org/10.1099/MGEN.0.000209>.

Glossary of Symbols

Symbol/Label	Meaning
	Manufacturer
IVD	In Vitro Diagnostic Medical Device
	Temperature limit
LOT	Batch Code
REF	Catalog Number
	Do not re-use
	Consult instructions for use or consult electronic instructions for use
	Contains sufficient for <n> tests
	Use-by date
	Do not use if package is damaged and Consult instructions for use
EC REP	Authorized representative in the European Community/ European Union
UDI	Unique device identifier
R_x only	USA: Caution: Federal law restricts this device to sale by or on order of a Physician
CE	European Conformity Mark
UK CA	UK Conformity Mark

© 2022 Thermo Fisher Scientific Inc. All rights reserved.
 ATCC® is a trademark of American Type Culture Collection.
 All other trademarks are the property of Thermo Fisher Scientific and its subsidiaries.
 This information is not intended to encourage use of these products in any manner that might infringe the intellectual property rights of others. Please contact your local sales representative for all enquiries.

Oxoid Limited, Wade Road, Basingstoke, Hampshire, RG24 8PW, England

For technical assistance please contact your local distributor.

Revision Information

Version	Date of issue and modifications introduced
1.0	2022-11-01. New document

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Podaci za kontakt

Međunarodni telefon: +44 (0) 1256 841144

E-pošta: microbiology@thermofisher.com

Internetska stranica: www.thermofisher.com

Literatura

- Public Health England. 2019. "UK Standards for Microbiology Investigations – X and F Factor Test." https://assets.publishing.service.gov.uk/government/uploads/system/uploads/attachment_data/file/800302/TP_38i4.pdf
- MacFaddin JF. X and V Factors. Biochemical Tests for Identification of Medical Bacteria. 3rd ed. Philadelphia: Lippincott Williams and Wilkins; 2000. p. 451-3. B, III.
- Jones AM. *Haemophilus influenzae* and *H. parainfluenzae*: the influence of media and CO₂ on differentiation using X, V and XV discs. *MedLab Sci* 1982;39:189-91. B, III.
- Kerr, G, K Forbes, A Williams, and T Pennington. 1993. "An Analysis of the Diversity of *Haemophilus Parainfluenzae* in the Adult Human Respiratory Tract by Genomic DNA Fingerprinting." *Epidemiology and Infection* 111 (1): 89–98. <https://doi.org/10.1017/S0950268800056715>.
- Cleary, D, V Devine, D Morris, K Oxman, R Gladstone, S Bentley, S Faust, and S Clarke. 2018. "Pneumococcal Vaccine Impacts on the Population Genomics of Non-Typeable *Haemophilus Influenzae*." *Microbial Genomics* 4 (9). <https://doi.org/10.1099/MGEN.0.000209>.

Rječnik simbola

Simbol/oznaka	Značenje
	Proizvođač
	Medicinski proizvod za dijagnostiku in vitro
	Gornja granica temperature
	Šifra serije
	Kataloški broj
	Nemojte ponovno upotrebljavati
	Proučite upute za upotrebu ili pogledajte elektroničke upute za upotrebu
	Sadrži dovoljnu količinu za <n> testova
	Rok upotrebe
	Nemojte upotrebljavati ako je pakiranje oštećeno i pogledajte upute za upotrebu
	Ovlašteni zastupnik u Europskoj zajednici / Europskoj uniji
	Jedinstvena identifikacijska oznaka proizvoda
	SAD: Opres: Saveznim zakonom određeno je da proizvod smiju prodavati samo liječnici / da se proizvod smije prodavati samo po nalogu liječnika
	Europska oznaka sukladnosti
	Oznaka sukladnosti za UK

© 2022. Thermo Fisher Scientific Inc. Sva su prava pridržana.

Oznaka ATCC® zaštitni je znak Američke zbirke tipskih kultura.

Svi drugi zaštitni znakovi vlasništvo su društva Thermo Fisher Scientific i njegovih podružnica. Ove informacije nisu namijenjene poticanju upotrebe ovih proizvoda na bilo koji način koji bi mogao dovesti do povrede prava intelektualnog vlasništva drugih. U slučaju bilo kakvih pitanja obratite se svom lokalnom prodajnom predstavniku.

Oxoid Limited, Wade Road, Basingstoke, Hampshire, RG24 8PW, Engleska

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Kontaktteave

Rahvusvaheline telefoninumber: +44 (0) 12 5684 1144

E-post: microbiology@thermofisher.com

Veebileht: www.thermofisher.com

Viited

- Public Health England. 2019. "UK Standards for Microbiology Investigations – X and F Factor Test." https://assets.publishing.service.gov.uk/government/uploads/system/uploads/attachment_data/file/800302/TP_38i4.pdf
- MacFaddin JF. X and V Factors. Biochemical Tests for Identification of Medical Bacteria. 3rd ed. Philadelphia: Lippincott Williams and Wilkins; 2000. p. 451–3. B, III.
- Jones AM. Haemophilus influenzae and H. parainfluenzae: the influence of media and CO₂ on differentiation using X, V and XV discs. *MedLab Sci* 1982;39:189–91. B, III.
- Kerr, G, K Forbes, A Williams, and T Pennington. 1993. "An Analysis of the Diversity of Haemophilus Parainfluenzae in the Adult Human Respiratory Tract by Genomic DNA Fingerprinting." *Epidemiology and Infection* 111 (1): 89–98. <https://doi.org/10.1017/S0950268800056715>.
- Cleary, D, V Devine, D Morris, K Oxman, R Gladstone, S Bentley, S Faust, and S Clarke. 2018. "Pneumococcal Vaccine Impacts on the Population Genomics of Non-Typeable Haemophilus Influenzae." *Microbial Genomics* 4 (9). <https://doi.org/10.1099/MGEN.0.000209>.

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Oxoid Limited, Wade Road, Basingstoke, Hampshire, RG24 8PW, Inglismaa

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E-mail : microbiology@thermofisher.com


Site Web : www.thermofisher.com

Références

- Public Health England. 2019. "UK Standards for Microbiology Investigations – X and F Factor Test." https://assets.publishing.service.gov.uk/government/uploads/system/uploads/attachment_data/file/800302/TP_38i4.pdf
- MacFaddin JF. Facteurs X et V. Biochemical Tests for Identification of Medical Bacteria. 3rd ed. Philadelphia : Lippincott Williams and Wilkins ; 2 000. p. 451-3. B, III.
- Jones AM. *Haemophilus influenzae* et *H. parainfluenzae*: the influence of media and CO₂ on differentiation using X, V and XV discs. *MedLab Sci* 1982;39:189-91. B, III.
- Kerr, G, K Forbes, A Williams et T Pennington. 1993. "An Analysis of the Diversity of *Haemophilus Parainfluenzae* in the Adult Human Respiratory Tract by Genomic DNA Fingerprinting." *Epidemiology and Infection* 111 (1): 89–98. <https://doi.org/10.1017/S0950268800056715>.
- Cleary, D, V Devine, D Morris, K Oxman, R Gladstone, S Bentley, S Faust et S Clarke. 2018. "Pneumococcal Vaccine Impacts on the Population Genomics of Non-Typeable *Haemophilus Influenzae*." *Microbial Genomics* 4 (9). <https://doi.org/10.1099/MGEN.0.000209>.

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Oxoid Limited, Wade Road, Basingstoke, Hampshire, RG24 8PW, Angleterre

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Kontaktinformationen

Internationales Telefon: +44 (0) 1256 841144

E-Mail: microbiology@thermofisher.com



Website: www.thermofisher.com

Literatur

1. Public Health England. 2019. „UK Standards for Microbiology Investigations – X and F Factor Test.“ https://assets.publishing.service.gov.uk/government/uploads/system/uploads/attachment_data/file/800302/TP_38i4.pdf
2. MacFaddin, J.F. X and V Factors. Biochemical Tests for Identification of Medical Bacteria. 3. Aufl. Philadelphia: Lippincott Williams und Wilkins; 2000. p. 451-3. B, III.
3. Jones AM. *Haemophilus influenzae* und *H. parainfluenzae*: the influence of media and CO₂ on differentiation using X, V and XV discs. *MedLab Sci* 1982;39:189-91. B, III.
4. Kerr, G, K Forbes, A Williams und T Pennington. 1993. „An Analysis of the Diversity of *Haemophilus parainfluenzae* in the Adult Human Respiratory Tract by Genomic DNA Fingerprinting.“ *Epidemiology and Infection* 111 (1): 89–98. <https://doi.org/10.1017/S0950268800056715>.
5. Cleary, D., V. Devine, D. Morris, K. Oxman, R. Gladstone, S. Bentley, S. Faust und S. Clarke. 2018. „Pneumococcal Vaccine Impacts on the Population Genomics of Non-Typeable *Haemophilus influenzae*.“ *Microbial Genomics* 4 (9). <https://doi.org/10.1099/MGEN.0.000209>.

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1. Public Health England. 2019. "UK Standards for Microbiology Investigations – X and F Factor Test." https://assets.publishing.service.gov.uk/government/uploads/system/uploads/attachment_data/file/800302/TP_38i4.pdf
2. MacFaddin JF. X and V Factors. Biochemical Tests for Identification of Medical Bacteria. 3η έκδ. Philadelphia: Lippincott Williams και Wilkins· 2000. σελ. 451-3. Β, ΙΙΙ.
3. Jones AM. *Haemophilus influenzae* and *H. parainfluenzae*: the influence of media and CO₂ on differentiation using X, V and XV discs. *MedLab Sci* 1982;39:189-91. Β, ΙΙΙ.
4. Kerr, G, K Forbes, A Williams και T Pennington. 1993. "An Analysis of the Diversity of *Haemophilus Parainfluenzae* in the Adult Human Respiratory Tract by Genomic DNA Fingerprinting." *Epidemiology and Infection* 111 (1): 89–98. <https://doi.org/10.1017/S0950268800056715>.

5. Cleary, D, V Devine, D Morris, K Oxman, R Gladstone, S Bentley, S Faust και S Clarke. 2018. "Pneumococcal Vaccine Impacts on the Population Genomics of Non-Typeable *Haemophilus Influenzae*." *Microbial Genomics* 4 (9). <https://doi.org/10.1099/MGEN.0.000209>.

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Sito web: www.thermofisher.com

Bibliografia

- Public Health England. 2019. "UK Standards for Microbiology Investigations – X and F Factor Test." https://assets.publishing.service.gov.uk/government/uploads/system/uploads/attachment_data/file/800302/TP_38i4.pdf
- MacFaddin JF. X and V Factors. Biochemical Tests for Identification of Medical Bacteria. 3a ed. Philadelphia: Lippincott Williams & Wilkins; 2000. p. 451-3. B, III.
- Jones AM. *Haemophilus influenzae* and *H. parainfluenzae*: the influence of media and CO₂ on differentiation using X, V and XV discs. *MedLab Sci* 1982;39:189-91. B, III.
- Kerr, G, K Forbes, A Williams e T Pennington. 1993. "An Analysis of the Diversity of *Haemophilus parainfluenzae* in the Adult Human Respiratory Tract by Genomic DNA Fingerprinting." *Epidemiology and Infection* 111 (1): 89-98. <https://doi.org/10.1017/S0950268800056715>.
- Cleary, D, V Devine, D Morris, K Oxman, R Gladstone, S Bentley, S Faust e S Clarke. 2018. "Pneumococcal Vaccine Impacts on the Population Genomics of Non-Typeable *Haemophilus influenzae*." *Microbial Genomics* 4 (9). <https://doi.org/10.1099/MGEN.0.000209>.

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Oxoid Limited, Wade Road, Basingstoke, Hampshire, RG24 8PW, Inghilterra

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- Public Health England. 2019. “UK Standards for Microbiology Investigations – X and F Factor Test.” https://assets.publishing.service.gov.uk/government/uploads/system/uploads/attachment_data/file/800302/TP_38i4.pdf
- MacFaddin JF. X and V Factors. Biochemical Tests for Identification of Medical Bacteria. 3rd ed. Philadelphia: Lippincott Williams and Wilkins; 2000. p. 451–3. B, III.
- Jones AM. Haemophilus influenzae and H. parainfluenzae: the influence of media and CO2 on differentiation using X, V and XV discs. MedLab Sci 1982;39:189–91. B, III.
- Kerr, G, K Forbes, A Williams, and T Pennington. 1993. “An Analysis of the Diversity of Haemophilus Parainfluenzae in the Adult Human Respiratory Tract by Genomic DNA Fingerprinting.” Epidemiology and Infection 111 (1): 89–98. <https://doi.org/10.1017/S0950268800056715>.
- Cleary, D, V Devine, D Morris, K Oxman, R Gladstone, S Bentley, S Faust, and S Clarke. 2018. “Pneumococcal Vaccine Impacts on the Population Genomics of Non-Typeable Haemophilus Influenzae.” Microbial Genomics 4 (9). <https://doi.org/10.1099/MGEN.0.000209>.

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Oxoid Limited, Wade Road, Basingstoke, Hampshire, RG24 8PW, Anglija

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Strona internetowa: www.thermofisher.com

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1. Public Health England. 2019. „UK Standards for Microbiology Investigations – X and F Factor Test.” https://assets.publishing.service.gov.uk/government/uploads/system/uploads/attachment_data/file/800302/TP_38i4.pdf
2. MacFaddin JF. X and V Factors. Biochemical Tests for Identification of Medical Bacteria. 3rd ed. Philadelphia: Lippincott Williams and Wilkins; 2000. p. 451–3. B, III.
3. Jones AM. *Haemophilus influenzae* and *H. parainfluenzae*: the influence of media and CO₂ on differentiation using X, V and XV discs. *MedLab Sci* 1982;39:189–91. B, III.
4. Kerr, G, K Forbes, A Williams i T Pennington. 1993. „An Analysis of the Diversity of *Haemophilus Parainfluenzae* in the Adult Human Respiratory Tract by Genomic DNA Fingerprinting.” *Epidemiology and Infection* 111 (1): 89–98. <https://doi.org/10.1017/S0950268800056715>
5. Cleary, D, V Devine, D Morris, K Oxman, R Gladstone, S Bentley, S Faust i S Clarke. 2018. „Pneumococcal Vaccine Impacts on the Population Genomics of Non-Typeable *Haemophilus Influenzae*.” *Microbial Genomics* 4 (9). <https://doi.org/10.1099/MGEN.0.000209>.

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Site: www.thermofisher.com

Referências

- Public Health England. 2019. "UK Standards for Microbiology Investigations – X and F Factor Test." https://assets.publishing.service.gov.uk/government/uploads/system/uploads/attachment_data/file/800302/TP_38i4.pdf
- MacFaddin JF. X and V Factors. Biochemical Tests for Identification of Medical Bacteria. 3rd ed. Philadelphia: Lippincott Williams and Wilkins; 2000. p. 451-3. B, III.
- Jones AM. *Haemophilus influenzae* and *H. parainfluenzae*: the influence of media and CO₂ on differentiation using X, V and XV discs. *MedLab Sci* 1982;39:189-91. B, III.

4. Kerr, G, K Forbes, A Williams, and T Pennington. 1993. "An Analysis of the Diversity of Haemophilus Parainfluenzae in the Adult Human Respiratory Tract by Genomic DNA Fingerprinting." *Epidemiology and Infection* 111 (1): 89–98. <https://doi.org/10.1017/S0950268800056715>.
5. Cleary, D, V Devine, D Morris, K Oxman, R Gladstone, S Bentley, S Faust, and S Clarke. 2018. "Pneumococcal Vaccine Impacts on the Population Genomics of Non-Typeable Haemophilus Influenzae." *Microbial Genomics* 4 (9). <https://doi.org/10.1099/MGEN.0.000209>.

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Telefon internațional: +44 (0) 1256 841144

E-mail: microbiology@thermofisher.com






Website: www.thermofisher.com

Referințe

1. Sănătatea Publică din Anglia. 2019. "UK Standards for Microbiology Investigations – X and F Factor Test." https://assets.publishing.service.gov.uk/government/uploads/system/uploads/attachment_data/file/800302/TP_38i4.pdf
2. MacFaddin JF. X and V Factors. Biochemical Tests for Identification of Medical Bacteria. 3rd ed. Philadelphia: Lippincott Williams and Wilkins; 2000. p. 451-3. B, III.
3. Jones AM. Haemophilus influenzae and H. parainfluenzae: the influence of media and CO2 on differentiation using X, V and XV discs. *MedLab Sci* 1982;39:189-91. B, III.
4. Kerr, G, K Forbes, A Williams, and T Pennington. 1993. "An Analysis of the Diversity of Haemophilus Parainfluenzae in the Adult Human Respiratory Tract by Genomic DNA Fingerprinting." *Epidemiology and Infection* 111 (1): 89–98. <https://doi.org/10.1017/S0950268800056715>.
5. Cleary, D, V Devine, D Morris, K Oxman, R Gladstone, S Bentley, S Faust, and S Clarke. 2018. "Pneumococcal Vaccine Impacts on the Population Genomics of Non-Typeable Haemophilus Influenzae." *Microbial Genomics* 4 (9). <https://doi.org/10.1099/MGEN.0.000209>.

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Oxoid Limited, Wade Road, Basingstoke, Hampshire, RG24 8PW, Anglia

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EU
USA
CA
Zvyšok sveta

+800 135 79 135
1 855 236 0910
1 855 805 8539
+31 20 794 7071

Thermo Scientific™ Oxoid™ X and V Factor Diagnostic Discs

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Odkazy

- Public Health England. 2019. „UK Standards for Microbiology Investigations – X and F Factor Test.“ https://assets.publishing.service.gov.uk/government/uploads/system/uploads/attachment_data/file/800302/TP_38i4.pdf
- MacFaddin JF. X and V Factors. Biochemical Tests for Identification of Medical Bacteria. 3. vyd.. Philadelphia: Lippincott Williams and Wilkins; 2000. str. 451-3. B, III.
- Jones AM. *Haemophilus influenzae* and *H. parainfluenzae*: the influence of media and CO₂ on differentiation using X, V and XV discs. *MedLab Sci* 1982;39:189-91. B, III.
- Kerr, G, K Forbes, A Williams, a T Pennington. 1993. „An Analysis of the Diversity of *Haemophilus Parainfluenzae* in the Adult Human Respiratory Tract by Genomic DNA Fingerprinting.“ *Epidemiology and Infection* 111 (1): 89–98. <https://doi.org/10.1017/S0950268800056715>.
- Cleary, D, V Devine, D Morris, K Oxman, R Gladstone, S Bentley, S Faust, a S Clarke. 2018. „Pneumococcal Vaccine Impacts on the Population Genomics of Non-Typeable *Haemophilus Influenzae*.“ *Microbial Genomics* 4 (9). <https://doi.org/10.1099/MGEN.0.000209>.

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Oxoid Limited, Wade Road, Basingstoke, Hampshire, RG24 8PW, Anglicko

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Website: www.thermofisher.com

Referenser

- Public Health England. 2019. "UK Standards for Microbiology Investigations – X and F Factor Test." https://assets.publishing.service.gov.uk/government/uploads/system/uploads/attachment_data/file/800302/TP_38i4.pdf
- MacFaddin JF. X and V Factors. Biochemical Tests for Identification of Medical Bacteria. 3rd ed. Philadelphia: Lippincott Williams and Wilkins; 2000. p. 451-3. B, III.
- Jones AM. Haemophilus influenzae and H. parainfluenzae: the influence of media and CO2 on differentiation using X, V and XV discs. *MedLab Sci* 1982;39:189-91. B, III.
- Kerr, G, K Forbes, A Williams och T Pennington. 1993. "An Analysis of the Diversity of Haemophilus Parainfluenzae in the Adult Human Respiratory Tract by Genomic DNA Fingerprinting." *Epidemiology and Infection* 111 (1): 89–98. <https://doi.org/10.1017/S0950268800056715>.
- Cleary, D, V Devine, D Morris, K Oxman, R Gladstone, S Bentley, S Faust och S Clarke. 2018. "Pneumococcal Vaccine Impacts on the Population Genomics of Non-Typeable Haemophilus Influenzae." *Microbial Genomics* 4 (9). <https://doi.org/10.1099/MGEN.0.000209>.

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©2022 Thermo Fisher Scientific Inc. Med ensamrätt.