[www.thermofisher.com](http://www.thermofisher.com)

## Columbia CNA Aesculin Selective Agar

**REF PB5115A, PB5267E\* and PB5224E\*\***

\* This IFU is intended to be read in conjunction with the IFU for *Brilliance™* UTI Clarity (PO5159A).

\*\*This IFU is intended to be read in conjunction with the IFU for MacConkey Agar No.3, Mod (PO5140A)

Additional IFUs available at [www.thermofisher.com](http://www.thermofisher.com)

### Intended Use

Columbia CNA Aesculin Selective Agar (PB5115A, PB5267E and PB5224E) devices are selective media for the growth of Gram-positive microorganisms with clearly visible haemolysis from clinical samples (e.g. wound, throat, genital, nose, and groin swabs). The devices are used in a diagnostic workflow to aid clinicians in determining potential treatment options for patients suspected of having bacterial infections.

The devices are for professional use only, are not automated, nor are they companion diagnostics.

### Summary and Explanation

Gram-positive bacteria are found as regular commensal flora in the environment. However, some Gram-positive bacteria are pathogenic to humans and can present as opportunistic pathogens. Such bacteria cause invasive infections of the skin, bloodstream, gastrointestinal tract and urinary tract infections (UTI). Such infections can result in severe and lethal infections, particularly those that are nosocomial following surgical procedures or hospital stays. Resistance to antibiotics is a problem among many Gram-positive bacterial infections, such as penicillin resistance, methicillin resistance and vancomycin resistance. The emergence of these resistances, or presence of, during infection can further complicate treatment protocols and impact patient outcomes.

Gram-positive bacteria such as *Staphylococcus aureus*, *Streptococcus agalactiae* and *Enterococcus faecalis* are such opportunistic pathogens that can readily acquire or have antibiotic resistance. *S. aureus* is a Gram-positive coccus commonly found as a commensal of the human upper respiratory tract and skin microbiomes<sup>1,2</sup>. *S. aureus* is an opportunistic pathogen, entering the body through cuts or abrasions where it can cause a spectrum of disease in both immunocompetent and immunocompromised hosts<sup>1</sup>. This can range from superficial skin infections such as cellulitis and abscesses to life-threatening illnesses such as pneumonia, toxic shock syndrome and sepsis<sup>1</sup>. *S. aureus* is one of the leading causes of hospital acquired infections (HAIs), further complicated by the rapid emergence of antibiotic resistance in this pathogen<sup>2</sup>. Strains of *S. aureus* that have developed multiple drug-resistance to beta-lactam antibiotics, termed methicillin-resistant *S. aureus* (MRSA), are a significant problem in hospitals as well as the community<sup>2</sup>. As such, the surveillance and reporting of methicillin-sensitive *S. aureus* (MSSA), as well as MRSA, is now mandatory in UK hospitals<sup>3</sup>.

*S. agalactiae* can be responsible for invasive disease in pregnant women (urinary tract infection, amnionitis,

# Thermo

## SCIENTIFIC

endometritis) and in new-borns, potentially causing both maternal and neonatal disease (i.e. infection occurring during the first four weeks of life)<sup>4,5</sup>. *S. agalactiae* can also cause invasive disease in non-pregnant adults with comorbidities causing immunodeficiency, with the elderly particularly at risk of infection<sup>6</sup>.

*Enterococcus faecalis* is thought to be a significant nosocomial pathogen due to its ability to colonise the gut without resulting in symptomatic infection<sup>7-9</sup>. Such asymptomatic carriage can facilitate rapid spread in healthcare settings, complicating infection control and treatment. Previous antibiotic treatment, prolonged hospitalisation, immunosuppression, and comorbidities such as diabetes are all thought to contribute to an increased likelihood of colonisation by, and potential infection with, vancomycin-resistant *Enterococcus*.

Consequently, it is important to distinguish between Gram-positive organisms in clinical samples to ensure that the correct diagnosis is achieved, and a suitable treatment protocol can be implemented for the patient. Gram-positive bacteria are easily identifiable following Gram stain where the crystal violet stain will bind to the thick peptidoglycan layer present in the bacterial cell wall and stain the bacteria a purple colour under the microscope. Following a Gram stain, a morphological analysis of the bacteria can help to further differentiate genus and species, such as whether the bacteria appear coccoid or bacillus in shape. Growth on blood agar can help to further differentiate species, *S. aureus* will produce  $\beta$ -hemolysis<sup>2</sup> on blood agar, *S. agalactiae* can also present as  $\beta$ -hemolysis<sup>10</sup> but would test negative for the catalase test; *Enterococcus faecalis* would not present any hemolysis on Blood agar, but would produce a dark halo around colonies under UV light<sup>10</sup>.

### Principle of Method

Special peptone supply nutrients for growth. Starch is added to absorb any toxic metabolites. Sodium chloride maintains the osmotic equilibrium and agar is a solidifying agent. Sheep blood is a nutrient source and allows the determination of the haemolytic properties of the isolate. Colistin sulphate is a polymyxin antibiotic that is active against Gram-negative bacteria including *Pseudomonas aeruginosa*; *Proteus* species however are not inhibited. Nalidixic acid is also active against many Gram-negative bacteria including *Proteus* species. Aesculin is added to differentiate between aesculin negative and aesculin positive bacteria. Positive organisms produce a dark halo around the colony in UV light at 366 nm.

### Typical Formula

	grams per litre
Special peptone	23.0
Starch	1.0
Sodium chloride	5.0
Aesculin	1.0
Agar	10.0
Nalidixic acid	0.005
Colistin	0.0075
Defibrinated sheep blood	50.0 ml

### Physical Appearance

Colour	Flame red
Clarity	Opaque
Fill weight	19 ± 2.0g
pH	7.1 – 7.5

### Materials Provided

PB5115A: 10 x 90mm Columbia CNA Aesculin Selective Agar plates

PB5267E: 10 x 90 mm Columbia CNA Aesculin Selective Agar / *Brilliance* UTI Clarity biplates

PB5224E: 10 x 90 mm Columbia CNA Aesculin Selective Agar / MacConkey Agar No.3, Mod biplates

Each plate should only be used once.  
Each pack contains enough plates for 10 individual tests.

## Materials Required but Not Supplied

- 1) Inoculating loops
- 2) Swabs
- 3) Collection containers
- 4) Incubators
- 5) Quality control organisms

## Storage

- Store product in its original packaging at 2–12°C until used.
- The product may be used until the expiry date stated on the label.
- Store away from light.
- Allow product to equilibrate to room temperature before use.
- Do not incubate prior to use.

## Warnings and Precautions

- For *in vitro* diagnostic use only.
- For professional use only.
- Inspect the product packaging before first use.
- Do not use the product if there is any visible damage to the packaging or plates.
- Do not use the product beyond the stated expiry date.
- Do not use the device if signs of contamination are present.
- Do not use the device if the colour has changed or there are other signs of deterioration.
- Each plate is single use and should only be used once.
- It is the responsibility of each laboratory to manage waste produced according to their nature and degree of hazard and to have them treated or disposed of in accordance with any federal, state and local applicable regulations. Directions should be read and followed carefully. This includes the disposal of used or unused reagents as well as any other contaminated disposable material following procedures for infectious or potentially infectious products.

Refer to the Safety Data Sheet (SDS) for safe handling and disposal of the product ([www.thermofisher.com](http://www.thermofisher.com)).

## Serious Incidents

Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the relevant regulatory authority in which the user and/or the patient is established.

## Specimen Collection, Handling and Storage

Specimen should be collected and handled following the recommended guidelines, such as the UK Standards for Microbiology Investigations (UK SMI) B 29 (Public Health England, 2020).

## Procedure

- (1) Allow product to equilibrate to room temperature.
- (2) Inoculate the specimen onto the medium and streak to obtain isolated colonies.
- (3) If the material is being cultured from a swab, roll the swab over a small surface area and streak for the isolation of single colonies.
- (4) Incubate plates aerobically for 24 hours at  $37 \pm 2^\circ\text{C}$ .
- (5) Inspect plates to assess colony growth and colour under good lighting.

## Interpretation

Most plates will have an area of confluent growth after incubation. On Columbia CNA Aesculin Selective Agar, typical colony morphology is as follows:

The presence of shiny, light grey colonies, aesculin positive indicate the sample is *Enterococcus faecalis*.

Shiny white colonies indicate *Staphylococcus aureus*.

Shiny grey colonies indicate *Streptococcus agalactiae*.

Identification is presumptive and should be confirmed.

## Quality Control

It is the responsibility of the user to perform Quality Control testing taking into account the intended use of the medium, and in accordance with any local applicable regulations (frequency, number of strains, incubation temperature etc.).

The performance of this medium can be verified by testing the following reference strains.

### Incubation Conditions: 24 h @ 37°C, aerobic

Positive Controls	
Inoculum $20 - 110$ colony forming units (cfu)	
<i>Enterococcus faecalis</i> ATCC® 29212	Light grey, shiny colonies without haemolysis, aesculin positive
Positive Controls	
Inoculum $10^3 - 10^4$ cfu	
<i>Staphylococcus aureus</i> ATCC® 25923	White, shiny colonies.
<i>Staphylococcus agalactiae</i> ATCC® 13813	Grey shiny colonies, aesculin negative.
Negative Controls	
Inoculum $10^4 - 10^5$ cfu	
<i>Escherichia coli</i> ATCC® 25922	Complete inhibition
<i>*Pseudomonas aeruginosa</i> ATCC® 27853	Complete inhibition

*\*Pseudomonas aeruginosa* is a negative control for PB5115A only.

A satisfactory result is represented by recovery of positive strains equal to or greater than 50% of the control medium. Negative strains are restricted in growth or are inhibited.

## Limitations

The haemolytic reaction will be strain dependent; it may also vary for some microorganisms with the type of blood used. All identifications are presumptive and should be confirmed using appropriate methods. Due to variation in nutritional requirements or sensitivity to selective agents, some strains of the target organisms may be encountered that poorly grow or fail to grow on this medium. Incubation in a carbon dioxide enriched atmosphere may cause inhibition of staphylococci. Microorganisms other than streptococci and staphylococci which are resistant to the antimicrobial agents used in the medium will be able to grow.

## Performance Characteristics

Accuracy has been demonstrated through review of the QC data. Correct growth of Gram-positive microorganisms is confirmed by the inclusion of well-characterised isolates in the QC processes performed as part of the manufacture of each batch of the devices. The precision of measurement

was demonstrated by an overall pass rate of 100% for all products (PB5115A: 10 batches between 15.10.2021 and 09.12.2021; PB5267E: 10 batches between 26.10.2021 to 15.12.2021; PB5224E: 10 batches between 03.01.2022 to 08.02.2022). This shows that the performance is reproducible.

The devices are tested in-house as part of the QC process since the products were launched in 2008 for PB5115A and PB5224E and in 2020 for PB5267E. For target organisms, when using 20-110 cfu (or  $10^3$  -  $10^4$  cfu for PB5267E and PB5224E) of *Enterococcus faecalis*, and  $10^3$  -  $10^4$  cfu of *Streptococcus agalactiae* and *Staphylococcus aureus* and incubating the device at  $36 \pm 1^\circ\text{C}$  for 24 hours in aerobic conditions, the user can recover organisms with colony size and morphology as detailed in this document. For non-target organisms, when using  $>10^4$  cfu for *Escherichia coli* and *Pseudomonas aeruginosa* (PB5115A only) and incubating the device at  $36 \pm 1^\circ\text{C}$  for 24 hours in aerobic conditions, the users can recover organisms with colony size and morphology as detailed in this document.

### Bibliography

1. Foster, Timothy J. 1996. "Staphylococcus." In Medical Microbiology. 4th Edition., 4th ed. Galveston.
2. Public Health England. 2020. 'Investigation of Specimens for Screening for MRSA'. UK Standards for Microbiology Investigations. B 29 (7). <https://www.gov.uk/government/publications/smi-b-29-investigation-of-specimens-for-screening-for-mrsa>.
3. National Institute for Health and Care Excellence. 2016. 'Healthcare-Associated Infections'. QS113. [www.nice.org.uk/guidance/qs113](http://www.nice.org.uk/guidance/qs113).
4. Melin, Pierrette. "Group B Streptococcal Disease in the Newborn: Maternal GBS-Screening Methods and Antimicrobial Prophylaxis." *European Obstetrics and Gynaecology-Touch Briefings* 3 (2008): 58-62.
5. Public Health England. 2018. 'Detection of Carriage of Group B Streptococci (Streptococcus Agalactiae).' UK Standards for Microbiology Investigations. Vol. B 58.
6. Parks, Tom, Lucinda Barrett, and Nicola Jones. 2015. 'Invasive Streptococcal Disease: A Review for Clinicians.' *British Medical Bulletin* 115 (1): 77-89. <https://doi.org/10.1093/bmb/ldv027>.
7. Agudelo Higuaita, Nelson I., and Mark M Huycke. 2014. 'Enterococcal Disease, Epidemiology, and Implications for Treatment.' In *Enterococci: From Commensals to Leading Causes of Drug Resistant Infection* [Internet], edited by MS Gilmore and DB Clewell, 1-35. Boston: Massachusetts Eye and Ear Infirmary; 2014.
8. Levitus, Matthew, Ayesan Rewane, and Thomas B. Perera. 2020. 'Vancomycin-Resistant Enterococci.' In *StatPearls* [Internet]. Treasure Island (FL): StatPearls Publishing.
9. Said, Mina S., Ekta Tirhani, and Emil Lesho. 2021. 'Enterococcus Infections.' In *StatPearls* [Internet]. Treasure Island (FL): StatPearls Publishing.
10. Public Health England. 2014. 'Identification of Streptococcus Species, Enterococcus Species and Morphologically Similar Organisms'. UK Standards for Microbiology Investigations. ID 4 (3). <https://www.gov.uk/government/publications/smi-id-4-identification-of-streptococcus-species->

11.

### Symbol Legend

Symbol	Definition
	Catalogue number
	In Vitro Diagnostic Medical Device
	Batch code
	Temperature limit
	Use-by date
	Keep away from sunlight
	Do not re-use
	Consult instructions for use or consult electronic instructions for use
	Contains sufficient for <n> tests
	Do not use if packaging damaged and consult instructions for use
	USA: Caution: Federal law restricts this device to sale by or on order of a Physician
	Manufacturer
	Authorized representative in the European Community/ European Union
	European Conformity Assessment
	UK Conformity Assessment
	Unique device identifier

© 2022 Thermo Fisher Scientific Inc. All rights reserved. ATCC® is a trademark of ATCC. All other trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries. This information is not intended to encourage use of these products in any manner that might infringe the intellectual property rights of others.

Oxoid GmbH, Am Lippeglacis 4-8, 46483 Wesel, Germany

For technical assistance please contact your local distributor.

#### Revision Information

Version	Date of issue and modifications introduced
1.0	2022-06-30. New document

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### Literatura

1. Foster, Timothy J. 1996. "Staphylococcus." In Medical Microbiology. 4th Edition., 4th ed. Galveston.
2. Public Health England. 2020b. 'Investigation of Specimens for Screening for MRSA'. UK Standards for Microbiology Investigations. B 29 (7). <https://www.gov.uk/government/publications/smi-b-29-investigation-of-specimens-for-screening-for-mrsa>.
3. National Institute for Health and Care Excellence. 2016. 'Healthcare-Associated Infections'. QS113. [www.nice.org.uk/guidance/qs113](http://www.nice.org.uk/guidance/qs113).
4. Melin, Pierrette. "Group B Streptococcal Disease in the Newborn: Maternal GBS-Screening Methods and Antimicrobial Prophylaxis." European Obstetrics and Gynaecology-Touch Briefings 3 (2008): 58-62.
5. Public Health England. 2018. 'Detection of Carriage of Group B Streptococci (Streptococcus Agalactiae)'. UK Standards for Microbiology Investigations. Vol. B 58.
6. Parks, Tom, Lucinda Barrett, and Nicola Jones. 2015. 'Invasive Streptococcal Disease: A Review for Clinicians.' British Medical Bulletin 115 (1): 77-89. <https://doi.org/10.1093/bmb/ldv027>.
7. Agudelo Higuera, Nelson I., and Mark M Huycke. 2014. 'Enterococcal Disease, Epidemiology, and Implications for Treatment.' In Enterococci: From Commensals to Leading Causes of Drug Resistant Infection [Internet], edited by MS Gilmore and DB Clewell, 1-35. Boston: Massachusetts Eye and Ear Infirmary; 2014.
8. Levitus, Matthew, Ayesan Rewane, and Thomas B. Perera. 2020. 'Vancomycin-Resistant Enterococci.' In StatPearls [Internet]. Treasure Island (FL): StatPearls Publishing.
9. Said, Mina S., Ekta Tirthani, and Emil Lesho. 2021. 'Enterococcus Infections.' In StatPearls [Internet]. Treasure Island (FL): StatPearls Publishing.
10. Public Health England. 2014a. 'Identification of Streptococcus Species, Enterococcus Species and Morphologically Similar Organisms'. UK Standards for Microbiology Investigations. ID 4 (3). <https://www.gov.uk/government/publications/smi-id-4-identification-of-streptococcus-species-enterococcus-species-and-morphologically-similar-organisms>.

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## Bibliographie

1. Foster, Timothy J. 1996. "Staphylococcus." In Medical Microbiology. 4th Edition., 4th ed. Galveston.
2. Public Health England. 2020b. 'Investigation of Specimens for Screening for MRSA'. UK Standards for Microbiology Investigations. B 29 (7). <https://www.gov.uk/government/publications/smi-b-29-investigation-of-specimens-for-screening-for-mrsa>.
3. National Institute for Health and Care Excellence. 2016. 'Healthcare-Associated Infections'. QS113. [www.nice.org.uk/guidance/qs113](http://www.nice.org.uk/guidance/qs113).
4. Melin, Pierrette. "Group B Streptococcal Disease in the Newborn: Maternal GBS-Screening Methods and Antimicrobial Prophylaxis." European Obstetrics and Gynaecology-Touch Briefings 3 (2008): 58-62.
5. Public Health England. 2018. 'Detection of Carriage of Group B Streptococci (Streptococcus Agalactiae)'. UK Standards for Microbiology Investigations. Vol. B 58.
6. Parks, Tom, Lucinda Barrett, and Nicola Jones. 2015. 'Invasive Streptococcal Disease: A Review for Clinicians.' British Medical Bulletin 115 (1): 77–89. <https://doi.org/10.1093/bmb/ldv027>.
7. Agudelo Higuera, Nelson I., and Mark M Huycke. 2014. 'Enterococcal Disease, Epidemiology, and Implications for Treatment.' In Enterococci: From Commensals to Leading Causes of Drug Resistant Infection [Internet], edited by MS Gilmore and DB Clewell, 1–35. Boston: Massachusetts Eye and Ear Infirmary; 2014.
8. Levitus, Matthew, Ayesan Rewane, and Thomas B. Perera. 2020. 'Vancomycin-Resistant Enterococci.' In StatPearls [Internet]. Treasure Island (FL): StatPearls Publishing.
9. Said, Mina S., Ekta Tirthani, and Emil Lesho. 2021. 'Enterococcus Infections.' In StatPearls [Internet]. Treasure Island (FL): StatPearls Publishing.
10. Public Health England. 2014a. 'Identification of Streptococcus Species, Enterococcus Species and Morphologically Similar Organisms'. UK Standards for Microbiology Investigations. ID 4 (3). <https://www.gov.uk/government/publications/smi-id-4-identification-of-streptococcus-species-enterococcus-species-and-morphologically-similar-organisms>.

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1. Foster, Timothy J. 1996. "Staphylococcus." In Medical Microbiology. 4th Edition., 4th ed. Galveston.
2. Public Health England. 2020b. 'Investigation of Specimens for Screening for MRSA'. UK Standards for Microbiology Investigations. B 29 (7). <https://www.gov.uk/government/publications/smi-b-29-investigation-of-specimens-for-screening-for-mrsa>.
3. National Institute for Health and Care Excellence. 2016. 'Healthcare-Associated Infections'. QS113. [www.nice.org.uk/guidance/qs113](http://www.nice.org.uk/guidance/qs113).
4. Melin, Pierrette. "Group B Streptococcal Disease in the Newborn: Maternal GBS-Screening Methods and Antimicrobial Prophylaxis." *European Obstetrics and Gynaecology-Touch Briefings* 3 (2008): 58-62.
5. Public Health England. 2018. 'Detection of Carriage of Group B Streptococci (Streptococcus Agalactiae)'. UK Standards for Microbiology Investigations. Vol. B 58.
6. Parks, Tom, Lucinda Barrett, and Nicola Jones. 2015. 'Invasive Streptococcal Disease: A Review for Clinicians.' *British Medical Bulletin* 115 (1): 77-89. <https://doi.org/10.1093/bmb/ldv027>.
7. Agudelo Higuera, Nelson I., and Mark M Huycke. 2014. 'Enterococcal Disease, Epidemiology, and Implications for Treatment.' In *Enterococci: From Commensals to Leading Causes of Drug Resistant Infection* [Internet], edited by MS Gilmore and DB Clewell, 1-35. Boston: Massachusetts Eye and Ear Infirmary; 2014.
8. Levitus, Matthew, Ayesan Rewane, and Thomas B. Perera. 2020. 'Vancomycin-Resistant Enterococci.' In *StatPearls* [Internet]. Treasure Island (FL): StatPearls Publishing.
9. Said, Mina S., Ekta Tirthani, and Emil Lesho. 2021. 'Enterococcus Infections.' In *StatPearls* [Internet]. Treasure Island (FL): StatPearls Publishing.
10. Public Health England. 2014a. 'Identification of Streptococcus Species, Enterococcus Species and Morphologically Similar Organisms'. UK Standards for Microbiology Investigations. ID 4 (3). <https://www.gov.uk/government/publications/smi-id-4-identification-of-streptococcus-species-enterococcus-species-and-morphologically-similar-organisms>.

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### Bibliografia

1. Foster, Timothy J. 1996. "Staphylococcus." In Medical Microbiology. 4th Edition., 4th ed. Galveston.
2. Public Health England. 2020b. 'Investigation of Specimens for Screening for MRSA'. UK Standards for Microbiology Investigations. B 29 (7). <https://www.gov.uk/government/publications/smi-b-29-investigation-of-specimens-for-screening-for-mrsa>.
3. National Institute for Health and Care Excellence. 2016. 'Healthcare-Associated Infections'. QS113. [www.nice.org.uk/guidance/qs113](http://www.nice.org.uk/guidance/qs113).
4. Melin, Pierrette. "Group B Streptococcal Disease in the Newborn: Maternal GBS-Screening Methods and Antimicrobial Prophylaxis." *European Obstetrics and Gynaecology-Touch Briefings* 3 (2008): 58-62.
5. Public Health England. 2018. 'Detection of Carriage of Group B Streptococci (Streptococcus Agalactiae).'UK Standards for Microbiology Investigations. Vol. B 58.

6. Parks, Tom, Lucinda Barrett, and Nicola Jones. 2015. 'Invasive Streptococcal Disease: A Review for Clinicians.' *British Medical Bulletin* 115 (1): 77–89. <https://doi.org/10.1093/bmb/ldv027>.
7. Agudelo Higuera, Nelson I., and Mark M Huycke. 2014. 'Enterococcal Disease, Epidemiology, and Implications for Treatment.' In *Enterococci: From Commensals to Leading Causes of Drug Resistant Infection* [Internet], edited by MS Gilmore and DB Clewell, 1–35. Boston: Massachusetts Eye and Ear Infirmary; 2014.
8. Levitus, Matthew, Ayesan Rewane, and Thomas B. Perera. 2020. 'Vancomycin-Resistant Enterococci.' In *StatPearls* [Internet]. Treasure Island (FL): StatPearls Publishing.
9. Said, Mina S., Ekta Tirthani, and Emil Lesho. 2021. 'Enterococcus Infections.' In *StatPearls* [Internet]. Treasure Island (FL): StatPearls Publishing.
10. Public Health England. 2014a. 'Identification of Streptococcus Species, Enterococcus Species and Morphologically Similar Organisms'. UK Standards for Microbiology Investigations. ID 4 (3). <https://www.gov.uk/government/publications/smi-id-4-identification-of-streptococcus-species-enterococcus-species-and-morphologically-similar-organisms>.

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Oxoid GmbH, Am Lippegelacis 4-8, 46483