Instructions for Use: *Brilliance*™ ESBL Agar

[REF] PO5302A

Intended Use

Brilliance™ ESBL Agar is a qualitative selective chromogenic medium screening for the presence of extended-spectrum beta-lactamase (ESBL) producing bacteria in faecal and rectal screening samples and for identifying clinical isolates. Used in a diagnostic workflow to aid clinicians in determining potential treatment options for patients suspected of having bacterial infections.

The device is for professional use only and is not intended to be used as a companion diagnostic device.

Summary and Explanation

In the early 1980s, the third generation cephalosporins were introduced in response to the increasing number of bacteria that were resistant to antibiotics due to the production of β -lactamases.¹ The new cephalosporins were resistant to most found lactamases such as SHV-1 and TEM-1 in *Escherichia coli* and *Klebsiella pneumoniae* and did not suffer from the nephrotoxic effects associated with aminoglycosides and polymyxins. Unfortunately, almost as soon as they were introduced, there were reports of a plasmid encoded β -lactamase (SHV-2) isolated from *Klebsiella ozaenae* that was capable of hydrolysing these new antibiotics.² The discovery of others, which were related to TEM-1 and TEM-2, soon followed.³ Because of their increased spectrum of activity, especially against the oxyimino-cephalosporins, these enzymes were called ESBL and were defined as β -lactamases that were capable of hydrolysing penicillin's first, second, and third generation cephalosporins.⁴ Because ESBLs are carried primarily on plasmids, they are easily transmitted among different bacterial species, and ESBLs are now the most common resistance mechanism of Gram-negative bacteria against β -lactam antibiotics.⁵ The spread of ESBL-producing bacteria has increased dramatically over the last two decades.^{6,7,8}

Principle of Method

Differentiation of ESBL-producing *E. coli*, members of the KESC group, and *Proteus*, *Morganella* and *Providencia* species is achieved through the inclusion of two chromogens that are targeted by specific enzymes: β -galactosidase and β -glucuronidase. The action of these enzymes on the chromogens causes release of the coloured component inside the bacterial cell resulting in coloured colonies. The colour produced depends on which enzymes the organisms produce. The presence of β -galactosidase enzymes in *E. coli* results in blue colonies, and rarely the presence of β -glucuronidase in some isolates results in pink colonies. The expression of β -galactosidase by KESC species produces green colonies. *Proteus*, *Morganella* and *Providencia* species do not utilise either chromogen but are able to deaminate tryptophan, which results in tan colonies with a brown halo. Cefpodoxime inhibits the growth of non ESBL-producing species and benzo[b]thiophene-2-boronic acid suppresses the growth of AmpC-producing bacteria facilitating a high degree of sensitivity for the detection of ESBL-producing bacteria. Additional antibiotics and antifungals in combination with sodium lauryl sulphate suppress the growth of other competing flora.

Typical Formula

	grams per litre
Peptone	12
Sodium chloride	5
Phosphate buffers	4
Agar	15
Chromogenic mix	4
Antibiotic mix	0.28

Physical Appearance

Colour	Oyster white
Clarity	Opaque
Fill weight	17g \pm 5%
pH	6.9 \pm 0.2

Materials Provided

- The pack contains 10 x 90mm agar plates, film wrapped.
- Each plate should only be used once.
- Each pack contains enough plates for 10 individual tests.

Materials Required but Not Supplied

- Inoculating loops
- Swabs
- Collection containers
- Incubators
- Quality control organisms

Storage

- Store product in its original packaging at 2–12°C until used.
- The product may be used until the expiry date stated on the label.
- Store away from light.
- Allow product to equilibrate to room temperature before use.
- Do not incubate prior to use.

Warnings and Precautions

- For *in vitro* diagnostic use only.
- For professional use only.
- Inspect the product packaging before first use
- Do not use the product if there is any visible damage to the packaging or plates.
- Do not use the product beyond the stated expiry date.
- Do not use the device if signs of contamination are present.
- Do not use the device if the colour has changed or there are other signs of deterioration.
- It is the responsibility of each laboratory to manage waste produced according to their nature and degree of hazard and to have them treated or disposed of in accordance with any federal, state and local applicable regulations. Directions should be read and followed carefully. This includes the disposal of used or unused reagents as well as any other contaminated disposable material following procedures for infectious or potentially infectious products.

Refer to the Material Safety Data Sheet (MSDS) for safe handling and disposal of the product at www.thermofisher.com

Materials of animal origin

Brilliance ESBL Agar contains yeast extract manufactured from microbial raw materials, and peptone manufactured from porcine, bovine and casein raw materials.

Specimen Collection, Handling and Storage

Requirements for specimen collection, handling and storage of samples are described in local procedures and guidelines, such as the UK Standards for Microbiology Investigations (UK SMI) B 59 (Public Health England, 2016).

Procedure

- Allow product to equilibrate to room temperature before use.
- Inoculate and streak the specimen onto the medium using a standard loop. *Brilliance* ESBL Agar can be inoculated direct from rectal screening swabs, faecal samples or from isolated colonies prepared as a liquid suspension approximately equivalent to 0.5 McFarland turbidity, according to local guidelines.
- Incubate plates aerobically for 18–24 hours at 35±2 °C.
- Visually inspect plates to assess colony growth and colour under good lighting.
- Negative plates should be incubated for an additional 24 hours and re-assessed.

Interpretation

The presence of blue, pink, green, colourless, or brown colonies indicate the sample is ESBL positive:

- Blue colonies indicate *E. coli*. In rare cases *E. coli* colonies may also be pink.
- Green colonies indicate *Klebsiella*, *Enterobacter*, *Serratia* and *Citrobacter* (KESC).
- Colourless colonies indicate *Salmonella*, *Acinetobacter* or other.
- Brown colonies with a halo indicate *Proteus*, *Morganella*, *Providencia*.

Quality Control

It is the responsibility of the user to perform Quality Control testing taking into account the intended use of the medium, and in accordance with any local applicable regulations (frequency, number of strains, incubation temperature etc.).

The performance of this medium can be verified by testing the following reference strains.

Incubation Conditions: 18 - 24 h @ 35° ± 2°C aerobic.

Positive Controls	
Colony count is ≥ 50% of the control medium count	
<i>Klebsiella pneumoniae</i> ATCC® 700603	1 - 2 mm, green colonies.
<i>Escherichia coli</i> TEM-3 NCTC 13351	1 - 2 mm, blue/turquoise colonies.
Negative Controls	
<i>Enterobacter cloacae</i> ATCC® 23355	No Growth
<i>Citrobacter freundii</i> NCTC 8581	No Growth
<i>Candida albicans</i> ATCC® 10231	No Growth to Inhibited Growth

Analytical Performance

A study of 58 ESBL pure cultures and 66 ESBL negative cultures (including non-ESBL *Enterobacteriaceae*) was undertaken,⁹ along with 10 un-spiked and spiked faecal samples. 0.5 McFarland standard suspension of each isolate was prepared and streaked on each of the test media. For positive cultures, 0.5 McFarland standard suspension was further serially diluted to 10⁻⁵, and 50 µl of the 10⁻³, 10⁻⁴ and 10⁻⁵ dilutions of each positive test organism was spread plated onto each of the test media. Plates were incubated according to the manufacturer's instructions for up to 48 hours. Plates were read by staff not involved in the development of the project. Sensitivity for *Brilliance* ESBL Agar was calculated based on the presence of correctly coloured blue, pink, green, or brown colonies (presumptive ESBL) and any other coloured or colourless colonies were reported. Specificity was calculated based on the number of true negative plates, i.e., the number of plates with correctly coloured non-ESBL colonies plus plates with no growth.

Performance	Incubation time (hr.)	<i>Brilliance</i> ESBL Agar (%)
Sensitivity	24	94.64
Specificity	24	94.03

Clinical Performance

Brilliance ESBL Agar has been evaluated through a series of external trials conducted at different European laboratories and hospitals, which compared and demonstrated the performance of the device in a clinical setting. The performance characteristics sensitivity, specificity, positive predictive value (PPV) and negative predictive value (NPV) have been assessed for this medium.

Brilliance ESBL Agar proved to be a consistently highly selective medium for the isolation of ESBL-producing organisms from clinical samples, detecting presumptive ESBL-producers within 24 hours. Additionally, it showed noteworthy inhibition of non-target organisms by yielding an excellent NPV across multiple studies, enabling the rapid exclusion of patients not carrying ESBL-producing bacteria.

In one performance study (Trial 1), a total of 528 clinical samples (including stool samples, rectal swabs) taken from patients at UCL Mont Godinne University Hospital, Yvoir, Belgium, were evaluated on *Brilliance* ESBL Agar.^{9, 10} A total of 200 ESBL confirmed isolates found in clinical samples taken from patients was included in this study. Each clinical specimen was homogenized in 1 mL sterile 0.85% saline and vortexed. A 50 µL sample of this suspension was inoculated onto the media. After incubation at 35°C for 24 hours, growth and colony morphology on the plates was inspected, according to an interpretation guide. Selection criteria for work-out (oxidase negative and coloured colonies) was applied.

In another performance study (Trial 2)⁹, a total of 504 rectal screening swab samples taken from patients at two hospitals in Belgium - AZ Sint Lucas and AZ Gezondheidszorg, were evaluated. 100 µL of the transport medium was inoculated onto *Brilliance* ESBL Agar. After incubation, growth and colony morphology on the plates was inspected, according to manufacturer's guidelines.

In another performance study (Trial 3)⁹, a total of 84 samples from clinical samples (including stool samples) were plated onto the media, and were evaluated as part of a randomised, investigator-blinded evolution of culture-based approaches study, which was conducted at

Department of Medical Microbiology, Vaccine & Infectious Disease Institute, Hasselt University, Belgium. After incubation, growth and colony morphology on the plates was inspected, according to manufacturer's guidelines.

Performance of *Brilliance* ESBL Agar.

Performance Characteristic	Brilliance ESBL Agar - Pre-Launch trials	Brilliance ESBL Agar – Post-Launch trials	
	Trial - 1 (Belgium)	Trial - 2 (Belgium)	Trial - 3 (Belgium)
Sensitivity (%)	94.9	88.0	99.4 ^a
Specificity (%)	95.7	87.0	99.2 ^a
PPV (%)	73.7	ND	ND
NPV (%)	99.3	ND	ND

ND – Not done.

^a – Mean value after 24-hour incubation.

Summary of results found in the studies evaluated in a literature review.

Study	Huang et al., 2010	Willems et al., 2013	Ongut et al., 2014	Blane et al., 2016
Sensitivity (%)	94.9	87.50	97.9	59
Specificity (%)	95.7	82.11	100	87
PPV (%)	73.7	38.89	100	ND
NPV (%)	99.3	98.06	96.9	ND

ND – Not done.

The Summary of Safety and Performance (SSP) for this device will be available in the European database on medical devices where it is linked to the device's Basic UDI-DI (5032384BrillianceESBL2V).

Refer to: Eudamed, <https://ec.europa.eu/tools/eudamed>

Limitations

Organisms with atypical enzyme patterns may give anomalous reactions on *Brilliance* ESBL Agar. Organisms with atypical MICs to the antibiotics contained in the medium may give anomalous reactions. Samples containing faecal material or blood may cause some localised discoloration within the medium. This discoloration should not be confused with a true chromogenic reaction, where coloured colonies are visible. In line with EUCAST guidelines, it is not recommended that antimicrobial susceptibility testing is carried out on colonies taken directly from *Brilliance* ESBL Agar; selective media can alter the MIC of an antimicrobial agent. Identifications are presumptive and should be confirmed.

Serious Incidents

Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the relevant regulatory authority in which the user and/or the patient is established

Bibliography

- Farber B, Moellering RC Jr. 1982. 'The third generation cephalosporins.' *Bull N Y Acad Med*. 58: 696–710. PMID: 6762896.
- Kliebe, C., B. A. Nies, J. F. Meyer, R. M. Tolxdorff-Neutzing, and B. Wiedemann. 1985. 'Evolution of Plasmid-Coded Resistance to Broad-Spectrum Cephalosporins'. *Antimicrobial Agents and*

- Brun-Buisson, C, A Philippon, M Ansquer, P Legrand, F Montravers, and J Duval. 1987. 'Transferable enzymatic resistance to third-generation cephalosporins during nosocomial outbreak of multiresistant *Klebsiella pneumoniae*.' *The Lancet* 330, no. 8554 302-306.
- Paterson, D L., and R A. Bonomo. 2005. 'Extended-spectrum β -lactamases: a clinical update.' *Clinical Microbiology Reviews* 18, no. 4 657-686.
- Al-Bayssari, C., Dabboussi, F., Hamze, M., and Rolain, J.-M. 2015. 'Detection of expanded-spectrum β -lactamases in Gram-negative bacteria in the 21st century.' *Expert Rev. Anti Infect. Ther.* 13, 1139–1158. <https://doi.org/10.1586/14787210.2015.1066247>.
- Kawamura, Kumiko, Noriyuki Nagano, Masahiro Suzuki, Jun-ichi Wachino, Kouji Kimura, and Yoshichika Arakawa. 2017. 'ESBL-Producing *Escherichia Coli* and Its Rapid Rise among Healthy People'. *Food Safety (Tokyo, Japan)* 5 (4): 122–50. <https://doi.org/10.14252/foodsafetyfscj.2017011>.
- Hu, Yanhong Jessika, Anju Ogyu, Benjamin J. Cowling, Keiji Fukuda, and Herbert H. Pang. 2019. 'Available Evidence of Antibiotic Resistance from Extended-Spectrum β -Lactamase-Producing Enterobacteriaceae in Paediatric Patients in 20 Countries: A Systematic Review and Meta-Analysis'. *Bulletin of the World Health Organization* 97 (7): 486-501B. <https://doi.org/10.2471/BLT.18.225698>.
- PHE. UK Standards for Microbiology Investigations: Detection of Enterobacteriaceae producing extended spectrum β -lactamases. Standards Unit, Microbiology Services, PHE Bacteriology. B 59. Issue no: 4.1, Issue date: 17.08.16.
- Data on file.
- Huang, Te-Din, Pierre Bogaerts, Catherine Berhin, Amelie Guisset, and Youri Glupczynski. 2010. 'Evaluation of *Brilliance*™ ESBL Agar, a Novel Chromogenic Medium for Detection of Extended-Spectrum-Beta- Lactamase-Producing Enterobacteriaceae'. *Journal of Clinical Microbiology* 48 (6): 2091–96. <https://doi.org/10.1128/JCM.02342-09>.

Glossary of symbols

Symbol/Label	Meaning
	Manufacturer
	In Vitro Diagnostic Medical Device
	Temperature limit
	Batch Code
	Catalog Number
	Do not re-use
	Consult instructions for use or consult electronic instructions for use

	Contains sufficient for <n> tests
	Use-by date
	Do not use if package is damaged and Consult instructions for use
	Authorized representative in the European Community/ European Union
	Unique device identifier
	USA: Caution: Federal law restricts this device to sale by or on order of a Physician
	European Conformity Mark
	UK Conformity Mark

The ATCC Licensed Derivative® Emblem, the ATCC Licensed Derivative® word mark, and the ATCC catalogue marks are trademarks of ATCC. Thermo Fisher Scientific Inc. is licensed to use these trademarks and to sell products derived from ATCC®

© 2021 Thermo Fisher Scientific Inc. All rights reserved. ATCC® is a trademark of ATCC. All other trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries. This information is not intended to encourage use of these products in any manner that might infringe the intellectual property rights of others.

Oxoid Limited, Wade Road, Basingstoke, Hampshire, RG24 8PW, England

2797
For technical assistance please contact your local distributor.

Version	Modifications introduced
1.0	2022-06-13 New document

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- Brun-Buisson, C., A. Philippon, M. Ansquer, P. Legrand, F. Montravers, and J. Duval. 1987. 'Transferable enzymatic resistance to third-generation cephalosporins during nosocomial outbreak of multiresistant *Klebsiella pneumoniae*.' *The Lancet* 330, no. 8554: 302–306.
- Paterson, D L., and R A. Bonomo. 2005. 'Extended-spectrum β -lactamases: a clinical update.' *Clinical Microbiology Reviews* 18, no. 4: 657–686.
- Al-Bayssari, C., Dabboussi, F., Hamze, M., and Rolain, J.-M. 2015. 'Detection of expanded-spectrum β -lactamases in Gram-negative bacteria in the 21st century.' *Expert Rev. Anti Infect. Ther.* 13, 1139–1158. <https://doi.org/10.1586/14787210.2015.1066247>.
- Kawamura, Kumiko, Noriyuki Nagano, Masahiro Suzuki, Jun-ichi Wachino, Kouji Kimura, and Yoshichika Arakawa. 2017. 'ESBL-Producing *Escherichia Coli* and Its Rapid Rise among Healthy People'. *Food Safety* (Tokyo, Japan) 5 (4): 122–50. <https://doi.org/10.14252/foodsafetyfscj.2017011>.
- Hu, Yanhong Jessika, Anju Ogyu, Benjamin J. Cowling, Keiji Fukuda, and Herbert H. Pang. 2019. 'Available Evidence of Antibiotic Resistance from Extended-Spectrum β -Lactamase-Producing Enterobacteriaceae in Paediatric Patients in 20 Countries: A Systematic Review and Meta-Analysis'. *Bulletin of the World Health Organization* 97 (7): 486–501B. <https://doi.org/10.2471/BLT.18.225698>.
- PHE. UK Standards for Microbiology Investigations: Detection of Enterobacteriaceae producing extended spectrum β -lactamases. Standards Unit, Microbiology Services, PHE Bacteriology. B 59. Issue no: 4.1, Issue date: 17.08.16.
- Data on file.
- Huang, Te-Din, Pierre Bogaerts, Catherine Berhin, Amelie Guisset, and Youri Glupczynski. 2010. 'Evaluation of Brilliance™ ESBL Agar, a Novel Chromogenic Medium for Detection of Extended-Spectrum-Beta- Lactamase-Producing Enterobacteriaceae'. *Journal of Clinical Microbiology* 48 (6): 2091–96. <https://doi.org/10.1128/JCM.02342-09>.

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- Paterson, D L., and R A. Bonomo. 2005. 'Extended-spectrum β -lactamases: a clinical update.' *Clinical Microbiology Reviews* 18, no. 4 657-686.
- Al-Bayssari, C., Dabboussi, F., Hamze, M., and Rolain, J.-M. 2015. 'Detection of expanded-spectrum β -lactamases in Gram-negative bacteria in the 21st century.' *Expert Rev. Anti Infect. Ther.* 13, 1139–1158. <https://doi.org/10.1586/14787210.2015.1066247>.
- Kawamura, Kumiko, Noriyuki Nagano, Masahiro Suzuki, Jun-ichi Wachino, Kouji Kimura, and Yoshichika Arakawa. 2017. 'ESBL-Producing *Escherichia Coli* and Its Rapid Rise among Healthy People'. *Food Safety* (Tokyo, Japan) 5 (4): 122–50. <https://doi.org/10.14252/foodsafetyfsci.2017011>.
- Hu, Yanhong Jessika, Anju Ogyu, Benjamin J. Cowling, Keiji Fukuda, and Herbert H. Pang. 2019. 'Available Evidence of Antibiotic Resistance from Extended-Spectrum β -Lactamase-Producing Enterobacteriaceae in Paediatric Patients in 20 Countries: A Systematic Review and Meta-Analysis'. *Bulletin of the World Health Organization* 97 (7): 486-501B. <https://doi.org/10.2471/BLT.18.225698>.
- PHE. UK Standards for Microbiology Investigations: Detection of Enterobacteriaceae producing extended spectrum β -lactamases. Standards Unit, Microbiology Services, PHE Bacteriology. B 59. Issue no: 4.1, Issue date: 17.08.16.
- Data on file.
- Huang, Te-Din, Pierre Bogaerts, Catherine Berhin, Amelie Guisset, and Youri Glupczynski. 2010. 'Evaluation of *Brilliance*™ ESBL Agar, a Novel Chromogenic Medium for Detection of Extended-Spectrum-Beta-Lactamase-Producing Enterobacteriaceae'. *Journal of Clinical Microbiology* 48 (6): 2091–96. <https://doi.org/10.1128/JCM.02342-09>.

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- Brun-Buisson, C., A Philippon, M Ansquer, P Legrand, F Montravers, and J Duval. 1987. 'Transferable enzymatic resistance to third-generation cephalosporins during nosocomial outbreak of multiresistant Klebsiella pneumoniae.' The Lancet 330, no. 8554 302-306.
- Paterson, D L., and R A. Bonomo. 2005. 'Extended-spectrum β -lactamases: a clinical update.' Clinical Microbiology Reviews 18, no. 4 657-686.
- Al-Bayssari, C., Dabboussi, F., Hamze, M., and Rolain, J.-M. 2015. 'Detection of expanded-spectrum β -lactamases in Gram-negative bacteria in the 21st century.' Expert Rev. Anti Infect. Ther. 13, 1139–1158. <https://doi.org/10.1586/14787210.2015.1066247>.
- Kawamura, Kumiko, Noriyuki Nagano, Masahiro Suzuki, Jun-ichi Wachino, Kouji Kimura, and Yoshichika Arakawa. 2017. 'ESBL-Producing Escherichia Coli and Its Rapid Rise among Healthy People'. Food Safety (Tokyo, Japan) 5 (4): 122–50. <https://doi.org/10.14252/foodsafetyfscj.2017011>.
- Hu, Yanhong Jessika, Anju Ogyu, Benjamin J. Cowling, Keiji Fukuda, and Herbert H. Pang. 2019. 'Available Evidence of Antibiotic Resistance from Extended-Spectrum β -Lactamase-Producing Enterobacteriaceae in Paediatric Patients in 20 Countries: A Systematic Review and Meta-Analysis'. Bulletin of the World Health Organization 97 (7): 486-501B. <https://doi.org/10.2471/BLT.18.225698>.
- PHE. UK Standards for Microbiology Investigations: Detection of Enterobacteriaceae producing extended spectrum β -lactamases. Standards Unit, Microbiology Services, PHE Bacteriology. B 59. Issue no: 4.1, Issue date: 17.08.16.
- Data on file.
- Huang, Te-Din, Pierre Bogaerts, Catherine Berhin, Amelie Guisset, and Youri Glupczynski. 2010. 'Evaluation of Brilliance™ ESBL Agar, a Novel Chromogenic Medium for Detection of Extended-Spectrum-Beta- Lactamase-Producing Enterobacteriaceae'. Journal of Clinical Microbiology 48 (6): 2091–96. <https://doi.org/10.1128/JCM.02342-09>.

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
Bibliographie

- Farber B, Moellering RC Jr. 1982. 'The third generation cephalosporins.' *Bull N Y Acad Med.* 58: 696–710. PMID: 6762896.
- Kliebe, C., B. A. Nies, J. F. Meyer, R. M. Tolxdorff-Neutzing, and B. Wiedemann. 1985. 'Evolution of Plasmid-Coded Resistance to Broad-Spectrum Cephalosporins'. *Antimicrobial Agents and Chemotherapy* 28 (2): 302–7. <https://doi.org/10.1128/aac.28.2.302>.
- Brun-Buisson, C, A Philippon, M Ansquer, P Legrand, F Montravers, and J Duval. 1987. 'Transferable enzymatic resistance to third-generation cephalosporins during nosocomial outbreak of multiresistant *Klebsiella pneumoniae*.' *The Lancet* 330, no. 8554 302-306.
- Paterson, D L., and R A. Bonomo. 2005. 'Extended-spectrum β-lactamases: a clinical update.' *Clinical Microbiology Reviews* 18, no. 4 657-686.
- Al-Bayssari, C., Dabboussi, F., Hamze, M., and Rolain, J.-M. 2015. 'Detection of expanded-spectrum β-lactamases in Gram-negative bacteria in the 21st century.' *Expert Rev. Anti Infect. Ther.* 13, 1139–1158. <https://doi.org/10.1586/14787210.2015.1066247>.
- Kawamura, Kumiko, Noriyuki Nagano, Masahiro Suzuki, Jun-ichi Wachino, Kouji Kimura, and Yoshichika Arakawa. 2017. 'ESBL-Producing *Escherichia Coli* and Its Rapid Rise among Healthy People'. *Food Safety (Tokyo, Japan)* 5 (4): 122–50. <https://doi.org/10.14252/foodsafetyfscj.2017011>.
- Hu, Yanhong Jessika, Anju Ogyu, Benjamin J. Cowling, Keiji Fukuda, and Herbert H. Pang. 2019. 'Available Evidence of Antibiotic Resistance from Extended-Spectrum β-Lactamase-Producing Enterobacteriaceae in Paediatric Patients in 20 Countries: A Systematic Review and Meta-Analysis'. *Bulletin of the World Health Organization* 97 (7): 486-501B. <https://doi.org/10.2471/BLT.18.225698>.
- PHE. UK Standards for Microbiology Investigations: Detection of Enterobacteriaceae producing extended spectrum β-lactamases. Standards Unit, Microbiology Services, PHE Bacteriology. B 59. Issue no: 4.1, Issue date: 17.08.16.
- Data on file.
- Huang, Te-Din, Pierre Bogaerts, Catherine Berhin, Amelie Guisset, and Youri Glupczynski. 2010. 'Evaluation of *Brilliance*™ ESB� Agar, a Novel Chromogenic Medium for Detection of Extended-Spectrum-Beta-Lactamase-Producing Enterobacteriaceae'. *Journal of Clinical Microbiology* 48 (6): 2091–96. <https://doi.org/10.1128/JCM.02342-09>.

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6. Kawamura, Kumiko, Noriyuki Nagano, Masahiro Suzuki, Jun-Ichi Wachino, Kouji Kimura, and Yoshichika Arakawa. 2017. 'ESBL-Producing Escherichia Coli and Its Rapid Rise among Healthy People'. *Food Safety* (Tokio, Japan) 5 (4): 122-50. <https://doi.org/10.14252/foodsafetyfscj.2017011>.
7. Hu, Yanhong Jessica, Anju Ogyu, Benjamin J. Cowling, Keiji Fukuda, and Herbert H. Pang. 2019. 'Available Evidence of Antibiotic Resistance from Extended-Spectrum β -Lactamase-Producing Enterobacteriaceae in Paediatric Patients in 20 Countries: A Systematic Review and Meta-Analysis'. *Bulletin of the World Health Organization* 97 (7): 486-501B. <https://doi.org/10.2471/BLT.18.225698>.
8. PHE. UK Standards for Microbiology Investigations: Detection of Enterobacteriaceae producing extended spectrum β -lactamases. Standards Unit, Microbiology Services, PHE Bacteriology. B 59. Ausgabe Nr.: 4.1, Ausgabedatum: 17.08.16.
9. Daten in den Akten.
10. Huang, Te-Din, Pierre Bogaerts, Catherine Berhin, Amelie Guisset, and Youri Glupczynski. 2010. 'Evaluation of *Brilliance*™ ESBL Agar, a Novel Chromogenic Medium for Detection of Extended-Spectrum-Beta- Lactamase-Producing Enterobacteriaceae'. *Journal of Clinical Microbiology* 48 (6): 2091–96. <https://doi.org/10.1128/JCM.02342-09>.

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- Farber B, Moellering RC Jr. 1982. 'The third generation cephalosporins.' *Bull N Y Acad Med*. 58: 696–710. PMID: 6762896.
- Kliebe, C., B. A. Nies, J. F. Meyer, R. M. Tolxdorff-Neutzing, and B. Wiedemann. 1985. 'Evolution of Plasmid-Coded Resistance to Broad-Spectrum Cephalosporins'. *Antimicrobial Agents and Chemotherapy* 28 (2): 302–7. <https://doi.org/10.1128/aac.28.2.302>.
- Brun-Buisson, C, A Philippon, M Ansquer, P Legrand, F Montravers, and J Duval. 1987. 'Transferable enzymatic resistance to third-generation cephalosporins during nosocomial outbreak of multiresistant *Klebsiella pneumoniae*.' *The Lancet* 330, no. 8554 302–306.
- Paterson, D L., and R A. Bonomo. 2005. 'Extended-spectrum β-lactamases: a clinical update.' *Clinical Microbiology Reviews* 18, no. 4 657–686.
- Al-Bayssari, C., Dabboussi, F., Hamze, M., and Rolain, J.-M. 2015. 'Detection of expanded-spectrum β-lactamases in Gram-negative bacteria in the 21st century.' *Expert Rev. Anti Infect. Ther.* 13, 1139–1158. <https://doi.org/10.1586/14787210.2015.1066247>.
- Kawamura, Kumiko, Noriaki Nagano, Masahiro Suzuki, Jun-ichi Wachino, Kouji Kimura, and Yoshichika Arakawa. 2017. 'ESBL-Producing *Escherichia Coli* and Its Rapid Rise among Healthy People'. *Food Safety (Tokyo, Japan)* 5 (4): 122–50. <https://doi.org/10.14252/foodsafetyfsci.2017011>.
- Hu, Yanhong Jessika, Anju Ogyu, Benjamin J. Cowling, Keiji Fukuda, and Herbert H. Pang. 2019. 'Available Evidence of Antibiotic Resistance from Extended-Spectrum β-Lactamase-Producing Enterobacteriaceae in Paediatric Patients in 20 Countries: A Systematic Review and Meta-Analysis'. *Bulletin of the World Health Organization* 97 (7): 486–501B. <https://doi.org/10.2471/BLT.18.225698>.
- PHE. UK Standards for Microbiology Investigations: Detection of Enterobacteriaceae producing extended spectrum β-lactamases. Standards Unit, Microbiology Services, PHE. Bacteriology. B 59. Issue no: 4.1, Issue date: 17.08.16.
- Data on file.
- Huang, Te-Din, Pierre Bogaerts, Catherine Berhin, Amelie Guisset, and Youri Glupczynski. 2010. 'Evaluation of *Brilliance*™ ESB� Agar, a Novel Chromogenic Medium for Detection of Extended-Spectrum-Beta- Lactamase-Producing Enterobacteriaceae'. *Journal of Clinical Microbiology* 48 (6): 2091–96. <https://doi.org/10.1128/JCM.02342-09>.

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Bibliográfia

- Farber B, Moellering RC Jr. 1982. 'The third generation cephalosporins.' *Bull N Y Acad Med.* 58: 696–710. PMID: 6762896.
- Kliebe, C., B. A. Nies, J. F. Meyer, R. M. Tolxdorff-Neutzling, and B. Wiedemann. 1985. 'Evolution of Plasmid-Coded Resistance to Broad-Spectrum Cephalosporins'. *Antimicrobial Agents and Chemotherapy* 28 (2): 302–7. <https://doi.org/10.1128/aac.28.2.302>.
- Brun-Buisson, C, A Philippon, M Ansquer, P Legrand, F Montravers, and J Duval. 1987. 'Transferable enzymatic resistance to third-generation cephalosporins during nosocomial outbreak of multiresistant *Klebsiella pneumoniae*.' *The Lancet* 330, no. 8554 302-306.
- Paterson, D L., and R A. Bonomo. 2005. 'Extended-spectrum β-lactamases: a clinical update.' *Clinical Microbiology Reviews* 18, no. 4 657-686.
- Al-Bayssari, C., Dabboussi, F., Hamze, M., and Rolain, J.-M. 2015. 'Detection of expanded-spectrum β-lactamases in Gram-negative bacteria in the 21st century.' *Expert Rev. Anti Infect. Ther.* 13, 1139–1158. <https://doi.org/10.1586/14787210.2015.1066247>.
- Kawamura, Kumiko, Noriyuki Nagano, Masahiro Suzuki, Jun-ichi Wachino, Kouji Kimura, and Yoshichika Arakawa. 2017. 'ESBL-Producing *Escherichia Coli* and Its Rapid Rise among Healthy People'. *Food Safety (Tokyo, Japan)* 5 (4): 122–50. <https://doi.org/10.14252/foodsafetyfscj.2017011>.
- Hu, Yanhong Jessika, Anju Ogyu, Benjamin J. Cowling, Keiji Fukuda, and Herbert H. Pang. 2019. 'Available Evidence of Antibiotic Resistance from Extended-Spectrum β-Lactamase-Producing Enterobacteriaceae in Paediatric Patients in 20 Countries: A Systematic Review and Meta-Analysis'. *Bulletin of the World Health Organization* 97 (7): 486-501B. <https://doi.org/10.2471/BLT.18.225698>.
- PHE. UK Standards for Microbiology Investigations: Detection of Enterobacteriaceae producing extended spectrum β-lactamases. Standards Unit, Microbiology Services, PHE Bacteriology. B 59. Issue no: 4.1, Issue date: 17.08.16.
- Data on file.
- Huang, Te-Din, Pierre Bogaerts, Catherine Berhin, Amelie Guisset, and Youri Glupczynski. 2010. 'Evaluation of *Brilliance*™ ESBL Agar, a Novel Chromogenic Medium for Detection of Extended-Spectrum-Beta- Lactamase-Producing Enterobacteriaceae'. *Journal of Clinical Microbiology* 48 (6): 2091–96. <https://doi.org/10.1128/JCM.02342-09>.

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- Kliebe, C., B. A. Nies, J. F. Meyer, R. M. Tolxdorff-Neutzling, and B. Wiedemann. 1985. 'Evolution of Plasmid-Coded Resistance to Broad-Spectrum Cephalosporins'. *Antimicrobial Agents and Chemotherapy* 28 (2): 302–7. <https://doi.org/10.1128/aac.28.2.302>.
- Brun-Buisson, C, A Philippon, M Ansquer, P Legrand, F Montravers, and J Duval. 1987. 'Transferable enzymatic resistance to third-generation cephalosporins during nosocomial outbreak of multiresistant *Klebsiella pneumoniae*.' *The Lancet* 330, no. 8554 302-306.
- Paterson, D L., and R A. Bonomo. 2005. 'Extended-spectrum β -lactamases: a clinical update.' *Clinical Microbiology Reviews* 18, no. 4 657-686.
- Al-Bayssari, C., Dabboussi, F., Hamze, M., and Rolain, J.-M. 2015. 'Detection of expanded-spectrum β -lactamases in Gram-negative bacteria in the 21st century.' *Expert Rev. Anti Infect. Ther.* 13, 1139–1158. <https://doi.org/10.1586/14787210.2015.1066247>.
- Kawamura, Kumiko, Noriyuki Nagano, Masahiro Suzuki, Jun-ichi Wachino, Kouji Kimura, and Yoshichika Arakawa. 2017. 'ESBL-Producing *Escherichia Coli* and Its Rapid Rise among Healthy People'. *Food Safety (Tokyo, Japan)* 5 (4): 122–50. <https://doi.org/10.14252/foodsafetyfscj.2017011>.
- Hu, Yanhong Jessika, Anju Ogyu, Benjamin J. Cowling, Keiji Fukuda, and Herbert H. Pang. 2019. 'Available Evidence of Antibiotic Resistance from Extended-Spectrum β -Lactamase-Producing Enterobacteriaceae in Paediatric Patients in 20 Countries: A Systematic Review and Meta-Analysis'. *Bulletin of the World Health Organization* 97 (7): 486-501B. <https://doi.org/10.2471/BLT.18.225698>.
- PHE. UK Standards for Microbiology Investigations: Detection of Enterobacteriaceae producing extended spectrum β -lactamases. Standards Unit, Microbiology Services, PHE Bacteriology. B 59. Issue no: 4.1, Issue date: 17.08.16.
- Data on file.
- Huang, Te-Din, Pierre Bogaerts, Catherine Berhin, Amelie Guisset, and Youri Glupczynski. 2010. 'Evaluation of Brilliance™ ESBL Agar, a Novel Chromogenic Medium for Detection of Extended-Spectrum-Beta-Lactamase-Producing Enterobacteriaceae'. *Journal of Clinical Microbiology* 48 (6): 2091–96. <https://doi.org/10.1128/JCM.02342-09>.

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Literatūra

- Farber B, Moellering RC Jr. 1982. 'The third generation cephalosporins.' *Bull N Y Acad Med*. 58: 696–710. PMID: 6762896.
- Kliebe, C., B. A. Nies, J. F. Meyer, R. M. Tolxdorff-Neutzing, and B. Wiedemann. 1985. 'Evolution of Plasmid-Coded Resistance to Broad-Spectrum Cephalosporins'. *Antimicrobial Agents and Chemotherapy* 28 (2): 302–7. <https://doi.org/10.1128/aac.28.2.302>.
- Brun-Buisson, C, A Philippon, M Ansquer, P Legrand, F Montravers, and J Duval. 1987. 'Transferable enzymatic resistance to third-generation cephalosporins during nosocomial outbreak of multiresistant *Klebsiella pneumoniae*.' *The Lancet* 330, no. 8554 302-306.
- Paterson, D L., and R A. Bonomo. 2005. 'Extended-spectrum β -lactamases: a clinical update.' *Clinical Microbiology Reviews* 18, no. 4 657-686.
- Al-Bayssari, C., Dabboussi, F., Hamze, M., and Rolain, J.-M. 2015. 'Detection of expanded-spectrum β -lactamases in Gram-negative bacteria in the 21st century.' *Expert Rev. Anti Infect. Ther.* 13, 1139–1158. <https://doi.org/10.1586/14787210.2015.1066247>.
- Kawamura, Kumiko, Noriyuki Nagano, Masahiro Suzuki, Jun-ichi Wachino, Kouji Kimura, and Yoshichika Arakawa. 2017. 'ESBL-Producing *Escherichia Coli* and Its Rapid Rise among Healthy People'. *Food Safety* (Tokyo, Japan) 5 (4): 122–50. <https://doi.org/10.14252/foodsafetyfscj.2017011>.
- Hu, Yanhong Jessika, Anju Ogyu, Benjamin J. Cowling, Keiji Fukuda, and Herbert H. Pang. 2019. 'Available Evidence of Antibiotic Resistance from Extended-Spectrum β -Lactamase-Producing Enterobacteriaceae in Paediatric Patients in 20

- PHE. UK Standards for Microbiology Investigations: Detection of Enterobacteriaceae producing extended spectrum β -lactamases. Standards Unit, Microbiology Services, PHE Bacteriology. B 59. Issue no: 4.1, Issue date: 17.08.16.
- Data on file.
- Huang, Te-Din, Pierre Bogaerts, Catherine Berhin, Amelie Guisset, and Youri Glupczynski. 2010. 'Evaluation of Brilliance™ ESBL Agar, a Novel Chromogenic Medium for Detection of Extended-Spectrum-Beta- Lactamase-Producing Enterobacteriaceae'. *Journal of Clinical Microbiology* 48 (6): 2091–96. <https://doi.org/10.1128/JCM.02342-09>.

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Bibliografia

- Farber B, Moellering RC Jr. 1982. 'The third generation cephalosporins.' *Bull N Y Acad Med*. 58: 696–710. PMID: 6762896.
- Kliebe, C., B. A. Nies, J. F. Meyer, R. M. Tolxdorff-Neutzing, and B. Wiedemann. 1985. 'Evolution of Plasmid-Coded Resistance to Broad-Spectrum Cephalosporins'. *Antimicrobial Agents and Chemotherapy* 28 (2): 302–7. <https://doi.org/10.1128/aac.28.2.302>.
- Brun-Buisson, C, A Philippon, M Ansquer, P Legrand, F Montravers, and J Duval. 1987. 'Transferable enzymatic resistance to third-generation cephalosporins during nosocomial outbreak of multiresistant *Klebsiella pneumoniae*.' *The Lancet* 330, no. 8554 302-306.
- Paterson, D L., and R A. Bonomo. 2005. 'Extended-spectrum β -lactamases: a clinical update.' *Clinical Microbiology Reviews* 18, no. 4 657-686.
- Al-Bayssari, C., Dabboussi, F., Hamze, M., and Rolain, J.-M. 2015. 'Detection of expanded-spectrum β -lactamases in Gram-negative bacteria in the 21st century.' *Expert Rev. Anti Infect. Ther.* 13, 1139–1158. <https://doi.org/10.1586/14787210.2015.1066247>.
- Kawamura, Kumiko, Noriyuki Nagano, Masahiro Suzuki, Jun-ichi Wachino, Kouji Kimura, and Yoshichika Arakawa. 2017. 'ESBL-Producing *Escherichia Coli* and Its Rapid Rise among Healthy People'. *Food Safety* (Tokyo, Japan) 5 (4): 122–50. <https://doi.org/10.14252/foodsafetyfsci.2017011>.
- Hu, Yanhong Jessika, Anju Ogyu, Benjamin J. Cowling, Keiji Fukuda, and Herbert H. Pang. 2019. 'Available Evidence of Antibiotic Resistance from Extended-Spectrum β -Lactamase-Producing Enterobacteriaceae in Paediatric Patients in 20 Countries: A Systematic Review and Meta-Analysis'. *Bulletin of the World Health Organization* 97 (7): 486-501B. <https://doi.org/10.2471/BLT.18.225698>.
- PHE. UK Standards for Microbiology Investigations: Detection of Enterobacteriaceae producing extended spectrum β -lactamases. Standards Unit, Microbiology Services, PHE Bacteriology. B 59. Issue no: 4.1, Issue date: 17.08.16.
- Data on file.
- Huang, Te-Din, Pierre Bogaerts, Catherine Berhin, Amelie Guisset, and Youri Glupczynski. 2010. 'Evaluation of *Brilliance*™ ESBL Agar, a Novel Chromogenic Medium for Detection of Extended-Spectrum-Beta- Lactamase-Producing Enterobacteriaceae'. *Journal of Clinical Microbiology* 48 (6): 2091–96. <https://doi.org/10.1128/JCM.02342-09>.

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
Bibliografia

- Farber B, Moellering RC Jr. 1982. 'The third generation cephalosporins.' Bull N Y Acad Med. 58: 696–710. PMID: 6762896.
- Kliebe, C., B. A. Nies, J. F. Meyer, R. M. Tolxdorff-Neutzling, and B. Wiedemann. 1985. 'Evolution of Plasmid-Coded Resistance to Broad-Spectrum Cephalosporins'. Antimicrobial Agents and Chemotherapy 28 (2): 302–7. <https://doi.org/10.1128/aac.28.2.302>.
- Brun-Buisson, C, A Philippon, M Ansquer, P Legrand, F Montravers, and J Duval. 1987. 'Transferable enzymatic resistance to third-generation cephalosporins during nosocomial outbreak of multiresistant *Klebsiella pneumoniae*.' The Lancet 330, no. 8554 302-306.
- Paterson, D L., and R A. Bonomo. 2005. 'Extended-spectrum β -lactamases: a clinical update.' Clinical Microbiology Reviews 18, no. 4 657-686.
- Al-Bayssari, C., Dabboussi, F., Hamze, M., and Rolain, J.-M. 2015. 'Detection of expanded-spectrum β -lactamases in Gram-negative bacteria in the 21st century.' Expert Rev. Anti Infect. Ther. 13, 1139–1158. <https://doi.org/10.1586/14787210.2015.1066247>.
- Kawamura, Kumiko, Noriyuki Nagano, Masahiro Suzuki, Jun-Ichi Wachino, Kouji Kimura, and Yoshichika Arakawa. 2017. 'ESBL-Producing *Escherichia Coli* and Its Rapid Rise among Healthy People'. Food Safety (Tokyo, Japan) 5 (4): 122–50. <https://doi.org/10.14252/foodsafetyfscj.2017011>.
- Hu, Yanhong Jessika, Anju Ogyu, Benjamin J. Cowling, Keiji Fukuda, and Herbert H. Pang. 2019. 'Available Evidence of Antibiotic Resistance from Extended-Spectrum β -Lactamase-Producing Enterobacteriaceae in Paediatric Patients in 20 Countries: A Systematic Review and Meta-Analysis'. Bulletin of the World Health Organization 97 (7): 486-501B. <https://doi.org/10.2471/BLT.18.225698>.
- PHE. UK Standards for Microbiology Investigations: Detection of Enterobacteriaceae producing extended spectrum β -lactamases. Standards Unit, Microbiology Services, PHE Bacteriology. B 59. Issue no: 4.1, Issue date: 17.08.16.
- Data on file.
- Huang, Te-Din, Pierre Bogaerts, Catherine Berhin, Amelie Guisset, and Youri Glupczynski. 2010. 'Evaluation of Brilliance™ ESBL Agar, a Novel Chromogenic Medium for Detection of Extended-Spectrum-Beta- Lactamase-Producing Enterobacteriaceae'. Journal of Clinical Microbiology 48 (6): 2091–96. <https://doi.org/10.1128/JCM.02342-09>.

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Bibliografie

- Farber B, Moellering RC Jr. 1982. 'The third generation cephalosporins.' *Bull N Y Acad Med*. 58: 696–710. PMID: 6762896.
- Kliebe, C., B. A. Nies, J. F. Meyer, R. M. Tolxdorff-Neutzling, and B. Wiedemann. 1985. 'Evolution of Plasmid-Coded Resistance to Broad-Spectrum Cephalosporins.' *Antimicrobial Agents and Chemotherapy* 28 (2): 302–7. <https://doi.org/10.1128/aac.28.2.302>.
- Brun-Buisson, C, A Philippon, M Ansquer, P Legrand, F Montravers, and J Duval. 1987. 'Transferable enzymatic resistance to third-generation cephalosporins during nosocomial outbreak of multiresistant *Klebsiella pneumoniae*.' *The Lancet* 330, no. 8554 302–306.
- Paterson, D L., and R A. Bonomo. 2005. 'Extended-spectrum β -lactamases: a clinical update.' *Clinical Microbiology Reviews* 18, no. 4 657–686.
- Al-Bayssari, C., Dabboussi, F., Hamze, M., and Rolain, J.-M. 2015. 'Detection of expanded-spectrum β -lactamases in Gram-negative bacteria in the 21st century.' *Expert Rev. Anti Infect. Ther.* 13, 1139–1158. <https://doi.org/10.1586/14787210.2015.1066247>.
- Kawamura, Kumiko, Noriyuki Nagano, Masahiro Suzuki, Jun-ichi Wachino, Kouji Kimura, and Yoshichika Arakawa. 2017. 'ESBL-Producing *Escherichia Coli* and Its Rapid Rise among Healthy People'. *Food Safety* (Tokyo, Japan) 5 (4): 122–50. <https://doi.org/10.14252/foodsafetyfscj.2017011>.
- Hu, Yanhong Jessika, Anju Ogyu, Benjamin J. Cowling, Keiji Fukuda, and Herbert H. Pang. 2019. 'Available Evidence of Antibiotic Resistance from Extended-Spectrum β -Lactamase-Producing Enterobacteriaceae in Paediatric Patients in 20 Countries: A Systematic Review and Meta-Analysis'. *Bulletin of the World Health Organization* 97 (7): 486–501B. <https://doi.org/10.2471/BLT.18.225698>.
- PHE. UK Standards for Microbiology Investigations: Detection of Enterobacteriaceae producing extended spectrum β -lactamases. Standards Unit, Microbiology Services, PHE Bacteriology. B 59. Issue no: 4.1, Issue date: 17.08.16.
- Data on file.
- Huang, Te-Din, Pierre Bogaerts, Catherine Berhin, Amelie Guisset, and Youri Glupczynski. 2010. 'Evaluation of *Brilliance*™ ESBL Agar, a Novel Chromogenic Medium for Detection of Extended-Spectrum-Beta- Lactamase-Producing

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- Farber B, Moellering RC Jr. 1982. 'The third generation cephalosporins.' Bull N Y Acad Med. 58: 696-710. PMID: 6762896.
- Kliebe, C., B. A. Nies, J. F. Meyer, R. M. Tolxdorff-Neutzling, and B. Wiedemann. 1985. 'Evolution of Plasmid-Coded Resistance to Broad-Spectrum Cephalosporins'. Antimicrobial Agents and Chemotherapy 28 (2): 302-7. <https://doi.org/10.1128/aac.28.2.302>.
- Brun-Buisson, C, A Philippon, M Ansquer, P Legrand, F Montravers, and J Duval. 1987. 'Transferable enzymatic resistance to third-generation cephalosporins during nosocomial outbreak of multiresistant Klebsiella pneumoniae.' The Lancet 330, no. 8554 302-306.
- Paterson, D L., and R A. Bonomo. 2005. 'Extended-spectrum β-lactamases: a clinical update.' Clinical Microbiology Reviews 18, no. 4 657-686.
- Al-Bayssari, C., Dabboussi, F., Hamze, M., and Rolain, J.-M. 2015. 'Detection of expanded-spectrum β-lactamases in Gram-negative bacteria in the 21st century.' Expert Rev. Anti Infect. Ther. 13, 1139-1158. <https://doi.org/10.1586/14787210.2015.1066247>.
- Kawamura, Kumiko, Noriyuki Nagano, Masahiro Suzuki, Jun-ichi Wachino, Kouji Kimura, and Yoshichika Arakawa. 2017. 'ESBL-Producing Escherichia Coli and Its Rapid Rise among Healthy People'. Food Safety (Tokyo, Japan) 5 (4): 122-50. <https://doi.org/10.14252/foodsafetyfsci.2017011>.
- Hu, Yanhong Jessika, Anju Ogyu, Benjamin J. Cowling, Keiji Fukuda, and Herbert H. Pang. 2019. 'Available Evidence of Antibiotic Resistance from Extended-Spectrum β-Lactamase-Producing Enterobacteriaceae in Paediatric Patients in 20 Countries: A Systematic Review and Meta-Analysis'. Bulletin of the World Health Organization 97 (7): 486-501B. <https://doi.org/10.2471/BLT.18.225698>.
- PHE. UK Standards for Microbiology Investigations: Detection of Enterobacteriaceae producing extended spectrum β-lactamases. Standards Unit, Microbiology Services, PHE Bacteriology. B 59. Issue no: 4.1, Issue date: 17.08.16.
- Data on file.
- Huang, Te-Din, Pierre Bogaerts, Catherine Berhin, Amelie Guisset, and Youri Glupczynski. 2010. 'Evaluation of Brilliance™ ESBL Agar, a Novel Chromogenic Medium for Detection of Extended-Spectrum-Beta-Lactamase-Producing

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Bibliografi

- Farber B, Moeller RC Jr. 1982. 'The third generation cephalosporins.' *Bull N Y Acad Med*. 58: 696–710. PMID: 6762896.
- Kliebe, C., B. A. Nies, J. F. Meyer, R. M. Tolxdorff-Neutzling, and B. Wiedemann. 1985. 'Evolution of Plasmid-Coded Resistance to Broad-Spectrum Cephalosporins.' *Antimicrobial Agents and Chemotherapy* 28 (2): 302–7. <https://doi.org/10.1128/aac.28.2.302>.
- Brun-Buisson, C, A Philippon, M Ansquer, P Legrand, F Montravers, and J Duval. 1987. 'Transferable enzymatic resistance to third-generation cephalosporins during nosocomial outbreak of multiresistant Klebsiella pneumoniae.' *The Lancet* 330, no. 8554 302–306.
- Paterson, D L., and R A. Bonomo. 2005. 'Extended-spectrum β -lactamases: a clinical update.' *Clinical Microbiology Reviews* 18, no. 4 657–686.
- Al-Bayssari, C., Dabboussi, F., Hamze, M., and Rolain, J.-M. 2015. 'Detection of expanded-spectrum β -lactamases in Gram-negative bacteria in the 21st century.' *Expert Rev. Anti Infect. Ther.* 13, 1139–1158. <https://doi.org/10.1586/14787210.2015.1066247>.
- Kawamura, Kumiko, Noriyuki Nagano, Masahiro Suzuki, Jun-ichi Wachino, Kouji Kimura, and Yoshichika Arakawa. 2017. 'ESBL-Producing Escherichia Coli and Its Rapid Rise among Healthy People'. *Food Safety (Tokyo, Japan)* 5 (4): 122–50. <https://doi.org/10.14252/foodsafetyfscj.2017011>.
- Hu, Yanhong Jessika, Anju Ogyu, Benjamin J. Cowling, Keiji Fukuda, and Herbert H. Pang. 2019. 'Available Evidence of Antibiotic Resistance from Extended-Spectrum β -Lactamase-Producing Enterobacteriaceae in Paediatric Patients in 20 Countries: A Systematic Review and Meta-Analysis'. *Bulletin of the World Health Organization* 97 (7): 486–501B. <https://doi.org/10.2471/BLT.18.225698>.
- PHE. UK Standards for Microbiology Investigations: Detection of Enterobacteriaceae producing extended spectrum β -lactamases. Standards Unit, Microbiology Services, PHE Bacteriology. B 59. Issue no: 4.1, Issue date: 17.08.16.
- Data on file.
- Huang, Te-Din, Pierre Bogaerts, Catherine Berhin, Amelie Guisset, and Youri Glupczynski. 2010. 'Evaluation of Brilliance™ ESBL Agar, a Novel Chromogenic Medium for Detection of Extended-Spectrum-Beta- Lactamase-Producing Enterobacteriaceae'. *Journal of Clinical Microbiology* 48 (6): 2091–96. <https://doi.org/10.1128/JCM.02342-09>.

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Oxoid GmbH,