<https://www.thermofisher.com>

EU  
US  
CA  
ROW

+800 135 79 135  
1 855 236 0910  
1 855 805 8539  
+31 20 794 7071

## Oxoid™ Optochin 5µg (OP5) Diagnostic Discs

**REF DD0001B & DD0001T**

### Intended Use

Impregnated with 'Optochin' (ethylhydrocuprein hydrochloride) and used in a diagnostic workflow to aid Microbiology professionals in identifying potential pathogenic bacteria suspected in patients having a *Streptococcus pneumoniae* infection. These discs provide a qualitative result that enables the differentiation of *Streptococcus pneumoniae* from other alpha-haemolytic streptococci. To be used with a pure, agar grown culture, the device is for professional use only and is neither automated, nor a companion diagnostic. This device is an in vitro diagnostic.

### Principle of the Test

Bowers and Jeffries<sup>1</sup> have shown that there is complete correlation between bile-solubility and full 'Optochin' susceptibility for the differentiation of *Streptococcus pneumoniae* from other streptococci.

Optochin (5 µg) OP5 Diagnostic Discs (DD0001) are paper discs impregnated with 'Optochin' (ethylhydrocuprein hydrochloride) which provide a convenient and reliable alternative to the bile-solubility test. *Pneumococci* are sensitive to 'Optochin' so that when grown on test medium, they exhibit a zone of inhibition around the impregnated disc, whilst other alpha-haemolytic streptococci either grow right up to the edge of the disc or, occasionally, show a very small zone of inhibition.

Pure cultures of clinical isolates of *Streptococcus pneumoniae* are inoculated onto test medium and the Optochin disc is placed on the surface. The agent within the disc diffuses through the agar to form a gradient. After incubation, the zone around the disc is measured.

A culture of *Streptococcus pneumoniae* shows a zone of inhibition around the impregnated disc of at least 14mm, whilst other alpha-haemolytic streptococci either grow right up to the edge of the disc or, occasionally, show a very small zone of inhibition of not more than 2mm.

### Components

Optochin OP5 Diagnostic Discs consist of 6mm diameter paper discs impregnated with a specific antimicrobial concentration. The discs are marked on both sides with an alpha-numeric code identifying the antimicrobial agent and concentration. Optochin OP5 Diagnostic Discs are supplied in cartridges of 50 discs, either as a single cartridge or as a pack of five. Cartridges are individually packed in a foil-sealed blister pack with a desiccant.

DD0001T.....	Σ 50
DD0001B.....	Σ 250

### Materials Provided

- 1 x cartridge with spring cap and plunger (DD0001T)
- 5 x cartridge with spring cap and plunger (DD0001B)
- Desiccant tablet x1 (DD0001T)
- Desiccant tablet x 5 (DD0008T)
- Foil individually sealing each cartridge with its desiccant
- Optochin discs x 50 (DD0001T)
- Optochin discs x 250 (DD0008T)

Optochin diagnostic discs contain the active ingredient 'Optochin' (ethylhydrocupreine hydrochloride) which is an Antimicrobial agent, typically white to light yellow when in dry form.

### Materials Required but Not Supplied

Agar plates with appropriate media, inoculum suspension medium, sterile loops and swabs, sterile forceps, McFarland turbidity standards, incubator, modified atmosphere environments, antibiotic disc dispensers, quality control strains, apparatus to measure zone sizes and interpretative criteria for local standard methods.

### Warnings and Precautions

- Not intended for susceptibility testing.
- This product is for *in vitro* diagnostic use only.
- Follow instructions for use.

- Observe aseptic techniques and established precautions against all microbiological hazards throughout all procedures.
- Cultures, containers and other contaminated materials must be sterilized after use in accordance with guidelines for the handling and disposal of biohazardous waste.
- Refer to the Material Safety Data Sheet (MSDS) for safe handling and disposal of the product (found on [www.thermofisher.com](http://www.thermofisher.com)).
- Inspect the product packaging before first use. Do not use product if there is any visible damage to the foil seal.
- Do not use beyond the stated expiry date.
- Once the cartridge is open, ensure it is stored in an opaque desiccated environment to prevent degradation.
- If the discs do not produce the expected reaction with recommended control microorganisms, review the entire procedure.
- In the event of malfunction do not use the device.
- The device does not contain any carcinogenic, mutagenic or toxic to reproduction, endocrine disrupting substances or materials likely to cause sensitisation or an allergic reaction under normal use.

### Storage and Handling

Unopened cartridges must be stored at -20 to +8°C until required. Allow cartridges to reach room temperature before removing them from the packaging to minimize condensation. The expiry date is valid only for the unopened blister packs stored under correct conditions. Once opened, cartridges should be stored within a dispenser (sold separately) in the container provided (with an unsaturated (orange) desiccant), or other suitable opaque airtight container with a desiccant to protect the discs from light and moisture. Once opened from their desiccant-containing packaging, discs should be used within 7 days and only if stored as described in this IFU. Dispensers should be stored within the container at 2 - 8°C and be allowed to come to room temperature before opening to prevent the formation of condensation.

### Method

1. Specimens should be collected and handled following recommended guidelines<sup>3</sup>.
2. Inoculate a blood agar plate using a loop or swab with a pure culture of the test organism exhibiting the appropriate haemolysis from the primary isolation plate.
3. Using either sterile forceps, single disc ejector or disc dispenser; apply the disc/s to the surface of the pre-inoculated agar.
4. Lightly touch the disc with the sterile forceps to ensure complete contact between the disc and agar surface.
5. Invert and incubate plates at 35-37°C in 5% CO<sub>2</sub> for 18-24 hours.
6. Inspect for a zone of inhibition around the discs. If suspect colonies of pneumococci are isolated from the primary plating medium, repeat the described procedure with a pure culture of the suspect colonies.

### Interpretation

A zone of inhibition greater than 14 mm from an alpha-hemolytic isolate surrounding the optochin disc is indicative of *S. pneumoniae*. The growth of alpha-haemolytic streptococci is not usually inhibited by optochin. A diameter between 6 and 14 mm is questionable for *pneumococci* and the strain should be presumptively identified as a *pneumococcus* only if it is bile soluble<sup>4</sup>.

### User Quality Control Procedures

The use of positive control organisms such as *S. pneumoniae* ATCC® 6305 is recommended. *Streptococcus pyogenes* ATCC® 19615 may be used as a negative control organism. Testing of control organisms should be performed in accordance with established laboratory quality control procedures. If aberrant quality control results are noted, patient results should not be reported.

### Limitations

For professional use only.

This product is for in vitro diagnostic use only. It is recommended that the tests are performed on colonies from pure culture. Alpha-haemolytic streptococci susceptibility to optochin should be reported as presumptive identification of *S. pneumoniae*. Final identification of *S. pneumoniae* is made by additional tests including, bile solubility, inulin fermentation, the Neufeld (Quellung) reaction, and Meuso Virulence<sup>4</sup>. Optochin resistance has been reported in a small percentage of strains<sup>5</sup>. Reported test results should form part of a holistic approach to treatment and the reported test results should

be assessed by the clinician in conjunction with the patient's history, clinical presentations and results of other clinical tests.

Failure to observe instructions for use may lead to erroneous results.

#### Analytical Performance Characteristics

Repeatability and reproducibility were measured internally by testing *Streptococcus pneumoniae* (ATCC 49619) and *Streptococcus pyogenes* (ATCC 19615) on 495 and 189 batches of DD0001B and DD0001T respectively. Results were deemed comparable when either no deviation, or a deviation no greater than 2 mm, between the test discs was observed. Across all testing, >99% of results were deemed comparable. (Data on file.)

#### Clinical Performance Characteristics

A national multicentric surveillance study<sup>6</sup> was undertaken to collect isolates from 30 clinical microbiology laboratories across Portugal. Optochin susceptibility was tested via the use of the Oxoid Optochin (5 µg) OP5 Diagnostic Discs (DD0001) on blood agar plates. 100% of isolates collected from 1999 to 2003 demonstrated optochin susceptibility. Isolates collected in 2004 and 2005 showed a small increase in resistance to 0.3% and 3.2% respectively.

Another study<sup>7</sup> collected nasopharyngeal samples from asymptomatic children aged 4 months to 6 years in 2001, 2002, 2003 and in 2006. Of the 1,973 *S. pneumoniae* isolates in the study, resistance (by year) to optochin ranged from 1.3 to 3.2%.

A total of 266 strains were isolated in 966 healthy children under five years old<sup>8</sup>. Three (1.1%) of the 266 isolates of *S. pneumoniae* were found to be optochin resistant.

#### Serious Incidents

Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the relevant regulatory authority in which the user and/or the patient is established.

#### Contact Information

International Telephone: +44 (0) 1256 841144  
Email: [microbiology@thermofisher.com](mailto:microbiology@thermofisher.com)  
Website: [www.thermofisher.com](http://www.thermofisher.com)

#### Glossary of Symbols

Symbol/Label	Meaning
	Manufacturer
	In Vitro Diagnostic Medical Device
	Temperature limit
	Batch Code
	Catalog Number
	Do not re-use
	Consult instructions for use or consult electronic instructions for use
	Contains sufficient for <n> tests
	Use-by date

	Do not use if package is damaged and Consult instructions for use
	Authorized representative in the European Community/ European Union
	Unique device identifier
	USA: Caution: Federal law restricts this device to sale by or on order of a Physician
	European Conformity Mark
	UK Conformity Mark

#### References

1. Bowers, E.F. and L.R. Jeffries. 1955. *J. Clin. Pathol.* 8:58.
2. Bowen, M.K., L.C. Thiele, B.D. Stearman, and I.G. Schaub. 1957. *J. Lab. Clin. Med.* 49:641.
3. Murray, P.R., E.J. Baron, J.H. Jorgensen, M.A. Pfaller, and R.H. Tenover. 2003. *Manual of Clinical Microbiology*, 8th ed. ASM, Washington, D.C.
4. Lennette, E.H., A. Balows, W.J. Hausler, and J.P. Truant. 1980. *Manual of Clinical Microbiology*, 3rd ed. ASM, Washington, D.C.
5. Nunes S, Sá-Leão R, de Lencastre H. 2008. *J. Clin. Microbiol.*;46(1):321-324
6. Aguiar et al. 2006. *Microbial Drug Resistance* 12 (4): 239-245
7. Nunes, Sá-Leão, and De Lencastre.2008. *Journal of Clinical Microbiology*, 46 (1): 321-324
8. N Jindal and U Arora. 2010. *Indian Journal of Medical Microbiology*. 10.4103/0255-0857.58741.

The ATCC Licensed Derivative® Emblem, the ATCC Licensed Derivative® word mark, and the ATCC catalogue marks are trademarks of ATCC. Thermo Fisher Scientific Inc. is licensed to use these trademarks and to sell products derived from ATCC® cultures.

© 2023 Thermo Fisher Scientific Inc. All rights reserved.  
ATCC® is a trademark of American Type Culture Collection.  
All other trademarks are the property of Thermo Fisher Scientific and its subsidiaries.  
This information is not intended to encourage use of these products in any manner that might infringe the intellectual property rights of others. Please contact your local sales representative for all enquiries.

#### Revision Information

Version	Date of issue and modifications introduced
1.0	2023-01-18

Oxoid Limited, Wade Road, Basingstoke, Hampshire, RG24 8PW, UK

[www.thermofisher.com](http://www.thermofisher.com)

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© 2023 Thermo Fisher Scientific Inc. Med ensamrätt.