Steri-Shield® Personal Protection System Flyte® Hood

Without Peel Away Lens **REF** 0408-800-000
0408-800-400

With Peel Away Lens **REF** 0408-800-100
0408-800-500

Instructions For Use

R_x ONLY

CE 0197

STERILE EO

ENGLISH (EN)
ESPAÑOL (ES)
DEUTSCH (DE)
FRANÇAIS (FR)
ITALIANO (IT)
NEDERLANDS (NL)
SVENSKA (SV)
DANSK (DA)
SUOMI (FI)
PORTUGUÊS (PT)
NORSK (NO)
POLSKI (PL)
ΕΛΛΗΝΙΚΑ (EL)
日本語 (JA)
中文 (ZH)
한국어 (KO)

Introduction

This instructions for use manual is the most comprehensive source of information for the safe and effective use of your product. This manual may be used by in-service trainers, physicians, nurses, and surgical technologists. Keep and consult this reference manual during the life of the product.

The following conventions are used in this manual:

- A **WARNING** highlights a safety-related issue. ALWAYS comply with this information to prevent patient and/or healthcare staff injury.
- A **CAUTION** highlights a product reliability issue. ALWAYS comply with this information to prevent product damage.
- A **NOTE** supplements and/or clarifies procedural information.

For additional information, especially safety information, or in-service training, contact your Stryker sales representative or call Stryker customer service. Outside the US, contact your nearest Stryker subsidiary.

NOTE: The user and/or patient should report any serious product-related incident to both the manufacturer and the Competent Authority of the European Member State where the user and/or patient is established.

Compliance Statements

The CE marking affixed to the equipment indicates compliance with the following European Community Regulations:

- Medical Device Regulation (EU) 2017/745
- Regulation (EU) 2016/425 - Personal Protective Equipment (PPE)

This PPE has been EC-type examined by the following notified body:

BSI Group The Netherlands B.V.
John M. Keynesplein 9
1066 EP Amsterdam, The Netherlands
BSI 2797

The lenses were developed and tested in accordance with the requirements of BS EN 166: 2002.

Marking Definitions

EN166 S1F3

CODE	DEFINITION
EN166	European Norm - Personal Eye Protection Specification
S	Code of Manufacturer, Stryker
1	Optical Class
F	Mechanical Strength, Low Energy Impact
3	Protection against droplets and splashes of liquids

Indications for Use

The Stryker Flyte Hoods are components of the Stryker Personal Protection System and are intended to protect the patient, healthcare personnel and operating room personnel against contamination, exposure of infectious bodily fluids, and the transfer of microorganism and particulate material. The Flyte Hoods are sterile, single use only.

Contraindications

None known.

Limitations of Use

WARNINGS:

- This personal protection system is not intended to be used as a respiratory protection system.
- The lenses of this system provide limited protection to the face from flying particles or the splash of hazardous liquids. If these lenses are worn over ophthalmic spectacles, impacts may be transmitted. The lenses are not shatterproof.
- The hoods are not intended for use against oily particulate such as paint mist, oil mist or detergents.
- Avoid contact with open flame. The hood fabric is not intended to be exposed to open flame and may burn.
- The hoods are not intended for use in atmospheres immediately dangerous to life or health.
- The materials contained within this system are not known allergens. However, if an allergic reaction occurs, contact a physician.

For Use With

Flyte Hoods are components of the Stryker Flyte System and are intended to be worn over any Stryker Flyte Helmet.

DESCRIPTION	REF
Flyte Helmet	0408-600-000
Flyte Helmet with Fiber Optic Headlight	0408-635-000
Flyte Helmet with Battery Powered Headlight	0408-645-000

User/Patient Safety

WARNINGS:

- Only trained and experienced healthcare professionals should use this equipment.
- Before using any system component, or any component compatible with this system, read and understand the instructions. Pay particular attention to WARNING information. Become familiar with the system components prior to use.
- The healthcare professional performing any procedure is responsible for determining the appropriateness of this equipment and the specific technique used for each patient. Stryker, as a manufacturer, does not recommend surgical procedure or technique.
- Upon initial receipt and before use, visually inspect the package for damage to confirm the integrity of the sterile barrier. Do not use the product if damage is apparent, the sterile barrier is compromised, or the package is unintentionally opened. If contamination is suspected, discard the contaminated hood and use a new hood.
- Use only Stryker-approved components and accessories, unless otherwise specified. DO NOT modify any component or accessory.

- DO NOT use in the presence of flammable anesthesia gases.
- DO NOT reuse, reprocess, or repack this device. This device is intended for a single use only. DO NOT clean or re-sterilize this device. Cleaning or re-sterilization may significantly affect the filtration efficiency and result in tears in the garment or other loss of mechanical integrity. Failure to comply may lead to infection or cross-infection and result in patient and/or healthcare staff injury.

Description

The hoods with Peel Away lenses feature a three-layer protective lens system. When the outer two lenses become soiled, they can be peeled away to reveal another sterile lens surface.

Symbol Definitions

The symbols located on the equipment and/or labeling are defined in this section or in the *Symbol Definition Chart*. See the *Symbol Definition Chart* supplied with the equipment.

SYMBOL	DEFINITION
	General warning sign
	Do not re-use
	Lens Slot - Helmet Hook
	Lens Magnet - Helmet Magnet
	This system DOES NOT provide respiratory protection.
TOP	Top of product within wrap



Instructions

To Don the Hood (figures 1 to 6)

WARNINGS:

- ALWAYS wear the helmet and hood or toga as a total system to obtain the required personal protection.
- Before and during use, ALWAYS inspect the equipment for rips, tears, cuts, or holes. ALWAYS make sure seams and closing systems are secure and function properly. Scratched or damaged lenses should be replaced. DO NOT use the equipment if damage is apparent.
- Per hospital protocol, ALWAYS use proper sterile technique, especially when positioning the user's arms to permit donning without contamination.
- DO NOT contaminate the exterior surface of the hood during the removal of the user's face mask.
- ALWAYS follow current hospital protocol governing the safe removal, handling, and disposal of potential biohazard waste.

**NOTES:**

- The user must don a surgical helmet and power pack then scrub before donning a hood. See the instructions for use supplied with the helmet.
- Hoods are worn in conjunction with surgical gowns where the hood material extends beyond the user's upper back, shoulders, and front, while the surgical gown covers the lower portion of the hood up to the user's neck.
- The hood must be donned with the help of a circulating assistant.
- The hood is folded and packaged to permit handling without causing contamination to the exterior surfaces of the hood.
- Be sure to adjust helmet settings and secure cable connections to the helmet prior to donning the hood.

- Make sure the white filter material of the hood or toga is positioned correctly over the fan intake of the helmet.
- When the lens is soiled, see the *To Remove the Peel Away Lens* section.

To Remove the Peel Away Lens (figures A and B)

WARNING: ALWAYS remove the Peel Away lens away from the sterile field.

NOTES:

- The lens may be removed by a circulating assistant or the user.
- Each layer of lens and exposed surfaces are sterile.

Storage and Handling**CAUTIONS:**

- ALWAYS use the original packaging container for transportation of the equipment.
- Store the equipment within the specified environmental condition value(s) throughout its useful life. See the *Specifications* section.

Specifications**MATERIAL STANDARDS COMPLIANCE**

Flammability:	Flammability of Clothing Textiles under United States Consumer Product Safety Commission (US CPSC) 16 Code of Federal Regulations (CFR) Part 1610 Class 1
---------------	---

HOOD MATERIAL

Hood Materials:	Polypropylene, polyester, acrylic
Peel Away Lens Material:	Clear polyester laminate
Standard Lens Material:	Clear polycarbonate
Color:	Blue

ENVIRONMENTAL CONDITIONS

Operation, Storage, and Transportation Temperature Limit:
Temperature exposure upper limit of 40 °C (104 °F).

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
Declaration of Conformity (DOC)
Application of Regulation (EU) 2016/425

Manufacturer: Stryker Instruments
Address: 4100 E. Milham Avenue
Kalamazoo, MI 49001, USA
Phone: 269-323-7700
Fax: 800-253-3210
269-324-5412
European Representative: Stryker European Operations B.V.
Address: Herikerbergweg 110
1101 CM Amsterdam
The Netherlands
Phone: 011-33472453651

Notified Body: BSI Group The Netherlands B.V.
Address: John M. Keynesplein 9
1066 EP Amsterdam
The Netherlands

Product Group Name: Stryker Steri-Shield Personal Protection Systems
Stryker Flyte and Flyte SurgiCool Steri-Shield Personal Protection System

Product Class and Rule: Optical Class 1
Impact Grade F – High Speed particles with low energy impact
Options 3 Protection against droplets and splashes of liquid
Category 2
(See Product List in Attachment to DoC below)

Product Number: We hereby declare under our sole responsibility that these products conform with the relevant provisions of the Personal Protective Equipment Regulation (EU) 2016/425, and where such is the case, with the national standard transposing harmonized standard No EN 166:2001.

BSI Netherlands, 2797, performed the EU type examination (Module B) and issued the EU Type Examination Certificate CE 710628.

I, the undersigned, declare on the basis of the above information that the system described above is in compliance with the requirements of Regulation (EU) 2016/425.

This declaration hereby authorizes the CE Mark to be affixed to the above mentioned product.

Reviewed and Approved By:  **Brittny Larsen**
Manager, Regulatory Affairs
Date: 2020-04-23 15:39:24.09 Date: 23-Apr-2020

Declaration of Conformity Number: DOC-044
Revision Level: AB

Attachment to DOC

T4/T5 Helmet Personal Protection System

Product Number	Product Name	Grade of Lens Material
0400-800-000	T4 Hood	1
0400-800-100	T5 Hood w/ Peel Away Face Shield	3
0400-810-000	T4 Zippered Toga (S/M)	1
0400-820-000	T4 Zippered Toga (L/XL)	1
0400-830-000	T4 Zippered Toga (Regular)	1
0400-840-000	T4 Zippered Toga (Large)	1
0400-850-000	T4 Zippered Toga (XL, Tall)	1
0400-710-000	T4 Pullover Toga (S/M)	1
0400-720-000	T4 Pullover Toga (L/XL)	1
0400-760-000	T4 Pullover Toga (Large)	1
0400-770-000	T4 Pullover Toga (XL)	1

Flyte Helmet Personal Protection System

Product Number	Product Name	Grade of Lens Material
0408-800-000	Flyte Hood	1
0408-800-100	Flyte Hood, Peel-Away	3
0408-800-400	Flyte Hood	1
0408-800-500	Flyte Hood, Peel-Away	3
0408-801-400	Flyte SurgiCool Hood	1
0408-801-500	Flyte SurgiCool Peel Away Hood	3
0408-820-000	Flyte Zippered Toga, Large	1
0408-830-000	Flyte Zippered Toga, X-Large	1
0408-840-000	Flyte Zippered Toga, 2X Large	1
0408-850-000	Flyte Zippered Toga, 3X Large	1
0408-821-000	Flyte SurgiCool Zippered Toga, Large	1
0408-831-000	Flyte SurgiCool Zippered Toga, X Large	1
0408-841-000	Flyte SurgiCool Zippered Toga, 2X Large	1
0408-851-000	Flyte SurgiCool Zippered Toga, 3X Large	1
0408-820-100	Flyte Toga, Peel-Away, Large	3
0408-830-100	Flyte PeelAway Toga, X-Large	3
0408-840-100	Flyte PeelAway Toga, 2X Large	3
0408-850-100	Flyte PeelAway Toga, 3X Large	3
0408-821-100	Flyte SurgiCool Zippered Toga, Peel Away, Large	3
0408-831-100	Flyte SurgiCool Zippered Toga, Peel Away, X Large	3
0408-841-100	Flyte SurgiCool Zippered Toga, Peel Away, 2X Large	3
0408-851-100	Flyte SurgiCool Zippered Toga, Peel Away, 3X Large	3
0408-720-000	Flyte Pullover Toga, Large	1
0408-730-000	Flyte Pullover Toga, X-Large	1
0408-740-000	Flyte Pullover Toga, 2X Large	1
0408-750-000	Flyte Pullover Toga, 3X Large	1
0408-721-000	Flyte SurgiCool Pullover Toga, Large	1
0408-731-000	Flyte SurgiCool Pullover Toga, X Large	1
0408-741-000	Flyte SurgiCool Pullover Toga, 2X Large	1
0408-751-000	Flyte SurgiCool Pullover Toga, 3X Large	1

Stryker Instruments
1941 Stryker Way
Portage, Michigan
(USA) 49002
1-269-323-7700
1-800-253-3210

stryker[®]