(240)900-029(10)0505

ENGLISH

INSTRUCTIONS FOR USE
Accessories for electrosurgical
units

Monopolar electrosurgical
electrodes, reusable

520-100 to 527-100

Monopolar
Electrosurgical Electrode (reusable)

528-650 to 529-050

Monopolar
Electrosurgical Electrode (reusable)

530-000 to 532-400

Monopolar
Electrosurgical Electrode (reusable)

534-000 to 548-200

Monopolar
Electrosurgical Electrode (reusable)

580-210

Needle electrode, angled 90°,
0,6 x 1.5mm, L: 115mm, handle 4mm ,

580-220

Needle electrode, angled 90°,
0,6 x 4mm, L: 115mm, handle 4mm ,

580-240

Knife electrode, angled 45°,
0,4 x 3mm, L: 115mm, handle 4mm ,

580-250

Keying electrode , angled 90°,
1,5 x 3mm, L: 115mm, handle 4mm ,

580-260

Round electrode, brush, L:
115mm, handle 4mm

580-270

Ball electrode straight, 3mm,
L: 115mm, handle 4mm ,

580-280

Round electrode, hollow, L:
115mm, handle 4mm

580-290

Round electrode, flat, L:
115mm, handle 4mm

Adaptors for electrodes, reusable

529-350

Adaptor for 2.4 mm electrodes
to handle 4mm

529-360

Adaptor for 1.6mm electrodes
to handle 4mm

Extensions electrodes, reusable

527-500

Extension for electrodes 4mm,
L: 100mm, handle 4mm

527-550

Extension for electrodes 4mm,
L: 150mm, handle 4mm

533-900

Extension for electrodes
2.4mm, L: 100mm, handle 2.4mm

533-950

Extension for electrodes
2.4mm, L: 150mm, handle 2.4mm

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tel. +48 22 455 66 66
support@emed.pl

Reports of medical incidents:
incidents@emed.pl

- Consult instruction for use
- Protect from moisture
- Protect from sunlight
- LOT

LOT
- REF

Catalogue number
- Caution
- Use to

The electrosurgical accessories from EMED comply with the General safety and performance requirements defined by Council Directive 93/42/EEC and amendment 2007/47/EC.

CAUTION
The present instructions for use do not substitute for those for the electrosurgical unit. Always follow the instructions for use for the electrosurgical unit provided by its producer and the guidance given in it, in particular the contraindications for the performance of electrosurgical procedures. This also applies to the instructions for the accessories used, including electrosurgical neutral electrodes to be used during monopolar application with the electrosurgical unit. Specifications, safety instructions and warnings contained in the respective instructions for use must always be observed and followed.

1. Intended use
Reusable monopolar electrodes are intended for cutting and coagulation in open surgery. They are inserted into suitable handle and connected to EMED electrosurgical units and other compatible units.

The adaptors allow to connect an electrode with a different diameter than the handle's. Electrode extensions allow to extend the monopolar electrode, which is used for cutting and coagulation.

2. Contraindications
Electrosurgery is not recommended for pregnant women or people with: implanted electronic devices (e.g. heart pacemakers, nerve stimulators, cardioverters or hearing implants), metal implants, hypertension, diabetes, blood clotting disorders. The above factors bring a risk of side effects. Accessories can not be used if, in the opinion of an experienced doctor or in

accordance with the current specialist literature, this could create a risk to the patient, for example because of his general state of health or other existing contraindications, and the level of risk to the health or life of the patient would exceed the benefits of their use.

CAUTION
When electrosurgical procedures are carried out there is a risk of interference with the operation of or damage to metal implants or implantable electronic devices (e.g. heart pacemakers, nerve stimulators, cardioverters or hearing implants). In these cases, the use of monopolar techniques is not recommended.
It is not permitted to use electrosurgical procedures on patients with implanted electronic devices in outpatient clinic conditions.

CAUTION
Risk of ignition of the gases in the upper and lower parts of the gastrointestinal tract
Prior to the use of the electrodes, make sure that endogenous gases have not collected on the application site.

CAUTION
Risk of ignition of the gases in the respiratory system
Given the risk of ignition, the following guidance must be followed:
- oxygen or other flammable gases and substances must not be fed into the respiratory system,
- when the patient needs ventilation, stop the coagulation.

CAUTION
Exogenous burns
They are caused by the ignition of flammable liquids or gases. Burns of this type can also occur as a result of an explosion of flammables/explosives. The possible causes of this type of burns include e.g.: the ignition of agents used for cleaning and disinfecting the operating field and anaesthetic gases etc.

CAUTION
Endogenous burns
They are caused by a high current intensity in the patient's tissue. The possible causes of this type of burns include:
- the incorrect positioning of the patient in such a manner that his/her skin touches the electrically conductive parts,
- the direct contact between the conducting cables and the patient's skin.

CAUTION
The electrodes are not intended for direct contact with the heart or central circulatory system.

CAUTION
Destroy all the accessories which may have been contaminated by prions (CJD). It is forbidden to use them again! Only disposable products should be used for treatment of patients with Creutzfeldt-Jakob disease.

3. Technical data
Maximum rated voltage of the accessories:

REF		Vp (max)
520-100 to 524-600, 524-700, 527-500 to 534-450, 540-100 to 542-200, 580-260, 580-290	524-460, 524-800, 534-450, 542-200, 580-280,	4 300 Vp
524-600, 525-200, 526-400, 536-400	524-900 to 525-800 to 534-460 to	500 Vp
527-000, 548-200 to 580-250, 580-270	527-100, 580-250,	1 300 Vp

CAUTION
Power output of the electrosurgical unit
Always set the correct power for a given accessory.
Avoid too high output settings. Check the maximum voltage rating of the accessory in the instructions for use or on the label. If in doubt, contact the manufacturer. The maximum output voltage of the MODE used in electrosurgical unit must not exceed the maximum permissible voltage of the of the connected accessory and must not exceed 4 MHz! Too low or too high settings of the electrosurgical unit may cause an incorrect operation of the device and unintended burns.

CAUTION
To reduce the risk of tissue overheating it is recommended that the current is applied briefly at short intervals.

The electrodes are connected by the handle with a monopolar cable to the monopolar output of the electrosurgical unit.

2.4mm electrodes can be connected directly to handles REF 218-250, REF 327-14S or comparable. They also can be connected to 4 mm handles with the adaptor REF 529-350.

4mm electrodes can be connected directly to handles REF 215-250, REF 325-14S or comparable product.

1.6mm electrodes can be connected with 4mm handles REF 215-250, REF 325-14S or a comparable product using an adapter REF 529-360.

The electrodes are activated by means of a footswitch or electrosurgical handle. The product can be operated with HF generators with monopolar operation mode of the following manufacturers:

ERBE	KLS Martin
Covidien	ValleyLab
EMED	Tekno
BOWA	Berchold

The manufacturer is not responsible for any other use which does not comply with the recommendations.

The type "BF"/"CF" applied part of the electrosurgical unit is extended with the accessory connected to it.

4. Warnings

CAUTION
Electrosurgical accessories can only be used by trained and qualified staff.

CAUTION
WARNING: Do not make any modifications to the product and do not use products which have undergone any modifications. EMED is not responsible for the use of products that have been modified in any way.

CAUTION
Refer to the instructions for use or the label for the maximum rated voltage of the accessory. Do not use settings and modes where the maximum voltage exceeds the permissible maximum voltage of the accessory. In case of doubt, contact the manufacturer.

CAUTION
Before starting a procedure, check if the operation of the unit does not interfere with the operation of the systems to which the patient is connected. The test is carried out by pressing first the yellow (cutting), then the blue (coagulation) footswitch button, while watching, at the same time, the operation of the monitoring systems.

CAUTION
Before each use, carry out a visual inspection for damage or defects on the surface of the electrodes, adaptors and extensions. Do not use products that show any damage/defects.

CAUTION
It is unacceptable to activate the electrodes if they are not in the operator's field of vision and if they are in contact with metal objects and/or optics. Procedures should only be performed if sufficient visibility of the application area is ensured.

CAUTION
It is very important to ensure that no flammable substances (anesthetics, oxidizing gases, endogenous gases, etc.) are present in the immediate vicinity during the complete electrosurgical application; otherwise there is a danger of explosion. Use non-flammable disinfectants agents; do not use e.g. alcohol-based tinctures or similar. All oxygen compounds must be tight and leak-proof during the procedure.

CAUTION
The tip of the electrode should not be touched during use. After the procedure, the electrode tip may still be hot and may cause burns.

CAUTION
Do not put the accessory directly on the patient. Inadvertent activation of the accessory may cause burns to the patient.

CAUTION
Perform procedures only when there is adequate visibility of the application area. Do not activate the accessory without visual control!

CAUTION
There is a risk of physical injury from sharp electrode tips.

CAUTION
Do not use spray coagulation (SPRAY) mode when using a needle electrode.

CAUTION
Regularly remove blood and tissue residues from the accessory during use e.g. with a disposable electrosurgical instrument cleaner (e.g. REF F7520).

CAUTION
When the cutting or coagulation quality deteriorates for a given setting do not increase the power without prior:
- checking the application of the neutral electrode,
- checking of the correct connection of the cable plugs and accessories,
- making sure that the function has been activated with the correct button in the footswitch (yellow-cut, blue-coagulation),
- checking if the insulation of the cables and electrodes is not damaged,
- checking if the active electrode is clean. When it is not clean this can cause a deterioration of the operational quality of the electrosurgical unit. If the coagulation quality does not improve although the recommended actions have been taken, the accessory should be replaced by a new one.

CAUTION
Accessories can potentially cause muscle stimulation, which may cause unexpected movement of the of the patient in the surgical field.

5. Performance of procedures
New accessories are reusable and non-sterile devices.
Prior their first and each subsequent application, electrodes, adaptors and extensions have to be subjected to cleaning, disinfection and sterilization according to a validated procedure (EN ISO 17665) (see section "Processing - Cleaning, Disinfection and Sterilization").

CAUTION
Immediately after purchasing a new product, before and after the procedure, carefully examine the condition of the elctrode, adaptor and extension paying particular attention to any damage/breakage on their surface. Before each use, test the operation of the product after connecting to the unit, press the blue button of the handle or footswitch

and then the yellow button and see if it activates. If visible defects or malfunctions are found, contact the manufacturer or its authorised representative immediately.

CAUTION

Do not use electrodes, adptors and extensions with visibly damaged/bent insulation as they may cause user/patient burns or even a fire!

CAUTION

Tissue may adhere to the distal end of the active electrode or it may discolour after use. Such adherent residue or discolouration is not considered a cause for complaint. In addition, due to the presence of an electric arc and mechanical force, deformation of the electrode or wear of the insulating material may occur, which are also not grounds for complaint.

Preparation for the procedure

- 1) Insulate the patient from any conductive metal elements. The operating table should be grounded. It is forbidden to place the connection cables directly on the patient, as this could lead to burns. Place the patient on a dry, electrically-insulated sheet. Make sure that the different parts of the patient's body do not touch one another. Dry gauze may be used for insulation. Fluids accumulated, e.g. in body cavities, should be removed before starting the procedure.
- 2) Fix the neutral electrode on the patient's body.

CAUTION

Refer to the instructions for use of the neutral electrode and follow the the information contained therein, especially those concerning application of the neutral electrode.

- 3) Place the electrodes in their appropriate handels (electrode adaptors can be used to match the diameter or electrode extensions can be used if applicable). Make sure that the electrode /adaptor/extension is firmly placed in the handle.
- Do this carefully to avoid damage to the accessory, handels/cables and injury to the patient or surgical staff. The electrode tip may be damaged if excessive force is applied.
- Then connect the handles cables, the neutral electrode and the footswitch to the corresponding sockets of the electrosurgical unit.
- 4) Before starting the application, make sure that the handle or cable you are using is properly connected to the electrosurgical unit and that the correct power setting is selected and displayed.
- 5) Monopolar electrodes are activated with the footswitch or by the button in a monopolar handle.

It is recommended to use a system smoke evacuation system.

6. Processing - cleaning, disinfection and sterilisation
Before their first and each subsequent use, the electrodes, adaptors and extentions should be subjected to the cleaning, disinfection and sterilisation procedure.

EMED recommends described below validated reprocessing procedures . Equivalent deviant processes are possible. It is the sole responsibility of the user to safeguard the suitability of the actual applied procedure by suitable means, e.g. validation, routine examination, verification of material compatibility, etc. In view of the design, the materials used and the intended use of the product, a maximum limit cannot be defined for the number of cleaning, disinfection or sterilization cycles that may be possible. During proper use of the electrosurgical electrodes, the products are subject to natural wear and tear depending on type and duration of their application. Damage to electrodes caused by technical factors resulting from intensive use cannot be avoided, therefore a visual inspection has to be done prior to each application (see paragraph „Visual Inspection“).

CAUTION

The electrodes, adaptors and extensions should be safely stored and transported to the area where contaminated devices are kept, in hermetic transport packaging protecting them against mechanical damage and physical, chemical and microbiological contamination.

CAUTION

Electrodes, adaptors and extensions should be processed separately.

CAUTION

Do not immerse the electrodes, adaptors and extensions in solutions containing more than 120mg/l chlorides (e.g. isotonic NaCL solution). Longer contact with the solution may result in pitting and stress corrosion.

CAUTION

When cleaning manually it is forbidden to use metal brushes, mineral wool, sanding blocks or any other coarse material.

Validated procedure for cleaning, disinfection and sterilization

If the cleaning, disinfection, sterilisation process is different from the recommended one, it must be validated for conformity to the relevant standards. Automatic cleaning and disinfection are recommended, given the low effectiveness of manual cleaning.

Preparation for cleaning:

The electrodes, adaptors and extensions should be removed from the packaging and placed in a container / device intended for cleaning / sterilization.

Manual pre-cleaning

Manual pre-cleaning of the electrodes, adpators and extensions should be carried out immediately after the end of the procedure.

Clean the product under tap water with drinking water quality, use only soft brushes (<40°C), as otherwise particles or dried discharges may adhere to them, making subsequent cleaning and sterilisation difficult or impossible. It must be ensured that areas with difficult access are thoroughly cleaned then rinsed tap water several times (<40°C). Cavities have to be rinsed intensively using 3 x 20 ml cold mains water (<40°C) with a syringe. This pre-cleaning step has to be done prior to the further manual cleaning or prior to cleaning with the cleaning and disinfection machine.

After each completed procedure, check the electrodes, adapters and electrode extenders for visible damage. Broken or damaged parts of the devices may pose a risk to the patient.

8. Automated cleaning and disinfection

CAUTION

A washer-disinfector must in principle ensure proven performance (e.g. bear the CE marking in accordance with the Standard ISO 15883).

CAUTION

Automated thermal disinfection must be carried out in accordance with national regulations in accordance with the ISO 15883-1 standard, taking into account the value of A₀: A₀-value> 3000).

CAUTION

The information provided by the manufacturer of the cleaner or disinfectant shall be followed.

Only those intended for use with medical devices made of metal and steel, with CE sign and having a pH value between 5.5 and 12.3 shall be used. A thermal disinfection programme must be chosen. The instructions for use of the agent manufacturer must be followed for the course of the programme. Make sure that in the course of cleaning accessories do not touch each other.

CAUTION

Do not clean electrodes, adaptors and extensions with sharp-edged or pointed objects.

Electrodes/adaptors/extenions should be placed in a suitable device in the basket/container intended for accessories and then Start the program with the following steps:

- 1. Pre-rinse with cold (<40°C) tap water with drinking water quality for 1 min.
- 2. Discharge
- 3. Second rinsing with cold (<40°C) water with drinking water quality for 3 min.
- 4. Discharge
- 5. Washing at 55 °C with 0.5% alkaline agent for 5 min (0.5% Neodisher® MediClean Forte is recommended).
- 6. Discharge

- 7. Neutralization with warm water (>40°C) and neutralizer (0.1% NeodisherZ) for 3 minutes.
- 8. Discharge
- 9. Rinsing with warm (>40°C) demineralised water for 2 minutes.
- 10. Discharge

Automatic disinfection:

Automatic thermal disinfection must be carried out taking into account requirements for A₀ values (see EN ISO 15883-1, A₀ value >3000) for 5 minutes at >90°C.

Automatic drying:

Takes place for 30 minutes at 90°C.

Visual inspection

At the end of the programme, a visual check for cleanliness and residual contamination should be carried out.

CAUTION

Make sure that any dirt is removed from the electrodes/ adaptors/ extensions. Improper cleaning may cause a build-up of protein deposits, and protein coagulation on disinfected instruments may interfere with sterilization.

If necessary, repeat the cleaning and disinfection process until the accesory is optically clean. Any cavities and insufficiently dried surfaces should be dried with sterile compressed air.

After drying the electrodes, adaptors and extensions shall be placed immediately in disposable paper-foil sterilisation packaging (according to EN ISO 11607-1).

9. Sterilisation

The sterilisation process compliant with the Standard ISO 17665.

Only accessories that have been previously cleaned and disinfected may be sterilized.

Electrodes, adaptors and extensions are intended exclusively for steam sterilisation in autoclaves.

EMED is not responsible for the use of other sterilization methods, e.g. in ethylene oxide.

FOR STERILISATION USE THE FRACTIONAL VACUUM METHOD:
Temperature: at least 134°C
Time: at least 5 min
Drying time – at least 10 min.

The maximum temperature must not exceed 137°C.
Do not sterilise with hot air, ethylene oxide
STERRAD and gamma radiation.

The user has sole responsibility for protecting and maintaining the sterile condition of the electrodes, adaptors and extensions after the sterilization process. If the above described and recommended chemicals and devices for manual pre-cleaning and/or automatic cleaning and disinfection are not available, the user is responsible for validating his procedure.

10. Repair and modification

It is not permissible to repair defective electrodes, adaptors and extensions. They must be replaced by new ones. Unauthorized modifications and repairs are strictly prohibited and will entail invalidation of the manufacturer's warranty.

11. Packing and storage conditions

The electrodes, adaptors and extensions must be stored in a clean and dry environment. The accessories must always be handled with the utmost care during transportation, cleaning, upkeep, sterilization and storage. This applies in particular for fine tips and other sensitive areas. It is the operator's responsibility to ensure that the sterile condition is preserved after the sterilization process.

12. Returns

Returns will be accepted only if they are marked as "hygienically safe" or "not contaminated" and have been securely packaged for shipping.

13. Disposal

The electrodes, adaptors, extensions and the packaging material of the accessories must be disposed of in accordance with the regulations and laws specific to the country in which they are used.

13. Medical incidents

Any serious incident relating to the product must be reported to the manufacturer and to the competent authority of the Member State in which the user or patient is established.

14. Update of the instructions for use

Keep the instructions throughout the working life of accessories. The instruction for use are also available in electronic form after scanning the QR code situated below the code on the packaging.

Copyright© EMED.

All rights reserved. Any copying, distribution, publishing in whole or in part without written EMED permission is prohibited.

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Copyright©EMED.