IFU: **900-023**
Version: **06.01**
Last update: **2022-10-19**

(240)900-023(10)0601

ENGLISH	
INSTRUCTIONS FOR USE Accessories for electrosurgical units Electrode handles for electrosurgery, reusable	
104-230	- Electrode handle 4mm, small, 2 switches, cable 3m, ERBE plug
104-250	- Electrode handle 4mm, small, 2 switches, cable 5m, ERBE plug
105-230	- Electrode handle 2.4mm, small, 2 switches, cable 3m, ERBE plug
105-250	- Electrode handle 2.4mm, small, 2 switches, cable 5m, ERBE plug
110-145	- Electrode handle 4mm, large, 2 switches, cable 4m, ERBE plug
111-145	- Electrode handle 2.4mm, large, 2 switches, cable 4m, ERBE plug
215-230	- Electrode handle 4mm, small, 2 switches, cable 3m, 3-pin plug
215-235	- Electrode handle 4mm, small, 2 switches, cable 3m, SDS plug
215-250	- Electrode handle 4mm, small, 2 switches, cable 5m, 3-pinplug
215-255	- Electrode handle 4mm, small, 2 switches, cable 5m, SDS plug
218-230	- Electrode handle 2.4mm, small, 2 switches, cable 3m, 3-pin plug
218-235	- Electrode handle 2.4mm, small, 2 switches, cable 3m, SDSplug
218-250	- Electrode handle 2.4mm, small, 2 switches, cable 5m, 3-pin plug
218-255	- Electrode handle 2.4mm, small, 2 switches, cable 5m, SDS plug
227-145	- Electrode handle 4mm, large, 2 switches, cable 4m, MARTIN plug
228-145	- Electrode handle 2.4mm, large, 2 switches, cable 4m, MARTIN plug
322-145	- Electrode handle 4mm, large, 2 switches, cable 4m, 3-pin plug
322-145	- Electrode handle 4mm, large, 2 switches, cable 4m, SDS plug
327-145	- Electrode handle 2.4mm, large, 2 switches, cable 4m, 3-pin plug
327-145	- Electrode handle 2.4mm, large, 2 switches, cable 4m, SDS plug

**EMED SP. Z O. O. SP. K.**
ul. Ryżowa 69A,
05-816 Opacz-Kolonia, Poland
SRN PL-MF-000012648

**2274**

Tel. +48 22 455 66 66
support@emed.pl

Reports of medical incidents:
Tel. +48 22 455 66 77
incidents@emed.pl

Used symbols:	
	Manufacturer
 2274	Conforms to Directive 93/42/EEC certified by Notified Body No. 2274
	Consult instruction for use
	Caution
	Protect from moisture
	Protect from sunlight
	Batch number
	Catalogue number
 (yellow background)	Warning of risk arising from a situation endangering the safety of the patient or user or interfering with the performance of the device
	Use to
	Safe temperature range (min. and max. values)
	Safe humidity range (min. and max. values)

The electrosurgical accessories from EMED comply with the Essential Requirements defined by Council Directive 93/42/EEC and subsequent amendments.

WARNING
Before starting to use the electrode handles read these Instructions for Use and keep them for the whole period of their use in a place accessible to the medical personnel and service staff.
The present Instructions for Use only apply to the handles and do not substitute for the instructions for use of the electrosurgical unit.
Follow the instructions for use for the electrosurgical unit provided by its manufacturer and the guidance given in it, in particular the contraindications, information and guidance for the performance of electrosurgical procedures and the safety rules for monopolar procedures.
Also read the instructions for use of the monopolar electrode to be used with the handle.

1. Intended use
They electrode handles are intended to control and directly influence the operation of monopolar modes in open electrosurgery.

1.1 Indications
In light of their general intended use, the electrode handles have no specifically defined purpose, with consideration given to clinical indications.

1.2 Contraindications
Electrosurgery is not recommended for pregnant women and persons with: implantable electronic devices (e.g. cardiac pacemakers, nerve stimulators, cardioverters or hearing implants), metal implants, arterial hypertension, diabetes and blood coagulation disorders.

These factors pose a risk of adverse events.

The contraindications include:
- contraindications in accordance with the current state of knowledge, specific of the operating procedure to be carried out,
- contraindications related to the patient's general condition,
- any other situations where the level of the risk posed by hazards to the patient's health or life would outweigh the benefits of the use of the device.

The person performing a surgical procedure is responsible for a decision to apply electrosurgery, while assessing, at the same time, the possible hazards.

Provide the patient with any information concerning the risks, adverse events or side effects related to the use of the medical device.

1.3 Target patient group
The electrode handles have no specifically defined purpose, with consideration given to the patients' clinical conditions, gender or age.

1.4 Intended users
The electrode handles may only be used by a qualified healthcare professional with formal education in the relevant field of medicine and knowledge in the field of electrosurgery. They are not intended for use by lay persons.

1.5 Use environment
The device may only be used in rooms intended for medical purposes.

2. Clinical benefits
The use of the cutting and coagulation techniques using electrode handles in open surgery procedures brings clinical benefits which have been confirmed in the literature, including, among others: shorter operative time, lesser intraoperative blood loss, shorter wound incision time, less severe postoperative pain and fewer complications.

3. Adverse events
There is a risk of burns. In order to minimise the risks occurring when electrosurgery is used, absolutely comply with the safety rules which are described in detail in the Instructions for Use.

WARNING
When procedures are carried out there is a risk of interference with or damage to implantable electronic devices (e.g. cardiac pacemakers, nerve stimulators, cardioverters or hearing implants).
In these cases, the use of the monopolar technique is not recommended. Before starting a procedure, consult a cardiologist or a relevant specialist.
It is not permitted to apply electrosurgical procedures on patients with implanted electronic devices in outpatient clinic conditions.

WARNING
Risk of ignition of the gases in the upper and lower parts of the gastrointestinal tract
Prior to activating the electrosurgical current, make sure that endogenous gases have not accumulated on the application site. The endogenous gas in the gastrointestinal tract can be flammable or explosive.
Risk of ignition of the gases in the respiratory system
Because of the risk of ignition, the following guidance must be followed:
- oxygen or other flammable gases and substances must not be fed into the respiratory system,
- if the patient needs ventilation, stop the cutting and/or coagulation (apply alternating ventilation)
Exogenous burns
They are caused by the ignition of flammable liquids or gases. The possible causes of this type of burns include, among others, the ignition of cleaning and disinfecting agents, anaesthetic gases etc.

Endogenous burns
High-frequency leakage currents can cause burns in the patient's tissues far from the electrode application sites, e.g. close to metal implants or if these sites are in contact with conductive elements, generating high current intensity. These burns can also be caused by direct contact between the conducting cables and the patient's skin.

WARNING
Do not activate the handles in the immediate vicinity of nerve structures. This may cause the unintentional stimulation of these structures.

4. Characteristics of the devices and their mode of action

New electrode handles are reusable and non-sterile devices. Before their first and each subsequent use, the handles should be subjected to the validated cleaning, disinfection and sterilisation process.

The sales unit is an individually packaged electrode handle marked with a separate REF number. The unit packaging contains the instructions for use of the device.

The operation of the device is achieved by using high-frequency electric current generated by an electrosurgical unit. The handles deliver high-frequency current from the electrosurgical unit to the active electrodes inserted into them. This current is delivered by the active electrode at the point of its contact with the patient's tissue. As a result of the current flow, heat is generated to produce a desirable electrosurgical effect. The desirable effect on the tissue depends of the type of the electrode tip and the selected mode and operating settings of the electrosurgical unit.

5. Technical specifications
The electrode handles are intended to work with high-frequency electrosurgical units from EMED equipped with monopolar cutting and coagulation modes and sockets compatible with the plugs of the handle cables, e.g. the SPECTRUM REF 100-013 electrosurgical unit.

The cutting and coagulation functions are activated, after the connection of the handle to the appropriate monopolar electrode, by means of buttons placed on the body of the handle (yellow for cutting, blue for coagulation). They require the application of the neutral electrode on the patient. The electrode will receive HF current and remove it to the electrosurgical unit, thus closing, in addition, the monopolar current circuit between the unit and the patient.

The manufacturer is not responsible for any different uses which are not consistent with the recommendations.

The handles can be connected to electrosurgical units with a non-detachable cable.

The BF/CF applied part of the electrosurgical unit used is expanded with the accessory connected to it.

WARNING
In order to reduce the risk of tissue overheating, brief current activation at short time intervals is recommended.

6. Maximum voltage
The maximum rated voltage of the devices is specified in the table below:

REF	Vp (max)
104-230, 104-250, 105-230, 105-250, 215-230, 215-235, 215-250, 215-255, 218-230, 218-235, 218-250, 218-255	4,300 Vp
110-145, 111-145, 227-145, 228-145, 322-145, 322-145, 327-145, 327-145	6,000 Vp

WARNING
Avoid too high output settings. The maximum output voltage of the used operating mode of the electrosurgical unit must not exceed the maximum permissible voltage of the connected electrode handle! The selected output power or output level should be adequate to enable the performance of the intended procedure. Both too low and too high settings of the electrosurgical unit may cause a malfunction of the device, damage to it and unintended burns. The settings should be selected according to the surgeon's experience by referring to the clinical recommendations for a given procedure and the medical practice pursued.

7. Performance of procedures

WARNING
Immediately after the purchase of a new handle, before starting a procedure and after its completion, carefully examine the condition of the handle, cable and plug checking them for any visible damage (scratched, cracked or broken elements or damaged insulation). Before each use test the operation of the handle (see the point "Operation test"). When used as intended, the electrode handles are subject to greater or lesser wear depending on the intensity and duration of their use. It is impossible to avoid damage to the handles and cable caused by the physical factors resulting from the frequency of their use.
In case visible defects or a malfunction are found, immediately replace the device by a new one and contact the manufacturer or the authorised representative of the manufacturer. Do not use damaged devices!

WARNING
Do not use cables with visible bends or damaged insulation. Cables with damaged/bent insulation or with broken conductors can cause user/patient burns or even a fire.

WARNING
When performing procedures on patients who are connected to monitoring systems (ECG), the monitoring electrodes should be placed as far away as possible from the application sites of monopolar electrodes in light of the risk of interference with these systems. Electrosurgical procedures may also cause interference with other systems. The probability of interference is increased by entangled accessory or unit cables laid close to the systems and their cables. Before starting a procedure, check if the operation of the electrosurgical unit does not interfere with the systems to which the patient is connected. This test should be carried out in accordance with the instructions for use of the electrosurgical unit, while watching the operation of the monitoring systems. In addition, the use of monitoring systems equipped with high-frequency current protections is recommended. Needle electrodes should not be used with monitoring systems in light of the risk of current concentration at the site of their application, the thermal effect of which may lead to burns.

WARNING
During a procedure and after its completion, put the electrode handle away carefully so as to avoid possible mechanical damage to the devices. Do not throw it! Do not wrap the cable around the device!

WARNING
Disconnect the handle cable by holding it by the plug. Never pull the cable directly by its cord since it may cause damage to the device.

WARNING
Cables must not be coiled, entangled with other cables or placed directly on the patient's skin, since this can cause burns. Always keep the connected device isolated from the patient and user, given the risk of inadvertent activations.

Operation test
In order to test the buttons in the handle, connect the handle to the electrosurgical unit and then switch on the unit. Check if the unit is not active in any way before any of the cutting or coagulation buttons is pressed.
If the activation of the electrosurgical unit occurs without pressing the handle button, it must not be used any more. Return the device to the supplier.
Subsequently, activate successively the yellow and blue buttons. Each time check if the correct operating mode (cutting/coagulation) is activated when a given button is pressed.
Each time an electrode is replaced in the handle make sure that it has been fully inserted into the handle socket.

Connection of the handles to the electrosurgical unit
1. First fix the neutral electrode on the patient's body.

Read the instructions for use of the neutral electrode and strictly follow the guidance contained in them, especially on the application of the neutral electrode.

2. The patient should be isolated from any conductive elements. The operating table should be grounded.
The patient should be placed on a dry, electrically insulated sheet. Prevent contact between different parts of the patient's body. Dry gauze can be used as an insulator.
Fluids accumulated, e.g. in body cavities, should be removed before starting the procedure.

3. Subsequently, connect the electrode handle to the appropriate monopolar electrode. Make sure that the electrode is sufficiently firmly seated in the handle and the hexagonal part of the electrode shaft fits the handle socket office. This must be done carefully to avoid damage to the device and injury to the patient or surgical personnel.

4. Then connect the neutral electrode and the cables of the application electrode handles to the corresponding sockets of the electrosurgical unit. Make sure that the plug of the cable is sufficiently firmly seated in the socket of the electrosurgical unit.

The connection cables must not be laid directly on the patient since this can lead to burns. The cables must not be coiled or entangled.

5. The cutting function is activated with the **yellow button in the electrode handle or the yellow pedal in the footswitch**. The coagulation function is activated with the **blue button in the electrode handle or the blue pedal in the footswitch**.

8. Warnings
General
WARNING
WARNING: Do not make any modifications to the device and do not use devices which have undergone any modifications.
EMED bears no responsibility for the use of modified devices.

WARNING
The handles are not intended for direct contact with the heart, central circulatory system and central nervous system.

WARNING
Exercise special caution when activating a handle with a connected electrode close to other metal objects (instruments and optics) in light of the risk of unintended tissue coagulation.

The recommended safe distance is at least 10mm between the electrode being activated and the other metal instruments in the operative field. It is forbidden to activate the handle when the electrode connected to it is in contact with metal objects and/or optics.

WARNING
The electrode handles may be activated only and exclusively when the electrodes are visually controlled by the operator in order to avoid an incidental burn or activation at an unintended site. High-frequency current may only be activated when the active electrodes inserted into the handle are in contact with the tissue.

WARNING
When the cutting and coagulation quality deteriorates for a given setting do not increase the power without first:
- checking the application of the neutral electrode,
- checking the correct connection of the handle cable plugs,
- making sure that the function has been activated with the correct button/pedal in the handle or the footswitch (yellow for cutting, blue for coagulation),
- checking if the insulation of the handle, cable and electrode is not damaged,
- checking if the active electrode is clean; when it is not clean this can cause a deterioration of the operational quality of the electrosurgical unit.
If the cutting or coagulation quality does not improve although the recommended steps have been taken, replace the device with a new one.

WARNING
It is recommended that flammable anaesthetics or oxidising gases, such as nitrous oxide (N₂O) or oxygen, should be avoided during procedures performed in the head or chest region unless these agents are removed by suction.

WARNING
Use non-flammable disinfectants. If the use of flammable disinfectants cannot be avoided, they should be left to evaporate before the procedure. The use of the handles close to flammable and explosive materials poses the risk of their ignition and the start of a fire.

WARNING
After each completed procedure check if the electrode is inserted into the handle and if the handles show no visible damage. Broken off or damaged parts of the handles may pose a risk to the patient.

WARNING
Research has demonstrated that the smoke arising during electrosurgical procedures can be potentially harmful to patients and the surgical team. In order to protect people and their health, EMED recommends the use of a smoke detector to remove the smoke emitted during surgical operations.

WARNING
The use of EMED accessories with accessories, cables or medical equipment from other manufacturers can result in higher electromagnetic emissions or lower interference resistance. The accessories should be safely combined with each other or with other devices in accordance with their intended use, the interface specification and compatibility information supplied by the manufacturer in the instructions for use. The manufacturer bears no responsibility for any other use/combination of medical devices which does not comply with the recommendations.

9. Reprocessing – cleaning, disinfection and sterilisation
EMED recommends the validated reprocessing procedures described below. The user bears the sole responsibility for applying a different procedure and ensuring its correctness by appropriate means, e.g. validation, routine testing, verification of material compatibility, etc.

Restrictions on reprocessing	
REFs of electrode handles	Number of reprocessing cycles
104-230, 104-250, 105-230, 105-250, 215-230, 215-235, 215-250, 215-255, 218-230, 218-235, 218-250, 218-255	150
110-145, 111-145, 227-145, 228-145, 322-145, 322-145, 327-145, 327-145	300

WARNING
The number of reprocessing cycles specified for the device must not be exceeded. EMED is not responsible for a device for which the number of reprocessing cycles has been exceeded.

Also remember that the devices are subject to natural wear depending on the type and duration of their use. For this reason check the condition of the handles and cable for possible damage and/or bent, cracked and broken conductors.

WARNING
Use only single-use devices when performing procedures on patients suffering from any disease of the group of spongiform encephalopathies (e.g. Creutzfeldt-Jakob disease). All the accessories that may have been contaminated with prions should be destroyed because of the risk of the infection of other patients. Do not reprocess or reuse them.

WARNING
Before their first and each subsequent use, the handles should be subjected to the reprocessing procedure (cleaning, disinfection and sterilisation). Incorrectly reprocessed instruments can be a source of infection.

WARNING
Before the reprocessing procedure disconnect the electrode from the handle. It is forbidden to clean, disinfect and sterilise the reusable electrode handles with connected electrodes. Before the cleaning, disinfection and sterilisation procedure disconnect the handle cable from the electrosurgical unit.

WARNING
Do not immerse the handles in solutions containing chlorides in a concentration of more than 120mg/l (e.g. isotonic NaCl solution) and hydrogen peroxide. Longer contact with the solution may result in pitting and/or stress corrosion.

WARNING
In order to avoid the coagulation of mucus, blood or other body fluids, do not soak the handles in water with a temperature of more than 40°C, alcohol, aldehyde disinfectants or antiseptics.

WARNING
Do not use ultrasonic washers to clean the handles.

WARNING
When cleaning the handles manually, do not use metal brushes, mineral wool, abrasive pads and any other coarse materials.

WARNING
Given the low effectiveness of manual cleaning, automatic cleaning and disinfection are recommended.

WARNING
If the reprocessing procedure is different from the one recommended one (as described in Section 9), it must be validated in accordance with the relevant ISO standards.

The validated methods must be strictly adhered to in order to minimise the risk of the transfer of infectious agents and other adverse impacts on the devices.

Use only CE marked agents to clean the handles. Always follow the recommendations of the manufacturer of the agent contained in its instructions. Use mildly alkaline cleaning agents designed for medical devices from plastics and metals.

The handles should be safely stored and transported to the area where contaminated devices are kept, in tight transport packaging protecting them against mechanical damage and physical, chemical and microbiological contamination.

The drying of contaminating residues (coagulation residues), such as blood or tissue fragments, makes cleaning difficult and reduces the effectiveness of the cleaning of the device. The cleaning procedure should be started not later than 1 hour after use.

Manual pre-cleaning
Manual pre-cleaning should be carried out immediately after each use, as it prevents tissues from sticking to the accessories. This pre-cleaning step should be taken prior to the further automatic cleaning and disinfection.

Manual pre-cleaning procedure
Remove any visible dirt with a soft plastic brush under running tap water. Failure to remove contaminants (e.g. blood, tissue, microbes) and cleaning agents may pose a risk for all the successive disinfection/sterilisation processes or the correct operation of the device.

Automatic cleaning and disinfection
WARNING
The washer-disinfector must, in principle, demonstrate essential performance (the CE mark in accordance with the ISO 15883 multi-part standard).

WARNING
Automatic thermal disinfection should be performed in compliance with the national regulations on the A₀ value (see EN ISO 15883-1, the A₀ value > 3,000).

In the course of the automatic cleaning and disinfection process, follow the guidance supplied by the manufacturer of the cleaning agent. Do not clean the handles together with objects having sharp edges and tips. Make sure that in the course of cleaning and disinfection the accessories do not touch each other.

The course of the process

Carefully place the handle on the instrument tray of the washer-disinfector, make sure that the cables are not entangled/bent or squeezed and do not touch other devices, and then start the programme with the following parameters:
- pre-rinsing for 2 min. with tap water at a temperature below 40°C,
- washing for 10 min. with tap water at a temperature of 55°C and the 0.5% Neodisher Mediclean Forte cleaning agent,
- neutralisation for 3 min. with tap water at a temperature of 60°C using the 0.2% Neodisher2 neutraliser,
- final rinsing for 2 minutes with tap water at a temperature below 40°C,
- thermal disinfection for 5 minutes at a temperature of 90°C using deionised water,
- drying until the device is completely dry.

During the automatic cleaning and disinfection process do not exceed a temperature of 95°C.

WARNING
After the cleaning and disinfection process is completed, carry out a visual check of the devices for cleanliness. Make sure that any contaminants have been removed from the handles. Inaccurate cleaning can cause an accumulation of protein precipitate and the process of protein coagulation on the devices subjected to sterilisation can disturb the process. If the handle still shows visible pollution, repeat the pre-cleaning and automatic cleaning/disinfection until the device is optically clean. If, after the cleaning and disinfection process, the handles show visible damage, cracks or damage to the insulation, replace them with new ones and do not reuse them.

The device does not require maintenance.

Packaging
After the cleaning and disinfection process is completed, immediately put the handle in single-use paper foil packaging intended for sterilisation, making sure that the handle cable is not entangled/bent or squeezed. The packaging must be carried out in accordance with the ISO 11607 standard on packaging for terminally sterilised medical devices.

Sterilisation
WARNING
In order to ensure the effectiveness of the process, only devices that have been previously cleaned and disinfected may be sterilised. The sterilisation process must meet the requirements of the ISO 17665 standard.

WARNING
Steam sterilisation is the only allowed validated sterilisation method. EMED bears no responsibility for the use of other sterilisation methods.
The electrode handles must be sterilised in an autoclave with a 3-times fractionated pre-vacuum procedure.

Recommended validated conditions for steam sterilisation
Temperature: 134°C
Duration: 3 min.
Drying time: 20 min

Comply with the recommendations laid down in the instructions supplied by the manufacturer of the sterilisation equipment concerning the input configuration, operation and duration of drying.
The maximum sterilisation temperature must not exceed 137°C.

Until use keep the sterilised device in sterilisation packaging. It is the user's sole responsibility to secure and maintain the sterile state of the handles after the sterilisation process. If the chemicals and equipment/tools described above and recommended for manual pre-cleaning and/or automatic cleaning and disinfection are not available, the user is responsible for validating their procedure. In addition, if the sterilisation procedure selected is different from the one described in this section, the user must validate it appropriately.

10. Conditions of use, transport and storage
Until their first use the devices should be stored in their original packaging.
Store and transport the handles in a clean and dry place, unexposed to the direct impact of sunlight. Always handle them with the utmost care during their transport, cleaning, sterilisation, use and storage. Do not exceed the use by date given on the packaging of the device.
Before a procedure wait for the device to reach room temperature.

11. Repairs
Damaged or defective handles must not be repaired and each time they have to be replaced with new ones. Any repairs must only be carried out by the authorised service of the manufacturer. Each modification, self-repair or non-compliance with these Instructions for Use exempt the manufacturer from the liability for the device.

12. Returns
Please disinfect the devices before returning them. Contaminated instruments which have been used in procedures pose a biological hazard to the service personnel. EMED has the right to refuse to accept accessories which have not been disinfected and safely packaged for shipment.
13. Disposal
In order to minimise the risk of transfer of infectious agents, immediately after their use ends, place the devices intended for disposal in marked containers for medical waste and dispose of them in compliance with the locally applicable provisions of law.

14. Guarantee
The guarantee for electrode handles covers manufacturing defects for 12 months. It does not cover the damage caused by the use of the device, since it is impossible to estimate the pace of the wear of the handles.

15. Medical incidents
Any serious incident relating to the device must be reported to the manufacturer and the competent authority of the Member State in which the user or patient reside. Report medical incidents to the e-mail address: incidents@emed.pl e-mail addresses of the national safety supervision bodies can be found on the Government's websites.
A serious incident means any incident that directly or indirectly led, might have led or might lead to any of the following:
(a) the death of a patient, user or other person,
(b) the temporary or permanent serious deterioration of a patient's, user's or other person's state of health,
(c) a serious public health threat.

16. Updates of the Instructions for Use
Keep the Instructions for Use throughout the period when the accessories are used. After each update the manufacturer publishes the new electronic version of the instructions on its website. The current Instructions for Use are also available in electronic form after the scanning of the QR code given below or the code on the packaging.

Copyright©EMED.
All rights reserved.
Any copying, distribution and publishing in whole or in part without EMED's written permission is prohibited.

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