ADP: **900-014**
Version: **04.03**
Last update:
2023-09-07

(240)900-014(10)0403

ENGLISH		
INSTRUCTIONS FOR USE: Accessories for electrosurgical units Electrosurgical cables and adaptors, reusable		
Monopolar cables		
REF	Compatible with socket	Compatible accessory
101-051	Ø 5mm, Erbe VIO / ICC / ACC	Monopolar laparoscopic instruments, Ø 4mm connector
106-04S	Ø 5mm, Erbe VIO / ICC / ACC	Wolf/Tonarra monopolar resectoscope
161-030	Ø 8mm	Monopolar endoscopic instruments, Ø 3mm connector
181-030	Ø 5mm, Erbe VIO / ICC / ACC	Monopolar endoscopic instruments 3mm
280-03S 280-05S	EMED SDS	Monopolar laparoscopic instruments 4mm
280-130 280-150	EMED 3-pin, 3-pin international	Monopolar laparoscopic instruments 4mm
281-03S	EMED SDS	Monopolar endoscopic instruments 3mm
281-130	EMED 3-pin, 3-pin international	Monopolar endoscopic instruments 3mm
360-030 360-050	Ø 8mm	Monopolar laparoscopic instruments 4mm
361-030 361-050	Olympus	Monopolar laparoscopic instruments 4mm
375-04S	Ø 8mm	Wolf/Tonarra monopolar resectoscope
405-04S	EMED SDS	Wolf/Tonarra monopolar resectoscope
405-14S	EMED 3-pin, 3-pin international	Wolf/Tonarra monopolar resectoscope
408-14S	EMED 3-pin, 3-pin international	Storz monopolar resectoscope
408-14S	EMED SDS	Storz monopolar resectoscope
409-04S	EMED SDS	Olympus monopolar resectoscope
409-14S	EMED 3-pin, 3-pin international	Olympus monopolar resectoscope
432-146	EMED 3-pin, 3-pin international	Flexible argon probe, oval connector/Luer Lock
432-46A	EMED SDSA	Flexible argon probe, oval connector
432-46S	EMED SDS	Flexible argon probe, oval connector/Luer Lock
Bipolar cables		
REF	Compatible with socket	Compatible accessory
101-000 101-001 101-100 101-101	Bipolar Erbe VIO / ICC / ACC	Bipolar forceps, bipolar laparoscopic instruments, bipolar electrodes – flat, EU type connector
287-030 287-050 287-130 287-150	Martin / Berchtold	Bipolar forceps, bipolar laparoscopic instruments, bipolar electrodes - flat, EU type connector
348-04S	EMED 2-pin, 2-pin 29mm	Storz bipolar resectoscope
348-04S	EMED SDS	Storz bipolar resectoscope
349-04S	EMED 2-pin, 2-pin 29mm	Olympus bipolar resectoscope
349-04S	EMED SDS	Olympus bipolar resectoscope
351-031 351-051 351-131 351-151	EMED 2-pin, 2-pin 29mm	Bipolar forceps, bipolar laparoscopic instruments, bipolar electrodes - flat, EU type connector
351-03S 351-05S	EMED SDS	Bipolar forceps, bipolar laparoscopic instruments, bipolar electrodes - flat, EU type connector

351-13S 351-15S	EMED SDS	Bipolar forceps, bipolar laparoscopic instruments, bipolar electrodes - flat, EU type connector
352-030 352-050	EMED 2-pin, 2-pin 29mm	Wolf bipolar laparoscopic instruments
352-03S	EMED SDS	Wolf bipolar laparoscopic instruments
354-04S	EMED 2-pin, 2-pin 29mm	Wolf/Tonarra bipolar resectoscope
354-04S	EMED SDS	Wolf/Tonarra bipolar resectoscope
357-04S	EMED 2-pin, 2-pin 29mm	Aesculap bipolar laparoscopic instruments
357-04S	EMED SDS	Aesculap bipolar laparoscopic instruments
358-030	EMED 2-pin, 2-pin 29mm	Bissinger BiTech bipolar scissors
358-03S	EMED SDS	Bissinger BiTech bipolar scissors
359-030 359-050	EMED 2-pin, 2-pin 29mm	Storz bipolar laparoscopic instruments
359-03S 359-05S	EMED SDS	Storz bipolar laparoscopic instruments
401-030 401-050	EMED 2-pin, 2-pin 29mm	ThermoStapler® bipolar vessel sealing clamps
401-03S 401-05S	EMED SDS	ThermoStapler® bipolar vessel sealing clamps
Neutral electrode cables		
REF	Compatible with socket	Compatible accessory
295-030	EU 6.3mm	Reusable neutral silicone electrode REF 242-003
295-050	EU 6.3mm	
385-030	EMED flat socket	
385-03S	US	
385-050	EMED flat socket	
385-05S	US	Single-use neutral electrodes REF 812-80H, REF 812-83H
294-030	EU 6.3mm	
294-050	EU 6.3mm	
380-030	EMED flat socket/USA	EMED flat socket/USA
380-050	EMED flat socket/USA	
Monopolar adaptors		
REF	Compatible with the socket of the unit	Compatible with the socket of the accessory
111-100	Ø 5mm, Erbe VIO / ICC / ACC	EMED 3-pin, 3-pin international
201-217	SDS	EMED 3-pin, 3-pin international
222-100	Martin / Berchtold	EMED 3-pin, 3-pin international
Bipolar adaptors		
REF	Compatible with the socket of the unit	Compatible with the socket of the accessory
201-21S	SDS	EMED 2-pin, 2-pin 29mm
281-60S	SDS	EMED 2-pin, 2-pin 29mm
340-100	Ø 5mm, Erbe VIO / ICC / ACC	EMED 2-pin, 2-pin 29mm
901-210	EMED 2-pin, 2-pin 29mm	Ø 5mm, Erbe VIO / ICC / ACC

EMED SP. Z O. O. SP. K.
ul. Ryżowa 69A,
05-816 Opacz-Kolonia,
Poland
SRN PL-MF-000012648

Tel. +48 22 455 66 66
support@emed.pl

Reports of medical incidents:
Tel. +48 22 455 66 77
incidents@emed.pl

Manufacturer

Batch code

The device conforms to the requirements of Regulation (EU) 2017/745

Catalogue number

(blue background)

Read the instructions for use

(yellow background)

Warning against a risk caused by a situation threatening the safety of the patient or user or disturbing the operation of the device

Date of manufacture

Follow the instructions for use

Protect against moisture

Protect against sunlight

Medical Device

Number of items in packaging

The electrosurgical accessories from EMED comply with the general safety and performance requirements as defined by Regulation (EU) 2017/745 of the European Parliament and of the Council (MDR).

WARNING

Before starting to use the cables read these Instructions for Use and keep them for the whole period of their use in a place accessible to the medical personnel and service staff. The present Instructions for Use do not substitute for the instructions for use of the electrosurgical unit.

Always follow the instructions for use of the electrosurgical unit provided by its manufacturer and the guidance given in it, in particular contraindications, for the performance of electrosurgical procedures, as well as safety rules for monopolar procedures using argon. This also applies to the instructions for use of accessories, including argon probes and neutral electrosurgical electrodes to be used during the monopolar operation with the the electrosurgical unit.

Specifications, safety instructions and warnings contained in the relevant instructions for use must always be complied with and followed.

1. Intended use

Electrosurgical cables are intended for the transmission of energy between an active device and a patient – without significant changes - for the purposes of monopolar or bipolar cutting or coagulation of tissues in electrosurgery. Electrosurgical cables REF 432-XXX are intended for the transmission of energy and argon gas between an active device and a patient – without significant changes - for the purposes of monopolar argon plasma coagulation of tissues in electrosurgery. Monopolar and bipolar adaptors are intended for the connection of electrosurgical units with different cable models of the same accessory and/or accessories from different manufacturers.

1.1 Indications

In light of their general intended use, the cables and adaptors have no specifically defined purpose, with consideration given to clinical indications.

1.2 Contraindications

Electrosurgery is not recommended for pregnant women and persons with: implantable electronic devices (e.g. cardiac pacemakers, nerve stimulators, cardioverters or hearing implants), metal implants, arterial hypertension, diabetes and blood coagulation disorders. These factors pose a risk of adverse events.

The contraindications include:

- contraindications in accordance with the current state of knowledge, specific of the operating procedure to be carried out,
- contraindications related to the patient's general condition,
- any other situations where the level of the risk posed by hazards to the patient's health or life would outweigh the benefits of the use of the device.

The person performing a surgical procedure is responsible for a decision to apply electrosurgery, while assessing, at the same time, the possible hazards.

Provide the patient with any information concerning the risks, adverse events or side effects related to the use of the medical device.

1.3 Target patient group

Electrosurgical cables and adaptors have no specifically defined purpose, with consideration given to the patients' clinical conditions, gender or age.

1.4 Intended users

Electrosurgical cables and adaptors may only be used by a qualified healthcare professional with formal education in the relevant field of medicine and knowledge in the field of electrosurgery. They are not intended for use by lay persons.

1.5 Use environment

The devices may only be used in rooms intended for medical purposes.

2. Clinical benefits

The clinical benefits brought by electrosurgical cables and their adaptors depend on the accessory connected to them, since these accessories may be used as intended only when they are connected to a cable.

3. Adverse events

There is a risk of burns. In order to minimise the risks occurring when electrosurgery is used, absolutely comply with the safety rules which are described in detail in the Instructions for Use.

4. Characteristics of the products and their mode of action

Electrosurgical cables are used to transmit electric current or argon gas along a specific path without significant changes in energy or substance. They are products designed to transmit energy between an active device such as an electrosurgical unit generating high-frequency current and the active electrode used in electrosurgery. Argon cables transmit the gas in addition to current. The neutral electrode cables are intended for the transmission of electrical energy from the neutral electrode to the electrosurgical unit. Adaptors enable the accessory cable plug to be fitted to the appropriate socket of the electrosurgical unit.

5. Technical specifications

Electrosurgical cables work with high-frequency electrosurgical units which are equipped with monopolar or bipolar cutting and coagulation modes and sockets compatible with the cable plugs. They operate as intended after a compatible active accessory is connected. Monopolar cables require the application of the neutral electrode on the patient, to receive HF current and return it to the unit, thus closing the monopolar current circuit between the unit and the patient. Active accessories connected by a monopolar/bipolar cable with the electrosurgical unit are activated by the footswitch pedals (the yellow one for cutting and the blue one for coagulation).

The manufacturer bears no responsibility for any other use inconsistent with the recommendations.

Maximum voltage of the products:	
Name	Vp (max)
Monopolar cables and monopolar adaptors	6000 Vp
Bipolar cables and bipolar adaptors	1200 Vp
Neutral electrode cables	500 Vp

The BF/CF applied part of the electrosurgical unit used is expanded with the accessory connected to it.

WARNING

Always set the correct power for a given accessory.

Avoid too high output settings. Check the maximum rated voltage of the accessory in the instructions for use or on the label. In case of doubt, contact the manufacturer. The maximum output voltage of the used operating mode of the electrosurgical unit must not exceed the maximum permissible voltage of the connected adaptor, the electrosurgical cable and the accessory connected to it. Both too low and too high settings of the electrosurgical unit may cause a malfunction of the product, damage to it and unintended burns. The settings should be selected on the basis of the surgeon's experience, by referring to the clinical recommendations for a given procedure and the medical practice pursued.

When the quality of cutting and coagulation deteriorates for a given setting do not increase the power without first:

- checking the application of the neutral electrode in the monopolar mode,
- checking the correct connection of the cable plugs, adaptors and accessories,
- making sure that the function has been activated with the appropriate footswitch pedal (the yellow one for cutting and the blue one for coagulation),
- checking if the insulation of the cables, adaptors and accessories is not damaged,
- checking if the active electrode is clean; if it is not clean this can cause a deterioration of the performance of the electrosurgical unit.

If the quality of cutting and coagulation does not improve although the recommended actions have been taken, immediately replace the product with a new one.

6. Warnings

General

WARNING

Do not make any modifications to the device and do not use devices which have undergone any modifications.

EMED bears no responsibility for the use of modified products.

WARNING

When procedures are carried out there is a risk of interference with or damage to metal implants and implantable electronic devices (e.g. cardiac pacemakers, nerve stimulators, cardioverters or hearing implants). In such cases, the use of the monopolar technique is not recommended. Before starting a procedure, consult a cardiologist or a relevant specialist. It is not permitted to use electrosurgical procedures on patients with implanted electronic devices in outpatient clinic conditions.

WARNING

Electrosurgical cables must not be placed close to the cables of monitoring devices, run together, parallel to the camera cables or be arranged in a loop. This can cause a distortion of the image of monitoring devices. Before starting a procedure, check if the operation of the unit does not interfere with the systems to which the patient is connected. This test should be carried out in accordance with the instructions for use of the electrosurgical unit, while watching the operation of the monitoring systems. In addition, the use of monitoring systems equipped with high-frequency current protections is recommended. Needle electrodes should not be used with monitoring systems in light of the risk of current concentration at the site of their application, the thermal effect of which may lead to burns.

WARNING

Do not place accessories and cables directly on the patient's skin. Always keep the connected product insulated from the patient and user. Unintended activation of the accessory can cause patient burns.

WARNING

Disconnect electrosurgical cables from the unit by holding them by the plug. Never pull directly by the cord of the cable, since this can damage it.

WARNING

Throughout the procedure make sure that there are no flammable substances (anesthetics, oxidising gases, endogenous gases, etc.) in the immediate vicinity of the product; otherwise there is a danger of explosion.

WARNING

Use non-flammable disinfectants; do not use e.g. alcohol-based or similar fluids. All oxygen compounds must be secured against leak during the procedure.

WARNING

Exercise special caution to ensure that the cable plug/adaptor is not flooded/wetted! In case the plug is flooded, immediately stop the activation, disconnect the cable/adaptor from the electrosurgical unit and replace the accessory by a new one. Flooded/wetted cables/adaptors must not be used again until they dry up!

WARNING

In order to avoid the risk of injury and electric shock to the patient and medical staff, make sure that the power is switched off before connecting or disconnecting the adapter/cable to/from the electrosurgical unit.

WARNING

The use of EMED accessories with accessories, cables or medical devices from other manufacturers may result in increased electromagnetic emissions or reduced immunity to interference. Accessories should be safely combined with each other and with other devices according to their intended use, interface specifications and compatibility information provided by the manufacturer in the Instructions for Use. The manufacturer bears no responsibility for any other use/combination of medical devices inconsistent with the recommendations.

⚠ WARNING
Use protective gloves when using and cleaning cables and cable adaptors.

7. Performance of procedures

⚠ WARNING
Immediately after the purchase of a new product, before starting a procedure and after its completion, carefully examine the condition of the product, its insulation and the plug, paying special attention to any damage/cracks on its surface. Damaged conductors can cause user/patient burns or even start a fire! Do not bend the cables! Do not throw them! Protect cables/adaptors against mechanical damage.

Before each use, test the operation of the product. Monopolar/bipolar cables should be tested after the electrochirurgical cable and accessory have been connected to the electrosurgical unit. Press the blue pedal of the footswitch and then its yellow pedal and check if it activates the corresponding function of the unit.

After the connection of the neutral electrode cable to the unit, connected to the corresponding neutral electrode correctly placed on the patient's body, check if the unit displays information that the connection is effective (the green LED at the neutral electrode socket or a message on the unit screen – with the neutral electrode icon illuminated in green).

If visible defects or a malfunction of the cable/adaptor are found, immediately replace the product with a new one and contact the manufacturer or its authorised representative. Damaged products must not be used!

When used as intended, the accessories are subject to greater or lesser wear depending on the intensity of their use. It is impossible to avoid damage to the cables/adaptors caused by technical factors resulting from their intensive use.

Connecting the cables to the electrosurgical unit

1) Place the patient on a dry, electrically-insulated sheet. Make sure that the different parts of the patient's body do not touch one another. Dry gauze may be used for insulation. Fluids accumulated, e.g. in body cavities, should be removed before starting the procedure. Insulate the patient from any conductive elements. The operating table should be grounded. For a monopolar procedure, first attach the neutral (passive) electrode to the patient's body. Read the instructions for use of the neutral electrode and strictly follow the guidance contained in it, especially the guidance on the application of the neutral electrode.

2) In order to connect the single-use neutral electrode cable to the appropriate neutral electrode:
- open the cable clip lever by moving it up,
- insert the tip of the single-use neutral electrode into the clip slot up to the base of the electrode,
- close the cable clip lever by pushing it down as far as it will go (the lever will be at the same level as the surface of the cable socket).
In order to connect the reusable neutral electrode to the cable, insert the electrode plug into the cable socket.

3) Monopolar/bipolar cables should be connected to a compatible accessory. Make sure that the plug of the accessory is sufficiently firmly seated in the socket of the cable. This should be done carefully to avoid damage to the accessory and cables. The accessory may be damaged if excessive force is applied.

4) After the cable REF 432-XXX is connected to an argon probe or the electrosurgical unit, check if all the components carrying the gas are firmly and tightly connected with each other. In the case of cables with 3-pin and SDS plugs, the LuerLock connector should be screwed clockwise to the argon socket of the unit. Make sure that the plug of the cable handle is sufficiently firmly seated in the socket of the unit. Before starting a procedure activate the Purge function to fill the tubes with argon. Use high purity, at least 4.8 grade (99.999%) argon.

5) Then connect the neutral electrode cable and the accessory cables to the appropriate sockets of the electrosurgical unit.

6) If necessary, attach adaptors to fit the diameter of the cable plug to the socket of the electrosurgical unit. Make sure that the cable plug/adaptor plug is sufficiently firmly seated in the socket of the electrosurgical unit. Activate the selected function - cutting or coagulation - using the footswitch.

⚠ WARNING
Do not use the neutral electrode cable (by pulling on it) to remove the neutral electrode from the patient's body.

8. Reprocessing – cleaning, disinfection and sterilisation
New electrosurgical cables and adaptors are reusable and non-sterile products.
EMED recommends the validated reprocessing procedures as described below. The user bears the sole responsibility for applying a different procedure and making sure that it is correct by suitable means, e.g. validation, routine examination, verification of material compatibility etc.

⚠ WARNING
Neutral electrode cables and adaptors are not sterile and do not require sterilization. Do not clean and disinfect them automatically or sterilise them because of the risk of damage to the product. They should be cleaned and disinfected manually using a soft, dust- and lint-free cloth.

In light of the materials used and the intended use, the maximum permissible number of cleaning and disinfection cycles for the neutral electrode cables and adaptors cannot be accurately determined. Monopolar/bipolar cables are intended for 300 sterilisation cycles; however, in the course of their correct use, they are also subject to natural wear depending on the manner and duration of their use.

Therefore, a visual inspection has to be carried out prior to each use of electrosurgical cables and adaptors. If visible defects (see "Product degradation signs") or a malfunction are found, the product should not be used and the manufacturer or its authorised representative should be contacted immediately.

Product degradation signs

- visible defects/cracks/breaks on the surface of the cable/adaptor, damage to the insulation (the cables should be carefully inspected with a magnifying glass),
- loose cable plugs,
- corrosion.

Cables should be reprocessed separately, after their disconnection from the electrosurgical unit, accessories and adaptors.

Before cleaning, adaptors should be disconnected from the electrosurgical unit and the cable.

Use only agents intended for cleaning and disinfection of medical devices made from silicone, metals and plastics. For cleaning and disinfection, use only CE marked agents which are compatible with each other. Always follow the instructions and recommendations of the manufacturer of the agent.

⚠ WARNING
Destroy all the accessories which may have been contaminated by prions (CJD). It is forbidden to use and reprocess them again! Only single-use products should be used in procedures on patients suffering from the Creutzfeldt-Jakob disease.

⚠ WARNING
In order to avoid the coagulation of mucus, blood or other body fluids, do not soak electrosurgical cables/adaptors in water with a temperature of more than 40°C, alcohol, aldehyde disinfectants or antiseptics.

⚠ WARNING
Do not immerse electrosurgical cables/adaptors in solutions containing more than 120mg/l of chlorides (e.g. isotonic NaCl solution). Longer contact with the solution may result in corrosion.

Validated procedure for reprocessing monopolar/bipolar cables

⚠ WARNING
Before first and each subsequent use, monopolar/bipolar cables should be subjected to the validated reprocessing procedure (cleaning, disinfection and sterilisation), in accordance with the approved procedure (EN ISO 17664). Incorrectly reprocessed accessories are a source of infection.

⚠ WARNING
If the procedure for reprocessing monopolar/bipolar cables is different from the recommended one (as described in Section 8), it must be validated in accordance with the relevant standards.

Monopolar/bipolar cables should be transported to the area where contaminated devices are kept, in tight transport packaging protecting them against mechanical damage and physical, chemical and microbiological contamination.

⚠ WARNING
Do not use ultrasonic washers in the process of cleaning monopolar/bipolar cables.

⚠ WARNING
Validated methods must be strictly adhered to in order to minimise the risk of transmission of infectious agents and other adverse impacts onto the products.

⚠ WARNING
Failure to remove contaminants (e.g. blood, tissues, microbes and cleaning agents) may pose a risk for all the successive disinfection/sterilisation processes or the correct operation of the product.

Automatic cleaning and disinfection are recommended, given the low effectiveness of manual cleaning.

Dried contaminating residues (coagulation deposits), such as blood or tissue fragments, hamper the cleaning of the product and reduce its effectiveness. The cleaning procedure should be started not later than within 1 hour of its use.

⚠ WARNING
When manually pre-cleaning monopolar/bipolar cables, do not use metal brushes, mineral wool, abrasive pads and any other coarse materials.

Manual pre-cleaning

Manual pre-cleaning should be carried out immediately after each use, as it prevents tissues from sticking to accessories. This pre-cleaning step should be taken prior to the further automatic cleaning and disinfection.

Remove any visible contamination with a soft, plastic brush under running tap water. In the case of cables REF 432-XXX, also rinse using a water gun all the hard-to-reach areas of cables or the inside of the argon tube.

Failure to remove contaminants (e.g. blood, tissues and microbes) may pose a risk for all the successive disinfection/sterilisation processes or the correct operation of the product.

Automatic cleaning and disinfection

⚠ WARNING
The washer-disinfector must, in principle, demonstrate essential performance (e.g. the CE mark in accordance with the PN-EN ISO 15883 standard).

⚠ WARNING
Automatic thermal disinfection should be performed in compliance with the national regulations (in accordance with the PN-EN ISO 15883 standard, taking into account the A₀ value: A₀ > 3000).

Choose a thermal disinfection programme. In the course of the programme, follow the guidance supplied by the manufacturer of the washer-disinfector and the manufacturer of the cleaning agent. Do not clean cables together with objects having sharp edges and tips.

Make sure that accessories do not come into contact with each other during cleaning.

The course of the process:

- Carefully place monopolar/bipolar cables in the instrument tray of the washer-disinfector. In the case of cables REF 432-XXX, also connect the argon tubes of cables to the rinse tubes of the device.
- Make sure the tubes are not entangled/folded or compressed and do not touch other products.
- Start the programme with the following parameters:
 - pre-rinsing** under tap water at a temperature below 40°C for 2 minutes
 - cleaning** at a temperature of 55°C for 10 minutes using tap water and a cleaning agent (0.5% Neodisher Mediclean Forte)
 - neutralisation** by rinsing for 3 minutes under tap water at a temperature of more than 60°C using a neutralising agent (0.2% Neodisher Z)
 - final rinsing** by rinsing for 2 minutes under tap water at a temperature below 40°C
 - thermal disinfection** for 5 minutes at a temperature of 90°C using demineralised water
 - drying** until they are completely dry.

In the course of the automatic cleaning/disinfection process, do not exceed a temperature of 95°C.

Visual inspection

After the programme is completed, carry out a visual inspection for cleanliness and residual contamination. Make sure that any contaminants have been removed from the products. Inaccurate cleaning can cause an accumulation of protein precipitate and the protein coagulation on the products subjected to disinfection can disturb the sterilisation process. If necessary, repeat the cleaning/disinfection process until any visible contaminants are removed.

The product does not require maintenance.

⚠ WARNING
Dry any cavities and insufficiently dried surfaces with sterile compressed air <2 bar. After drying, cables should be immediately placed in single-use paper and foil packaging for sterilisation, making sure that the cable is not entangled/bent or compressed (in accordance with EN ISO 11607-1).

Make sure that the products are dry after the cleaning and disinfection process, since water droplets promote microbial multiplication and rusting.

Sterilisation

⚠ WARNING
Only monopolar/bipolar cables which have been previously cleaned and disinfected may be sterilised. The sterilisation process must meet the requirements of the PN-EN ISO 17665-1 standard.

Please follow the recommendations laid down in the instructions of the manufacturer of the sterilisation device regarding the input configuration, handling and drying time. The maximum sterilisation temperature may not exceed 137°C.

⚠ WARNING
Steam sterilisation in autoclaves (in a fractionated pre-vacuum procedure) is the only allowed validated sterilisation method.
EMED is not responsible for the use of other sterilisation methods.

RECOMMENDED CONDITIONS FOR STEAM STERILISATION

Temperature: 134°C
Duration: 3 min.
Drying time: at least 20 min.

The sterilised device should be stored in a sterilisation package forming a sterile barrier system until use.

It is the user's sole responsibility to secure and maintain the sterile state of cables after the sterilisation process. If the chemicals and equipment described above and recommended for manual pre-cleaning and/or automatic cleaning and disinfection are not available, the user is responsible for validating its procedure. In addition, if the sterilisation procedure selected is different from the one described in the instructions, the user must validate it appropriately.

Procedure for cleaning and disinfection of neutral electrode cables and adaptors

⚠ WARNING
Follow the instructions laid down in the documentation of the manufacturer of the cleaning/disinfecting agent.

⚠ WARNING
Use a ready-to-use agent or a solution prepared according to the instructions on the leaflet of the preparation. Too high a concentration of the agent can damage the device, while too low a concentration will be ineffective.

The cleaning/disinfecting agents for neutral electrode cables should be those intended for rubber materials, too, and have a biocidal effect on bacteria, fungi and viruses.

Manual cleaning and disinfection should be carried out immediately after a procedure is completed.

The neutral electrode cable/adaptor should be disinfected again just before the next procedure. If alcohol-based disinfectants have been used, wait a moment for the alcohol to evaporate.

When manually cleaning neutral electrode cables/adaptors, do not use metal brushes, mineral wool, abrasive pads and any other coarse materials.

Cables and adapters should be cleaned using a clean, dust-free, absorbent cloth.

Operating instructions

- Prepare the agent in accordance with the manufacturer's guidance.
- Using a cloth apply the cleaning/disinfecting agent (e.g. Incidin or mikrozid® AF liquid) on the surface of the cable/adaptor. Make sure that the cable plug is not wetted/flooded.
- First clean visible dirt and make a visual inspection. If dirt/contamination is still visible, repeat the above steps until all contamination is removed, then wipe evenly the whole surface.
- Remove the excess agent from the whole surface.
- Thoroughly dry the surface of the cable/adaptor with a lint-free, absorbent and clean cloth.
- Allow the chemicals used to completely evaporate before use.

Products do not require maintenance.

Immediately after the cleaning and disinfection process, place the neutral electrode cables/adaptors in a clean, closed packaging to protect them from dirt and harmful environmental factors as recommended on the product label.

9. Conditions of use, transport and storage

Store and transport electrosurgical cables/adaptors in a clean and dry place, unexposed to the direct impact of sunlight. Always handle the products with the utmost care during their transport, cleaning, sterilisation, use and storage. The user is responsible for maintaining their sterile state after the sterilisation procedure. Until their first use the accessories should be stored in their original packaging. Before use the products must reach room temperature.

10. Repairs

Damaged or defective cables/adaptors must not be repaired and they must always be replaced with new ones. Any modification, self-repair or failure to follow the present Instructions for Use exempt the manufacturer from the liability for the product.

11. Returns

Please disinfect and sterilise the products before returning them. Contaminated products which have been used for procedures pose a biological hazard to service personnel. EMED has the right to refuse to accept accessories which have not been disinfected, sterilised and safely packaged for shipment.

12. Disposal

In order to minimise the risk of transfer of infectious agents, immediately after their use ends, place the products to be disposed of in marked containers for medical waste and dispose of them in compliance with the locally applicable provisions of law.

13. Guarantee

The guarantee for electrosurgical cables/adaptors covers manufacturing defects for a period of 12 months. It does not cover damage resulting from use of the product, since the rate of its wear cannot be estimated.

14. Medical incidents

Any serious incident relating to the device must be reported to the manufacturer and the competent authority of the Member State in which the user or patient reside. Report medical incidents to the e-mail address: incidents@emed.pl. The e-mail addresses of the national safety supervision bodies can be found on the Government's websites. A serious incident means any incident that directly or indirectly led, might have led or might lead to any of the following:
(a) the death of a patient, user or other person,
(b) the temporary or permanent serious deterioration of a patient's, user's or other person's state of health,
(c) a serious public health threat.

15. Updates of the Instructions for Use

Keep the Instructions for Use throughout the period when the accessories are used. The current Instructions for Use are also available in electronic form after the scanning of the QR code given below or the code on the packaging.

Copyright© EMED.
All rights reserved. Any copying, distribution, publication in whole or in part without written EMED permission is prohibited.

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