
Gebrauchsanweisung / Instructions for use / Istruzioni per l'uso / Mode d'emploi / Gebruiksaanwijzing /
 Instrucciones de uso / Käyttöohje / Bruksanvisning / Bruksanvisning / Инструкции за употреба / Οδηγίες χρήσης /
 Uputstvo za upotrebu / Návod k použití / Instrukcja obsługi / Instrucțiuni de utilizare / Návod na použitie /
 Navodila za uporabo / Használati utasítás / Upute za korištenje / Naudojimo instrukcijos / Lietošanas instrukcija /
 Kasutusjuhend / Инструкция по применению / Kullanım kılavuzu / Instruções de uso /取扱説明書 / 사용 설명서

DE XD 8151 Pumpenschlauch für ulrichINJECT CT motion
 EN XD 8151 Pump tubing for ulrichINJECT CT motion
 IT Tubo della pompa XD 8151 per iniettore CT motion ulrichINJECT
 FR Tubulure de pompe XD 8151 pour ulrichINJECT CT motion
 NL XD 8151 Pompslang voor ulrichINJECT CT motion
 ES Tubo de bomba XD 8151 para ulrichINJECT CT motion
 FI XD 8151 ulrichINJECT CT motion -injektorin pumpputetku
 NO XD 8151 Pumpeslange for ulrichINJECT CT motion
 SV XD 8151 Pompslang för ulrichINJECT CT motion
 BG XD 8151 Помпен маркуч за ulrichINJECT CT motion
 EL XD 8151 Σωλήνας αντλίας για τον ulrichINJECT CT motion
 SR XD 8151 Crevo za pumpu za ulrichINJECT CT motion
 CS XD 8151 Hadička čerpadla pro ulrichINJECT CT motion
 PL Wężyc pompy XD 8151 do wstrzykiwacza ulrichINJECT CT motion
 RO XD 8151 Tubulatură pompă pentru ulrichINJECT CT motion
 SK Hadička čerpadla XD 8151 pre injektor ulrichINJECT CT motion
 SL XD 8151 Cevke črpalke za pripomoček ulrichINJECT CT motion
 HU XD 8151 Szivattyúszerezék az ulrichINJECT CT motionhöz
 HR Cijev crpke XD 8151 za ulrichINJECT CT motion
 LT XD 8151 siurblio vamzdelis, skirtas „ulrichINJECT CT motion“
 LV XD 8151 sūkņa caurulīte ierīcei „ulrichINJECT CT motion“
 ET XD 8151 Pumbavoolik ulrichINJECT CT motioni jaoks
 RU XD 8151 Трубка насоса для ulrichINJECT CT motion
 TR XD 8151 ulrichINJECT CT motion için pompa hortumu
 BR Tubo flexível de bomba XD 8151 para ulrichINJECT CT motion
 JA ulrichINJECT CT motion専用XD 8151ポンプチュービング
 KO ulrichINJECT CT motion용 XD 8151 펌프 튜브

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EN Instructions for use

Abbreviations

CT motion	ulrichINJECT CT motion
CT	Computed tomography
MRI	Magnetic resonance imaging
PET	Positron Emission Tomography
CM	Contrast media
NaCl	Sodium chloride (0.9% physiological saline solution)
IFU	Instructions for use

About this document

These IFU are part of the product and they describe its safe and proper use.

- Read the instructions for use before using the product.
- Keep the instructions for use readily available with the product at all times.

Labeling and symbols

Symbol	Description
WARNING	Warning! Hazards to persons. Non-compliance can cause death or serious injuries to persons, or damage to the product.
→	Handling instructions
	CE mark with identification number of the notified body
	Caution
	Medical device
	Catalogue number
	Batch code
	Quantity
	Do not re-use! Intended for single use only for 24 hours.
	Use-by date
	Manufacturer
	Date of manufacture
	Consult instructions for use
	Sterilized using ethylene oxide
	Do not use if package is damaged and consult instructions for use
	Single sterile barrier system
	Do not use if protective caps are not properly attached
	Non-pyrogenic
	Humidity limitation, non-condensing
	Keep away from sunlight
	Temperature limit
	Pump tubing
	Keep dry

Purpose of use

The pump tubing for ulrich medical CT/MRI injectors is intended for the venous administration of contrast media (CM) and physiological saline solution (NaCl) during a CT/MRI examination. The pump tubing connects the injector (on the device side) to the patient tubing.

Indications

CT and PET-CT examinations with contrast media.

Contraindications

The pump tubing for ulrichINJECT CT motion is not intended for the administration of contrast media during high-pressure angiography or other applications that do not comply with the intended use.

Side effects

No side effects are known in association with the pump tubing. The residual risks defined in the risk management are covered by the safety notes in the instructions for use.

User profile

The pump tubing may only be used by qualified medical personnel (e.g. medical laboratory assistants, medical-technical radiology assistants, medical assistants and radiologists) and after having completed training in accordance with the instructions for use. Training is performed by ulrich medical or by a person authorized by ulrich medical only as part of the initial product training or of product updates requiring training. Users already trained are allowed to train other users of the same health care facility in compliance with the internal regulations of the respective facility as well as the national regulations.

Patient group

No restrictions. Patients receiving a contrast medium enhanced computed tomography scan.

Product description

The pump tubing on the patient side consists of a tube with a male Luer Lock connector with an integrated check valve (7) and a protective cap. The check valve opens only when medium flows in the defined direction. In the event of counter-pressure, the valve remains closed. The pressure sensor unit (6) contains an integrated particle filter. The particle filter holds back impurities and the pressure sensor unit controls the pressure in the pump tubing. On the device side, the pump tubing has three spikes (1–3) with an integrated air filter in the air supply duct.

- A Device side with pump tubing
- B Patient side with patient tubing
- 1, 3 Spike for contrast medium supply
- 2 Spike for media supply for NaCl
- 4 W-piece
- 5 Tubing segment for the roll pump
- 6 Pressure sensor unit with integrated particle filter
- 7 Male Luer lock with integrated check valve
- 8 Protective caps for tubing ends
- 9 Media flow direction

Warnings and precautionary measures

 **WARNING! Risk of injury by infection, contamination, material fatigue or use of tubing extensions and non-approved products!**

- The product is EO sterilized and delivered in sterile packaging.
- Handle all sterile components carefully in order to ensure sterility.
 - Do not re-sterilize or re-use the product.
 - Do not use any product from open or damaged sterile packaging.
 - Before use, visually check the packaging for damage and the sterile barrier for physical integrity. Do not use the product if the packaging or the sterile barrier is damaged.
 - Dispose of patient tubing upon disconnecting it from the patient in order to prevent contamination of the pump tubing and the environment.
 - Only use disposables approved for single use once and do not reprocess them. If disposables are reused, there is a risk of contamination and a hazard due to material fatigue. Improper handling jeopardizes user and patient safety.
 - Do not use if the protective caps are missing or are not properly attached.
 - Do not use the product after the expiry date.
 - Do not use product that was not stored under the recommended storage conditions.
 - Do not use tubing extensions (e.g. "Heidelberg extension").
 - Only use original products approved by ulrich medical.
 - Replace and dispose of the pump tubing at the latest after 24 hours.
 - After disconnecting the patient tubing, the system must immediately be equipped with new patient tubing, and the protective cap must not be removed from the Luer lock connector until it is connected to the patient.

 **WARNING! Risk of injury due to mixing of contrast media when using the Tandem software option!**

- When a different type of contrast medium is attached to the injector that is incompatible with the previously used medium: Change pump tubing.

Safe handling

- Before using the product, check that it is operable and in proper condition. Do not use the product if any parts are damaged.
- To avoid damage due to improper use, use the product in accordance with these instructions for use and the CT motion instructions for use.
- Product and accessories must only be used by persons who have the required knowledge, training or experience.
- Ensure that all connectors are correctly connected to each other. DO NOT OVERTIGHTEN, because this can break the threads of the Luer lock connector and lead to ruptures.

Reporting incidents

All complaints in connection with the safety, efficacy or performance of the product must be reported by the user to ulrich medical (complaint@ulrichmedical.com) or the local distributor within 72 hours. If one or more components of the system should have a malfunction (i. e. one or more of the performance specifications are not met or the intended performance is not achieved for other reasons) or if it is suspected that this may be the case, ulrich medical (complaint@ulrichmedical.com) or the local distributor must be notified immediately (within 24 hours). If one or more components of the system has ever failed and possibly led or contributed to the death or severe injury of a patient, user or third party, then ulrich medical or the local distributor must be notified by telephone immediately. Please use the Medical Device Vigilance Report Form from ulrich medical to report product complaints, malfunctions or product failure. It is available at: www.ulrichmedical.com/vigilancereport.

Use

- Remove the product from the packaging and check for damaged parts before use.
- Insert the pump tubing in the injector in accordance with the CT motion instructions for use.
- Use the product in accordance with these instructions for use and the CT motion instructions for use.
- Only use the pump tubing with ulrich medical original products:

Prod. number	Name
XD 8000	ulrichINJECT CT motion
XD 2035	Patient tubing 150 cm, 2 check valves
XD 2040	Patient tubing 250 cm, 2 check valves
XD 2045	Patient tubing 320 cm, 2 check valves

Technical information

Parameter	Value
Pressure tightness	Max. 20 bar
Material	PVC (DEHP free, latex free)
Max. service life	Any number of injections within a 24 hour period

Storage

- Store the product in the sterile packaging in a dry, dark and temperature-controlled room.
- Protect the sterile packaged product from environmental factors such as dirt, dust, humidity, heat, and sunlight.

Ambient conditions

	Use	Storage
Temperature	+5°C to +40°C	+15°C to +30°C
Relative humidity	10% to 90% non-condensing	30% to 60% non-condensing

Disposal

- Attach the protective caps to the spikes.
- After use, dispose of the product and packaging in accordance with the local applicable regulations for disposable products contaminated with blood. Of course, the system and system components can also be disposed of via the clinic's own disposal system. The disposal of medical waste in the European Union is subject to the European List of Waste Materials. If necessary, various levels of government enact their own waste regulations, which must be observed. Outside the European Union, the relevant country-specific regulations for the disposal of medical waste must be observed.

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Ulrich GmbH & Co. KG | Buchbrunnenweg 12 | 89081 Ulm | Germany