Ochrona zgodna z normami

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PL

OxyChem C310

Protection Coverall, Cat. III, Type 4, 5, 6. Medical product Class I

Certification

Type 6 - Protective clothing against liquid chemicals - low-pressure splashes	Type 5 - Protective clothing against particulates	Type 4 - Protective clothing against liquid chemicals - Spray liquid	Protective clothing against radioactive contamination	Electrostatic properties	Protection against infectious agents	Protective clothing - General requirements
EN 13034:2005 +A1:2009	EN ISO 13982-1:2004 +A1:2010	EN14605:2005 +A1:2009	EN 1073-2:2002	EN 1149-5:2018	EN 14126:2003 +AC:2004	EN ISO 13688:2013

Lightweight, resistant, with a breathable panel

The **OxyChem C310** chemical protective suit is designed for protection in areas exposed to infectious material and provides limited chemical protection type 4/5/6. It is also a component of the Individual Biological Protection Coverall (BOP), created in accordance with the WHO recommendation. It is made of fabric laminated with microporous film (**MPFL**) of 63g/m² +/- 2g/m². Outer layer: polyethylene film. Inner layer: polypropylene fibres, which is resistant to infectious agents penetration. Reinforced and additionally taped seams provide Type 4 protection. **OxyChem C310** also provides protection against radioactive contamination and has anti-electrostatic properties. Its construction is designed to provide the wearer with the highest level of safety, ergonomics and comfort during work.

Properties

- Material laminated with microporous film (**MPFL**) 63g/m² +/- 2g/m²
- Additionally, taped seams
- Three-panel hood
- Two-sided zipper with self-adhesive flap
- Two-piece gusset in crotch
- Elastic bands in cuffs, legs, waist and hood
- Durable seam overlock sewn on the inside
- Elastic thumb loop
- Antistatic properties

Aplication for jobs

- In the pharmaceutical industry
- In contact with toxic materials
- When securing crime scenes
- In veterinary medicine
- In hospitals
- When working with sick people, e.g. during a pandemic
- When working with asbestos
- When painting and spray-painting

Packing

The suit is individually wrapped in a sealed foil package.

Quantity in a carton - 50 pcs.

Quantity of overalls on a pallet - 1000 pcs.

Quantity of cartons on a pallet - 20 pcs.

Gross weight of cardboard - 13 kg

Dimensions of the cardboard - 48 cm x 28 cm x 60 cm

Technical parameters

Type of protection obtained according to EN 14126 (Protective clothing - Requirements and test methods for protective clothing against infectious agents) - Performance levels.

Standard	Title	Retrieved from type protection
EN 13034:2005+A1:2009	Protective clothing against liquid chemicals - Performance requirements for chemical protective clothing offering limited protective performance against liquid chemicals (Type 6 and Type PB[6] equipment)	Type 6B
EN ISO 13982-1:2004 + A1:2010	Protective clothing for use against solid particulates - Part 1: Performance requirements for chemical protective clothing providing protection to the full body against airborne solid particulates (Type 5 clothing)	Type 5B
EN 14605:2005+A1:2009	Protective clothing against liquid chemicals - Performance requirements for clothing with liquid-tight (Type 3) or spray-tight (Type 4) connections, including items providing protection to parts of the body only (Types PB [3] and PB [4])	Type 4B

Performed examinations and tests in accordance with EN 14126 (test for protection against infectious/infectious agents):

Standard test	Result obtained	Classification by result
Pressure test, resistance to penetration of blood and body fluids with synthetic blood (EN ISO 16603 + EN 14126 par. 4.1.4.1)	20kPa	Class 6
Resistance to pathogen penetration from blood using bacteriophage (Fag Fi-X174) (resistance to viral penetration)(EN ISO 16604 + EN 14126 par. 4.1.4.1)	> 75 min	Class 6
Resistance to bacterial penetration (EN ISO 22610 + 14126 par. 4.1.4.2)		
Resistance to penetration of biologically contaminated aerosols (EN ISO 22611 + EN 14126 par. 4.1.4.3)	> 5	Class 3
Resistance to penetration of contaminated dust (EN ISO 22612 + 14126 par. 4.1.4.4)	≤ 1	Class 3

Performed examinations and tests in accordance with EN 1073-2 (protection against contamination with radioactive particles):

Standard test	Result obtained	Classification by result
Total internal leakage (EN ISO 13982-2 + EN 1073-2)	TILE < 3% TILA < 2%	Class 2

Parameters relating to the strength of the garment (physical properties)

Standard test	Result obtained	Classification by result
Determination of resistance to abrasion of materials (EN 530 met. 2 + EN 14325 par. 4.4.1 + EN 13034 par. 4.1)	> 1000	Class 4
Determination of resistance to flexural fracture (EN ISO 7854 met. B + EN 14325 par. 4.5 + EN ISO 13982-1 par. 4.1)	> 100 000	Class 4
Determination of trapezoidal tearing strength (EN ISO 9073-4 + EN 14325 par. 4.7 + EN 13034 par. 4.1)	Longitudinal > 60 Transverse > 40	Class 4/3
Determination of resistance to push-out (EN ISO 13938-1 + EC 1 +EN 14126 par. 4.1.2)	> 320	Class 4
Tensile strength and elongation. Strip method (EN ISO 13934-1 + EN 14325 par. 4.9 + 13034 par. 4.1)	Longitudinal > 100 Transverse > 60	Class 3/2
Puncture resistance of materials (EN 863 + 1073-2 par. 4.2)	> 50	Class 3
Electrostatic properties test Surface resistivity (EN 1149-1 + EN 1149-5 par. 4.2.1)	< 2,5 x 10 ⁹ Ω	
Seam strength (max. breaking load) using the grab method (EN 13935-2 + EN 14325 para. 5.5 + EN ISO 13982-1 para. 4.2.2)	> 100	Class 3
Flammability tests (ignition and flammability) (13274-4 met 3 + EN 1073-2 par. 4.2)	Test passed	

Examinations and tests according to EN 13034

Standard test	Result obtained
Test of resistance against chemicals (EN ISO EN ISO 6530 + EN 14325 par. 4.12 and 4.13 + EN 13034 par 4.1)	H ₂ SO ₄ (30%) < 1 NaOH (10%) < 1 o-Ksylen < 1 N-butanol < 1

Examinations and tests according to EN 14605

Standard test	Result obtained
Resistance to penetration by sprayed liquid (Spray test) (EN ISO 17491-4 + EN 14605 par. 4.3.4)	Test passed
Requirements and test methods for non-ventilated protective clothing against particulate contamination (EN 1073-2)	Test passed

Additional research:

Standard test	Result obtained
Determination of pH of aqueous extracts (mean PH value of extracts 20 and 30) (EN ISO 3071 + EN ISO 13688 para. 4.2)	3,5 > pH > 9,5
Cytotoxicity tests (Mean % cell survival after incubation) (EN ISO 10993-5)	94.47% - is not cytotoxic
Dermatological test (Regulation (EC) No 1223/2009 of the European Parliament and of the Council of 30 November 2009 on cosmetic products) (irritation index X _{mn})	Non-irritant product X _{mn} = 0

Assessment of Declared Product Properties - Protective Suit OxyChem C210/C310 (by J.S. Hamilton Poland Sp. z O.O. - No. 99014/21/CGDA of 22.03.2021)

On the basis of the results of the conducted performance test, the properties of the product COMPOSITES OXYCHEM C210/C310 declared by the manufacturer were evaluated as follows:

The product allows for easy dressing and removal.	100% positive responses
The product is convenient to use.	96% positive responses
The product does not put undue pressure on the skin (welts).	92% positive responses
The product does not cause reddening at the contact with the skin.	92% positive responses
The product provides good thermal comfort while working.	100% positive responses
The product is airtight.	100% positive responses
The product adheres well at the cuffs, legs, hood.	96% positive responses
The product allows freedom of movement while working.	100% positive responses

In addition, it was assessed:

Proper application of the product.	100% positive responses
Suitable product texture.	96% positive responses
Overall product rating.	84% positive responses
Intention to purchase the product.	88% positive responses
Willingness to recommend the product to a friend.	88% positive responses

EN ISO 15223-1:2021	Medical devices – Symbols to be used with information to be supplied by the manufacturer – Part 1: General requirements
EN ISO 20417:2021	Medical devices – Information to be supplied by the manufacturer.
EN ISO 14971:2019	Medical devices. Application of risk management to medical devices.

The EU Declaration of Conformity is available at: <https://www.oxyline.eu/deklaracje-kombinezony.html>

The OxyChem C310 one-piece protective suit is a Class I, non-sterile medical device which complies with the requirements of Regulation 2017/745 of the European Parliament and of the Council of 5 April 2017 on medical devices, amending Directive 2001/83/EC, Regulation (EC) No 178/2002 and Regulation (EC) No 1223/2009 and repealing Council Directives 90/385/EEC and 93/42/EEC (Text with EEA relevance.) and complies with the Medical Devices Act of 20.05.2010. It is also a Category III PPE in accordance with Regulation EU 2016/425 of the European Parliament and of the Council of 9 March 2016 on personal protective equipment and repealing Council Directive 89/686/EEC (Text with EEA relevance). Product designed and placed on the market in accordance with Regulation (EU) 2017/745 of the European Parliament and of the Council of 5 April 2017 and Regulation EU 2016/425 of the European Parliament and of the Council of 9 March 2016.

Sizing	S	M	L	XL	XXL	XXXL
Height	164 - 170	170 - 176	176 - 182	182 - 188	189 - 194	194 - 200
Chest	84 - 92	92 - 100	100 - 108	108 - 117	116 - 124	124 - 132

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Eigenschaften

- Material laminated with microporous film (**MPFL**) $63\text{g/m}^2 \pm 2\text{g/m}^2$
- Additionally, taped seams
- Three-panel hood
- Two-sided zipper with self-adhesive flap
- Two-piece gusset in crotch
- Elastic bands in cuffs, legs, waist and hood
- Durable seam overlock sewn on the inside
- Elastic thumb loop
- Antistatic properties

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