

Single-use

StellaClip Hemoclips

Usage

For endoscopic hemostasis and clamping the wound in the GI tract.

Features

- Synchronous rotation supports accurate positioning.
- Repositionable clips.
- The handle's security lock system helps prevent unexpected pain for the patient.
- MR conditional.
- Packaging: boxes of 5.



Description	Model name	Item number	Maximum outer diameter [mm]	Working length [mm]	Maximum opening of clip [mm]	Compatible endoscopes channel diameter* [mm]
Uncoated	VDK-CD-11-230-O	500063	2.7 mm	2300	11	≥ Ø 2.8
	VDK-CD-11-230-P	500064	2.7 mm	2300	11	≥ Ø 2.8
Coated	VDK-CD-13-230-P	500012	2.7 mm	2300	13	≥ Ø 2.8
	VDK-CD-16-230-P	500065	2.7 mm	2300	16	≥ Ø 2.8

*Figures are based on the validation tests done by PENTAX Medical

Single-use

StellaNeedle Injection needles

Usage

Endoscopic injection treatment of esophageal and gastric varices.
Endoscopic injection of submucosa in the GI tract.

Features

- Classic handle designed to easily control the needle movement.
- 3-bevel tip designed for precise injection and prevents lumen from clogging.
- Needle shield designed to prevent damage to endoscope working channel.
- Packaging: boxes of 10.



Needle size	Model name	Item number	Outer diameter [mm]	Working length [mm]	Needle length [mm]	Compatible endoscopes channel diameter* [mm]
19G	VDK-IN-23-230-1904-A	500016	2.3	2300	4	≥ Ø 2.8
	VDK-IN-18-120-2205-A	500013	1.8	1200	5	≥ Ø 2.0
22G	VDK-IN-23-180-2204-A	500068	2.3	1800	4	≥ Ø 2.8
	VDK-IN-23-230-2205-A	500069	2.3	2300	5	≥ Ø 2.8
25G	VDK-IN-23-230-2504-A	500017	2.3	2300	4	≥ Ø 2.8
	VDK-IN-23-230-2506-A	500019	2.3	2300	6	≥ Ø 2.8

*Figures are based on the validation tests done by PENTAX Medical