

# Working Elements

for Two-Stem Electrodes with Stabilizers, 26 Fr.




For use with Resectoscope Sheaths 26040 SL, 26050 SC and 26050 SL and 4 mm HOPKINS® Telescope 12° 26105 FA




- Special Features:
- Resection in saline solution
  - High level of safety due to direct current return via the electrode
  - Very deep coagulation effect

Cutting by means of a spring  
Movable thumb support  
In rest position the electrode tip is inside the sheath.



26040 EBH    **Working Element Set, bipolar**  
including:  
**Working Element, bipolar**  
2x **Cutting Loop, bipolar**  
**Cutting Electrode, bipolar, pointed**  
**Coagulation Electrode HALF MOON®, bipolar, with ball end**  
**Bipolar High Frequency Cord**  
**Protection Tube**

**26040 GP1**    Pętla tnąca, bipolarna,  
26 Charr., średnica płaszczka 8 mm, kolor  
kodu: żółty




Working End	26 Fr., Sheath Diameter 8 mm color code: yellow	Description
	26040 GP1	Cutting Loop, bipolar
	26040 GD1	Cutting Loop, bipolar, small
	26040 BL1	Cutting Electrode, bipolar, pointed
	26040 NB1	Coagulation Electrode HALF MOON®, bipolar, ball end
	26040 JB1	Cutting Loop, bipolar, longitudinal, color code: yellow-orange

Units and Accessories for Intrauterine HF Electrosurgery see chapter 11, UNITS  
Components/Spare Parts see chapter 12

poz. 2

4-15

Reusable Instruments for Intrauterine HF Surgery (26 Fr.) – for use with AUTOCON® III 400



26040 EBH    **Working Element Set**, bipolar, including:  
                  **Working Element**, bipolar  
                  2x **Cutting Loops**, bipolar  
                  **Cutting Electrode**, bipolar, pointed  
                  **Coagulation Electrode HALF MOON®**, with ball end  
                  **Bipolar High Frequency Cord**  
                  **Protection Tube**

26105 FA    **HOPKINS® Telescope 12°**, enlarged view, diameter 4 mm, length 30 cm, **autoclavable**






26050 SC    **Resectoscope Sheath**, including connecting tube for in- and outflow, 26 Fr., oblique beak, rotatable inner sheath

26040 OC    **Standard Obturator**

The new AUTOCON® III 400



**26040NB** Elektroda koagulacyjna HALF MOON, bipolarna, kulka, 26 Fr.

	Order No.	Description
	26040 GP1	<b>Cutting Loop</b> , bipolar
	26040 NB	<b>Coagulation Electrode</b> , HALF MOON®, bipolar
	26040 BL1	<b>Cutting Electrode</b> , bipolar, pointed
	26040 JB1	<b>Cutting Loop</b> , bipolar, longitudinal
	280	<b>Protection Tube</b> , for sterilization and storage of electrodes, loops, curettes and knives

poz. 1

# Tubing Sets

for irrigation and suction



## Special Features:

- Automatic product recognition and mode selection via the pump

STERILE

For use with ENDOMAT® SELECT



031523-10



030647-10

poz. 3

031523-10

**Tubing Set for Irrigation, PC, sterile, for single use, package of 10, for use with HAMOU® ENDOMAT® SCB 26 3311 20-1 and ENDOMAT® SELECT**

030647-10

**Tubing Set for Suction, DS, sterile, for single use, package of 10, for use with ENDOMAT® SELECT**

**031523-10** Zestaw drenów płuczących, PC, jednorazowego użytku, sterylne, opakowanie 10 szt., do użytku z HAMOU® ENDOMAT® SCB 26331120-1 i ENDOMAT® SELECT.